



DENTAL INSTRUMENTS

CATALOGUE

✉ info@simzmedicalsupplies.com

🌐 www.simzmedicalsupplies.com

About Our Company



SIMZ medical supplies Co. is a trading company that is committed to delivering quality medical supplies. SIMZ prides itself in delivering single use and reusable Dental, Surgical and Orthopedic instruments. Our products are manufactured by using the highest quality stainless steel (AISI, ASTM and DIN standards, which is imported from Germany, Japan & France). We are committed to producing and serving the best quality instruments. We are devoted to delivering quality products at affordable prices.

1 Dental Instruments

Dental instruments are tools that dental professionals use to provide dental treatment. They include tools to examine, manipulate, treat, restore and remove teeth and surrounding oral structures

2 Surgical Instruments

A surgical instrument is a tool or device performing specific actions or carrying out desired effects during a surgery or operation, such as modifying biological tissue, or to provide access for viewing it

3 Orthopedic Instruments

Orthopedic surgery instruments are designed for surgeries related to muscles, joint, ligaments and bones. Our extensive range of orthopedic instruments accommodates a variety of medical approaches

01 EXTRACTING FORCEPS **01 TO 70**

02 ROOT ELEVATORS **71 TO 113**

03 DIAGNOSTIC, PERIODONTAL & RESTORATIVE INSTRUMENTS **114 TO 224**

04 MAXILLOFACIAL SURGERY INSTRUMENTS **225 TO 276**

05 ORTHODONTIC PLIERS & INSTRUMENTS **277 TO 304**

06 DENTAL IMPLANT INSTRUMENTS **305 TO 349**

07 IMPRESSION TRAYS **350 TO 365**

08 HOLLOWARE INSTRUMENTS **366 TO 375**

Handle Selection For Extracting Forceps

Up
Right Left
Down

English Pattern English Pattern English Pattern, Anatomical Handle American Pattern Deep Grip Forceps With Hole Pattern Hole Pattern

A **B** **C** **D** **E** **F**

Extraction Forceps

Dental forceps are pliers-like instruments that are used to grasp teeth.

A dentist will usually have a number of different ones on hand, each having a design that's tailored to:

The specific shape of the tooth it's intended to remove (such as incisors vs. molars having two roots vs. molars having three roots).

The location of the tooth in the mouth (such as front vs. back, or left lower molar vs. right upper molar.)

How they're used.

A dentist will grasp a tooth with their forceps and then firmly and deliberately rock it back and forth as much as it will.

Because the bone that surrounds the tooth is compressible, the socket will gradually start to expand.

As it does, the range of the dentist's side-to-side motions will increase.

In addition to this rocking motion, the dentist will also rotate the tooth back and forth.

This twisting action helps to rip and tear the tooth away from the ligament that binds it in place.

At some point, the socket will be enlarged enough, and the ligament torn away enough, that the tooth can be easily removed.



EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)

- Щипці екстракційні (англійська форма)
- Forcèps de Extracció (Forme Inglesa)
- استخراج الملقط (نمط اللغة الإنجليزية)



SIMZ-YM1-001
Fig. 1
upper laterals
and canines



SIMZ-YM1-002
Fig. 2
upper laterals
and canines



SIMZ-YM1-003
Fig. 3
upper laterals
and canines



SIMZ-YM1-004
Fig. 4
lower laterals
and canines



SIMZ-YM1-005
Fig. 7
upper premolars



SIMZ-YM1-006
Fig. 8
lower premolars



SIMZ-YM1-007
Fig. 13
lower premolars

EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملقط (نمط اللغة الإنجليزية)



SIMZ-YM1-008
Fig. 13S
lower premolars

SIMZ-YM1-009
Fig. 17
upper molars,
right



SIMZ-YM1-010
Fig. 18
upper molars,
left



SIMZ-YM1-011
Fig. 18A
upper molars



SIMZ-YM1-012
Fig. 18B
upper molars



SIMZ-YM1-013
Fig. 19
upper third molars

EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Dавiers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці скстракційні (англійська форма)
- Forceps de Extracci6n (Forme Inglesa)
- استخراج المنقط (نمط اللغة الإنجليزية)



SIMZ-YM1-014
Fig. 20
lower third molar



SIMZ-YM1-015
Fig. 21
lower molar



SIMZ-YM1-016
Fig. 22
lower molar



SIMZ-YM1-017
Fig. 22S
lower molar,
for children



SIMZ-YM1-018
Fig. 22G
lower molar



SIMZ-YM1-019
Fig. 23
lower molar, left



SIMZ-YM1-020
Fig. 24
lower molar, right

EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Щипці екстракційні (англійська форма)
- Daviers (Forme Anglaise)
- Forceps de Extracción (Forme Inglesa)
- Zahnzangen (Englische Form)
- استخراج الملقط (نمط اللغة الإنجليزية)



SIMZ-YM1-021
Fig. 29
upper roots



SIMZ-YM1-022
Fig. 29N
Roots of upper
incisors



SIMZ-YM1-023
Fig. 29S
upper roots,
for children



SIMZ-YM1-024
Fig. 30
upper roots



SIMZ-YM1-025
Fig. 31
lower roots



SIMZ-YM1-026
Fig. 32
lower roots



SIMZ-YM1-027
Fig. 30S
upper roots,
for children

EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)

- Щипці екстракційні (англійська форма)
- Forcèps de Extracció (Forme Inglesa)
- استخراج المنطق (نمط اللغة الإنجليزية)



SIMZ -YM1-028
Fig. 33
lower roots



SIMZ-YM1-029
33A
lower roots



SIMZ-YM1-030
Fig. 33B
lower roots



SIMZ-YM1-031
Fig. 33L
lower roots



SIMZ -YM1-032
Fig. 33S
lower roots,
for children



SIMZ -YM1-033
Fig. 33C
Lower roots



SIMZ -YM1-034
Fig. 33M
Lower roots
with long beak



SIMZ -YM1-035
Fig. 35N
Upper incisors
and premolars
deep-gripping



SIMZ-YM1-036
Fig. 36N
Lower premolars
deep-gripping



SIMZ-YM1-037
Fig. 37
upper incisors,
for children

EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملقط (نمط اللغة الإنجليزية)



SIMZ-YM1-038
Fig. 38
lower incisors and
canines, for children



SIMZ-YM1-039
Fig. 39
upper molars
for children



SIMZ-YM1-040
Fig. 39 R
upper molars right
for children



SIMZ-YM1-043
Fig. 39 L
upper molars left
for children



SIMZ-YM1-041
Fig. 39 A
upper molars
for children



SIMZ-YM1-042
Fig. 40
lower molars
for children

EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forcèps de Extracció (Forme Inglesa)
- استخراج المنطق (نمط اللغة الإنجليزية)



SIMZ -YM1-044
Fig. 41
upper roots



SIMZ-YM1-045
Fig. 44
upper roots



SIMZ-YM1-046
Fig. 45
lower roots



SIMZ-YM1-047
Fig. 46
Lower root



SIMZ-YM1-048
Fig. 46N
Lower premolars
deep-gripping



SIMZ-YM1-054
Fig. 50
upper roots



SIMZ-YM1-049
Fig. 48
lower roots



SIMZ-YM1-050
Fig. 49
upper roots



SIMZ-YM1-051
Fig. 49L
Very thin lower roots



SIMZ-YM1-052
Fig. 46M
lower roots,
gripping in depth



SIMZ-YM1-053
Fig. 46M
very fine lower roots

EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Щипці екстракційні (англійська форма)
- Daviers (Forme Anglaise)
- Forceps de Extracción (Forme Inglesa)
- Zahnzangen (Englische Form)
- استخراج الملقط (نمط اللغة الإنجليزية)




SIMZ-YM1-055
Fig. 51
upper roots



SIMZ-YM1-056
Fig. 51 A
upper roots



SIMZ-YM1-057
Fig. 51 L
upper roots



SIMZ-YM1-058
Fig. 51 C
Upper roots



SIMZ-YM1-059
Fig. 51 S
upper roots,
for children



SIMZ-YM1-060
Fig. 52
upper roots



SIMZ-YM1-061
Fig. 52 A
upper teeth
and resection

EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)


- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج المنقط (نمط اللغة الإنجليزية)



SIMZ-YM1-062
Fig. 54
for separating
upper molars



SIMZ-YM1-063
Fig. 55
for separating
upper molars



SIMZ-YM1-064
Fig. 56
for separating
lower molars



SIMZ-YM1-065
Fig. 59
Lower root



SIMZ-YM1-066
Fig. 65 R
upper molars,
right



SIMZ-YM1-067
Fig. 65 L
upper molars,
left



SIMZ-YM1-068
Fig. 67
upper third molars




SIMZ-YM1-069
Fig. 67 A
upper third molars

EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- (استخراج الملقط (نمط اللغة الإنجليزية))
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)



SIMZ-YM1-070
Fig. 67½ R
upper third molars,
right.



SIMZ-YM1-071
Fig. 67½ L
upper third molars,
left



SIMZ-YM1-072
Fig. 67N
Upper third molars



SIMZ-YM1-073
Fig. 68
lower roots



SIMZ-YM1-074
Fig. 68 A
lower roots



SIMZ-YM1-075
Fig. 73
lower molars



SIMZ-YM1-076
Fig. 73 R
lower molars,
right



SIMZ-YM1-077
Fig. 73 L
lower molars,
left

EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)

- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج المنقط (نمط اللغة الإنجليزية)



SIMZ-YM1-078
Fig. 73 A
lower molars



SIMZ-YM1-079
Fig. 73 S
lower molars



SIMZ-YM1-080
Fig. 73 B
lower molars
for children



SIMZ-YM1-081
Fig. 74
lower roots



SIMZ-YM1-082
Fig. 74 D
lower roots



SIMZ-YM1-083
Fig. 73 K
lower molars
for children

EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملقط (نمط اللغة الإنجليزية)



SIMZ-YM1-084
Fig. 74M
lower roots



SIMZ-YM1-085
Fig. 74N
lower roots



SIMZ-YM1-086
Fig. 74XN
lower roots and
crowded incisors



SIMZ-YM1-087
Fig. 75
lower premolars



SIMZ-YM1-088
Fig. 76
upper roots



SIMZ-YM1-089
Fig. 76N
upper roots

EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)

- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج المنطق (نمط اللغة الإنجليزية)



SIMZ-YM1-090
Fig. 76S
upper roots




SIMZ-YM1-091
Fig. 79
lower third molars



SIMZ-YM1-092
Fig. 79A
lower third molars



SIMZ-YM1-093
Fig. 79N
Lower third molars



SIMZ-YM1-094
Fig. 79C
Lower third molars
With Serrated
Tips



SIMZ-YM1-095
Fig. 81
lower premolars
and roots



SIMZ-YM1-096
Fig. 80
upper molar roots

EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Щипці екстракційні (англійська форма)
- Daviers (Forme Anglaise)
- Forceps de Extracción (Forme Inglesa)
- Zahnzangen (Englische Form)
- استخراج الملقط (نمط اللغة الإنجليزية)



SIMZ-YM1-097
Fig. 86A
lower molars



SIMZ-YM1-098
Fig. 86B
lower molars



SIMZ-YM1-099
Fig. 86C
lower molars



SIMZ-YM1-100
Fig. 87
lower molars



SIMZ-YM1-101
Fig. 89
upper molars,
right



SIMZ-YM1-104
Fig. 90
upper molars,
left



SIMZ-YM1-102
Fig. 89
upper molars, right,
gripping in depth



SIMZ-YM1-103
Fig. 90
upper molars, left,
gripping in depth

EXTRACTING FORCEPS (English Pattern)

- Шипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)


- Шипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج المنطق (نمط اللغة الإنجليزية)



SIMZ-YM1-105
Fig. 92
upper second molars,
right



SIMZ-YM1-106
Fig. 93
upper second molars,
left



SIMZ-YM1-107
Fig. 94
upper molars,
right



SIMZ-YM1-108
Fig. 91
for separating
lower molars



SIMZ-YM1-109
Fig. 99
lower molars



SIMZ-YM1-110
Fig. 97
Very thin
upper roots



SIMZ-YM1-111
Fig. 95
upper molars,
left

EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملقط (نمط اللغة الإنجليزية)



SIMZ-YM1-112
Fig. 99 1/4
lower molars



SIMZ-YM1-113
Fig. 99 1/2
lower molars




SIMZ-YM1-116
Fig. 100
upper roots



SIMZ-YM1-117
Fig. 101
upper premolars
and roots



SIMZ-YM1-114
Fig. 233
Lower root
(parallel beaks)



SIMZ-YM1-115
Fig. 235
Upper root
(parallel beaks)

EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)

- Щипці екстракційні (англійська форма)
- Forceps de Extracció (Forme Inglesa)
- استخراج المنقط (نمط اللغة الإنجليزية)



SIMZ-YM1-118
Fig. 104
upper roots



SIMZ-YM1-119
Fig. 102
upper molars,
right



SIMZ-YM1-120
Fig. 103
upper molars,
left



SIMZ-YM1-121
Fig. 105
lower incisors



SIMZ-YM1-122
Fig. 106
lower incisors



SIMZ-YM1-123
Fig. 107
upper canines

EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملقط (نمط اللغة الإنجليزية)



SIMZ-YM1-124
Fig. 112
upper roots



SIMZ-YM1-125
Fig. 113
upper roots



SIMZ-YM1-127
Fig. 136
upper premolars
and roots



SIMZ-YM1-128
Fig. 137
lower incisors
and roots



SIMZ-YM1-126
Fig. 114
lower roots



SIMZ-YM1-129
Fig. 143
lower roots


EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)

- Щипці екстракційні (англійська форма)
- Forceps de Extracció (Forme Inglesa)
- استخراج الملقط (نمط اللغة الإنجليزية)



SIMZ-YM1-130
Fig. MD1
MEAD
upper anteriors,
roots



SIMZ-YM1-131
Fig. MD2
MEAD
upper Molars



SIMZ-YM1-132
Fig. MD3
MEAD
lower roots



SIMZ-YM1-133
Fig. MD4
MEAD
lower Molars



SIMZ-YM1-134
FELSCH
upper third molars

EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Щипці екстракційні (англійська форма)
- Daviers (Forme Anglaise)
- Forceps de Extracción (Forme Inglesa)
- Zahnzangen (Englische Form)
- استخراج العنلق (نمط اللغة الإنجليزية)



SIMZ-YM1-135
STIKSA
lower molars
and wisdoms



SIMZ-YM1-136
STIKSA
lower molars and
wisdoms roots



SIMZ-YM1-137
HERTZ
lower molars
and wisdoms roots




SIMZ-YM1-138
HERTZ
lower incisors


- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)

- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملقط (نمط اللغة الإنجليزية)


Children Pattern




SIMZ-YM1-139
Fig. 158
Children's Upper Teeth Molars Universal




SIMZ-YM1-140
Fig. 159
Children's upper Molars Universal



SIMZ-YM1-141
Fig. 160
Children's Upper Lower Universal (Hawk's Bill)




SIMZ-YM1-142
Fig. 161
Children's Lower Centrals & Roots



SIMZ-YM1-143
Fig. 162
Children's Upper Centrals



SIMZ-YM1-144
Fig. 163
Children's Upper Centrals




SIMZ-YM1-145
Fig. 37
Children's Upper Incisors & Canines


EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Щипці екстракційні (англійська форма)
- Daviers (Forme Anglaise)
- Forceps de Extracción (Forme Inglesa)
- Zahnzangen (Englische Form)
- استخراج الملقط (نمط اللغة الإنجليزية)

Children Pattern



SIMZ-YM1-146
Fig. 123
Children's Lower Teeth Universal




SIMZ-YM1-147
Fig. 138
Children's Upper Universal & Roots





SIMZ-YM1-148
Fig. 157
Children's Upper Morals Universal



SIMZ-YM1-149
Fig. 39
Children's Upper Morals Universal



SIMZ-YM1-150
Fig. 39R
Children's Upper Morals Right



SIMZ-YM1-151
Fig. 39L
Children's Upper Morals Left

WISDOM EXTRACTING FORCEPS

- Щипцы экстракционные для зубов мудрости (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні для зубів мудрості (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- الحكمة استخراج الملقط



SIMZ-YM1-152
Routerier
Fig. 22 ½ L
Lower Molars and
Wisdoms, Left



SIMZ-YM1-153
Routerier
Fig. 22 ½ R
Lower Molars and
Wisdoms, Right



SIMZ-YM1-154
Routerier
Fig. 33 ½ L
Lower Molar and
Wisdom Roots, Left



SIMZ-YM1-155
Routerier
Fig. 33 ½ R
Lower Molar and
Wisdom Roots, Right



SIMZ-YM1-156
Fig. 186L
Lower Molars



SIMZ-YM1-157
Fig. 186R
Lower Molars

EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Щипці екстракційні (англійська форма)
- Daviers (Forme Anglaise)
- Forceps de Extracción (Forme Inglesa)
- Zahnzangen (Englische Form)
- استخراج الملقط (نمط اللغة الإنجليزية)



SIMZ-YM1-158
Busch
lower roots



SIMZ-YM1-159
Busch
lower roots



SIMZ-YM1-160
Witzel
lower roots



SIMZ-YM1-161
Fig. 146
Upper and lower roots
deep-gripping



SIMZ-YM1-162
Fig. 333S
Very thin lower roots



SIMZ-YM1-163
Fig. 151
upper roots



SIMZ-YM1-164
Fig. 351S
Very thin upper roots

EXTRACTING FORCEPS (English Pattern)


- Щипцы экстракционные (английская форма)
- Щипці екстракційні (англійська форма)
- Daviers (Forme Anglaise)
- Forceps de Extracción (Forme Inglesa)
- Zahnzangen (Englische Form)
- استخراج الملقط (نمط اللغة الإنجليزية)




SIMZ-YM1-165
Fig. 121
Upper Third
Molars



SIMZ-YM1-166
Fig. 136
Upper Premolars
and Roots



SIMZ-YM1-167
Fig. 137
Lower Incisors
and Roots



SIMZ-YM1-168
Fig. 167
Universal Patter
For Lower Teeth
Molars & Premolars



SIMZ-YM1-169
Fig. 168
Upper Canines



SIMZ-YM1-170
Fig. 169
Lower Canines &
Premolars



SIMZ-YM1-171
Fig. 233
Lower Root
(Parallel Beaks)



SIMZ-YM1-172
Fig. 235
Upper Root
(Parallel Beaks)

EXTRACTING FORCEPS (For Children)

- Щипцы экстракционные детские
- Щипці екстракційні дитячі
- Daviers - forme anglaise
- Fòrceps de extracció - forma inglesa
- Pinzper per estrazione per bambini
- استخراج المنطق (تأطقل)



SIMZ-YM1-173
English Pattern
Klein Fig. 4
Pedodontic Upper root

Pediatric Tooth Extracting Forceps



SIMZ-YM1-176
English Pattern
Klein Fig. 1
Pedodontic, Upper Incisors



SIMZ-YM1-177
English Pattern
Klein Fig. 5
Pedodontic Lower Incisors.



SIMZ-YM1-174
English Pattern
Klein Fig. 2
Pedodontic Upper Premolars



SIMZ-YM1-178
English Pattern
Klein Fig. 6
Pedodontic Lower Molars



SIMZ-YM1-175
English Pattern
Klein Fig. 3
Pedodontic Upper Molars Universal

Baby Tooth Extraction Forceps

SIMZ-YM1-180
Set of 7 Pieces



SIMZ-YM1-179
English Pattern
Klein Fig. 7
Pedodontic Lower Molars,

EXTRACTING FORCEPS (For Children)

- Щипцы экстракционные детские
- Daviers - forme anglaise
- Pinzaper per estrazione per bambini

- Щипці екстракційні дитячі
- Fòrceps de extracciòn - forma inglesa
- استخراج الملقط (للأطفال)



SIMZ-YM1-181
Fig. 151 SK
Lower Incisors,
Premolars, Roots



SIMZ-YM1-182
Fig. 222 SK
Lower Third Molars



SIMZ-YM1-183
Fig. 560
Upper Incisors

Deep Gripping Root Forceps




SIMZ-YM1-184
Fig. 146
Upper Roots



SIMZ-YM1-185
Fig. 149
Lower Roots




SIMZ-YM1-186
Fig. 153
Upper Roots



1/1

Obere Wurzeln, Universal
Upper roots, Universal
Radici superiori, Universale
Raices superiores, Universal
Racines superieures, Universel

Witzel
SIMZ-YM1-187
170mm



1/1

Untere Wurzeln, Universal
Lower roots, Universal
Radici inferiori, Universale
Raices inferiores, Universal
Racines inferieures, Universel

Witzel
SIMZ-YM1-188
165mm



Witzel 221
SIMZ-YM1-189
Universal Pattern For
Upper Teeth



1/1

Obere Wurzeln, Universal
Upper roots, Universal
Radici superiori, Universale
Raices superiores, Universal
Racines superieures, Universel

SIMZ-YM1-190
Fig. 400 (300)
165mm



1/1

Untere Wurzeln, Universal
Lower roots, Universal
Radici inferiori, Universale
Raices inferiores, Universal
Racines inferieures, Universel

SIMZ-YM1-191
Fig. 401 (301)
160mm



Witzel 223
SIMZ-YM1-192
Universal Pattern For
Lower Teeth



Diamantiert
With Diamond
Con diamante
Diamantada
Diamantées

SIMZ-YM1-193
130mm



Diamantiert
With Diamond
Con diamante
Diamantada
Diamantées

SIMZ-YM1-194
130mm



SIMZ-YM1-195
Ralk
Root Fragments



SIMZ-YM1-196
130mm



SIMZ-YM1-197
130mm



SIMZ-YM1-198
125mm

ROOTS SPLINTER FORCEPS

- Щипцы для удаления корневых осколков
- Pincas Fragmentos
- Wurzelsplitterzangen



- Щипці для видалення корневих осколків
- Pinzas Para Fragmentos
- جنور منشفة العظام




SIMZ-YM1-199
Splinter Forceps
Straight
180mm



SIMZ-YM1-200
Splinter Forceps
180mm
45°



SIMZ-YM1-202
Stieglitz
Root Fragments



SIMZ-YM1-201
Stieglitz
Root Fragments



SIMZ-YM1-203
Stieglitz Splinter Forceps
140mm 90°




SIMZ-YM1-204
Stieglitz Splinter Forceps
140mm Straight

EXTRACTING FORCEPS (For Children)


- Щипцы экстракционные детские
- Daviers - forme anglaise
- Pinzeper per estrazione per bambini
- Щипці екстракційні дитячі
- Fòrceps de extracció - forma inglesa
- استخراج الملقط (للأطفال)



SIMZ-YM1-205
Koehler
Luxation Forceps
Right & Left Molars



SIMZ-YM1-206
Upper Canines
& Premolars



SIMZ-YM1-207
Lower Molars
Either Side



SIMZ-YM1-208
Lower Third Molars



SIMZ-YM1-209
Albrecht
16cm



SIMZ-YM1-210
Albrecht
16cm


deciduous teeth

EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج المنطق (نمط اللغة الإنجليزية)

Deep Gripping Root Forceps

SIMZ-YM1-219
Upper Centrals
Fig. 34A



SIMZ-YM1-220
Upper Premolar
Fig. 35A



SIMZ-YM1-221
Upper Roots
Fig. 331



SIMZ-YM1-222
Lower Premolar
Fig. 336



SIMZ-YM1-223
Lower Premolar
Fig. 36A



EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Щипці екстракційні (англійська форма)
- Daviers (Forme Anglaise)
- Forceps de Extracción (Forme Inglesa)
- Zahnzangen (Englische Form)
- استخراج الملقط (نمط اللغة الإنجليزية)

Deep Gripping Root Forceps



SIMZ -YM1-224
Upper Third Molars
Fig. 67N



SIMZ -YM1-225
Lower Third Molars
Fig. 79N



SIMZ -YM1-226
Upper Premolars
Fig. 7N



SIMZ -YM1-227
Lower Premolars
Fig. 36N


EXTRACTING FORCEPS (English Pattern Hole Handle)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracció (Forme Inglesa)
- استخراج المنقط (نمط اللغة الإنجليزية)

(Deep Gripe Extracting Forceps)



- No Fracture
- No Tissue Damage
- Extraction With Easy



Thin Tip



SIMZ -YM1-228
Fig. 34
Upper Anterior

Wide Tip



SIMZ -YM1-229
Fig. 34A
Upper Anterior

Thin Tip



SIMZ -YM1-230
Fig. 35
Upper Universal

Wide Tip



SIMZ -YM1-231
Fig. 35A
Upper Universal

Thin Tip




SIMZ -YM1-232
Fig. 321
Lower Universal

A-TRAUMATIC EXTRACTION FORCEPS
Improved tapered beaks for complete root contact.
Cross serrations for a never slip grasp.
Easily enter the sub gingival.

EXTRACTING FORCEPS (English Pattern Hole Handle)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)

- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملقط (نمط اللغة الإنجليزية)




Deep Root Pick Forceps



EXTRACTING FORCEPS (English Pattern Hole Handle)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج المنقط (نقط اللغة الإنجليزية)



DEEP ROOT PICK FORCEPS



SIMZ-YM1-238
Fig. 79
Lower Third Molars



SIMZ-YM1-239
Fig. 79A
Apex Forceps Lower Universal
Deep Root Pick



SIMZ-YM1-240
Fig. 79 T
Lower Third Molars



SIMZ-YM1-241
Fig. 44
Apex Forceps Upper Root
Deep Root Pick




SIMZ-YM1-242
Fig. 45
Apex Forceps Lower Root
Deep Root Pick

EXTRACTING FORCEPS (English Pattern Hole Handle)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)


- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج المنطق (نمط الثقفة الإنجليزية)



DEEP ROOT PICK FORCEPS



- Щипцы экстракционные (английская форма)
- Щипці екстракційні (англійська форма)
- Daviers (Forme Anglaise)
- Forceps de Extracción (Forme Inglesa)
- Zahnzangen (Englische Form)
- استخراج الملقط (نمط الفتحة الإنجليزية)



DEEP ROOT PICK FORCEPS



SIMZ-YM1-243
Fig. 51
 Apex Forceps
 Upper Root
 Deep Root Pick



SIMZ-YM1-244
Fig. 46XL
 Apex Forceps Lower Root
 Deep Root Pick



SIMZ-YM1-245
 Upper Root
Fig. 44



SIMZ-YM1-246
 Lower Root
Fig. 45



SIMZ-YM1-247
 Upper Root
Fig. 51

EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Daviers - forme américaine
- Forceps de extracción - forma americana
- Щипці екстракційні (американська форма)
- Pinzere per estrazione - forma americana
- مستفراج الملقط (النمط الأمريكي)

Carbide Tipped Root Tip Forceps



Anatomical Handle



SIMZ -YM1-248
T.C.
Carbide Serrated Beak
Lower Root Tips & Molars



SIMZ-YM1-249
Root Tip-Lower
Carbide Serrated Beak
Lower Root Tips & Molars



SIMZ-YM1-250
Root Tip-Upper
Carbide Serrated Beak
Lower Root Tips & Molars

German Stainless Steel Extraction Forceps



SIMZ-YM1-251
Fig. 150AS
Serrated Beak
Upper Anteriors,
Split Beaks



SIMZ-YM1-252
Fig. 151AS
Serrated Beak
Lower Anteriors,
Split Beaks



SIMZ-YM1-253
Fig. 23
1st & 2nd Lower Molars
Universal, Cow Horn

Wrap Around Forcep



SIMZ-YM1-254
Wrap around extraction forcep
right angle "Wrap Around" Design
Provides Superior Access
Serrated Beak

EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма) • Щипці екстракційні (американська форма)
- Daviers - forme américaine • Pinzere per estrazione - forma americana
- Férceps de extracción - forma americana • استخراج الملقط (المنظ الأمريكي)

Diamond Coated Extraction forceps - Better Grip & Tactile Feel - Less Trauma



SIMZ-YM1-256
Fig. 151
Carbide Serrated Beak
Lower Incisors, Cuspids,
Bicuspid & Roots, Universal



SIMZ-YM1-257
Fig. 1STD
Carbide Serrated Beak
Upper Incisors, Cuspids,



SIMZ-YM1-258
Fig. 17
Carbide Serrated Beak
1st & 2nd Lower Molars,
Universal



SIMZ-YM1-259
Fig. 53L
Carbide Serrated Beak
1st & 2nd Upper Molars,
Left



SIMZ-YM1-260
Fig. 53R
Carbide Serrated Beak
1st & 2nd Upper Molars,
Right



SIMZ-YM1-261
Fig. 65
Carbide Serrated Beak
Upper Incisors, Roots &
Fragments



SIMZ-YM1-255
150
Carbide Serrated Beak
Upper Incisors, Cuspids,
Bicuspid & Roots,
Universal



SIMZ-YM1-262
Fig. 69
Carbide Serrated Beak
Upper & Lower Fragments
And Roots



SIMZ-YM1-263
Fig. 150S Pedro forcep
Carbide Serrated Beak
Pedo, Upper Primary Teeth
and Roots, Universal



SIMZ-YM1-264
Fig. 151S Pedro forcep
Carbide Serrated Beak
Pedo, Lower Primary Teeth
and Roots, Universal

EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Щипці екстракційні (американська форма)
- Daviers - forme américaine
- Pinzere per estrazione - forma americana
- Fórceps de extracción - forma americana
- استخراج المنطق (التمط الأمريكي)



SIMZ-YM1-265
45° Angled Root Forcep
"Fish-Hook" Serrations
Upper Atraumatic Extraction



SIMZ-YM1-266
90° Angled Root Forcep
"Fish-Hook" Serrations
Upper Atraumatic Extraction



SIMZ-YM1-267
Fig. 01
witzel
Radici inferiori
Lower roots



SIMZ-YM1-268
Fig. 02
witzel
Radici superiori
Upper roots



SIMZ-YM1-269
Fig. 01
Straight for pins
and fragments



SIMZ-YM1-270
Fig. 00 - Lindo Levien
Straight for roots and
fragments - TC coating



SIMZ-YM1-271
Fig. 01 - Lindo Levien
Curved for roots and
fragments - TC coating




SIMZ-YM1-272
Fig. 02
Curved for pins
and fragments

EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Щипці екстракційні (американська форма)
- Pinzper per estrazione - forma americana
- Stخراج المنطق (النمط الأمريكي)
- Daviers - forme américaine
- Fôrceps de extracción - forma americana



SIMZ-YM1-273
Fig. W1
Winter
upper centrals
and canines



SIMZ-YM1-274
Fig. 1
upper incisors
and roots



SIMZ-YM1-275
Fig. 2
upper and
lower roots



SIMZ-YM1-276
Fig. 4
upper and
lower premolars



SIMZ-YM1-277
Fig. 5
Physick
lower third molars



SIMZ-YM1-278
Fig. 6
separating forceps
spaltzange




EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Щипці екстракційні (американська форма)
- Daviers - forme américaine
- Pinzere per estrazione - forma americana
- Fórceps de extracción - forma americana
- استخراج الملقط (النمط الأمريكي)



SIMZ-YM1-279
Fig. 7
upper roots



SIMZ-YM1-280
Fig. 8
upper and lower
third molars



SIMZ-YM1-282
Fig. 10 H
upper molars



SIMZ-YM1-283
Fig. 13
upper incisors



SIMZ-YM1-281
Fig. 10S
upper roots




SIMZ-YM1-284
Fig. 14
lower incisors

EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Щипці екстракційні (американська форма)
- Pinzper per estrazione - forma americana
- Fórceps de extracción - forma americana
- استخراج المنطق (النمط الأمريكي)



SIMZ -YM1-285
Fig. 15
Harris
lower molars



SIMZ -YM1-287
Fig. 16
lower molars



SIMZ -YM1-288
Fig. 17
lower molars



SIMZ -YM1-286
Fig. 16S
lower molars
for children




SIMZ -YM1-289
Fig. 18 L
Harris
upper molars,
left



SIMZ -YM1-290
Fig. 18R
Harris
upper molars,
right

EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Щипці екстракційні (американська форма)
- Daviers - forme américaine
- Pinzeta per estrazione - forma americana
- Fórceps de extracción - forma americana
- استخراج المنطق (النمط الأمريكي)



SIMZ-YM1-291
Fig. 20R
upper molars,
Right

SIMZ-YM1-292
Fig. 20L
upper molars,
left

SIMZ-YM1-293
Fig. 23
lower molars

SIMZ-YM1-294
Fig. 24
upper molars

SIMZ-YM1-295
Fig. 27
lower molars

SIMZ-YM1-296
Fig. 29S
lower teeth,
for children

EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Щипці екстракційні (американська форма)
- Daviers - forme américaine
- Pinzere per estrazione - forma americana
- Fórceps de extracción - forma americana
- استخراج المقلد (الأمريكي)




SIMZ -YM1-297
Fig. 30S
upper teeth,
for children



SIMZ -YM1-298
Fig. 32
upper premolars
and molars




SIMZ -YM1-299
Fig. 32A
upper premolars
and roots



SIMZ -YM1-300
Fig. 30S
lower anteriors,
and roots




SIMZ -YM1-301
Fig. 35
upper roots



SIMZ -YM1-302
Fig. 36
lower anteriors,
and roots

EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Щипці екстракційні (американська форма)
- Daviers - forme américaine
- Pinzere per estrazione - forma americana
- Fórceps de extracción - forma americana
- استخراج العنقظ (النمط الأمريكي)



SIMZ -YM1-303
Fig. 37
lower anteriors
and roots




SIMZ -YM1-305
Fig. 53 R
upper molars,
right



SIMZ -YM1-306
Fig. 53 L
upper molars,
left



SIMZ -YM1-304
Fig. 38
upper anteriors
and roots



SIMZ -YM1-307
Fig. 62
incisors, premolars,
SIMZ deciduous teeth



SIMZ -YM1-308
Fig. 65
upper incisors
and roots

EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Щипці екстракційні (американська форма)
- Daviers - forme américaine
- Pinzere per estrazione - forma americana
- Férceps de extracción - forma americana
- استخراج الملقط (النمط الأمريكي)



SIMZ-YM1-309
Fig. 69 Tomas
upper and
lower roots



SIMZ-YM1-310
Fig. 85 Tomas
lower canines,
premolars and molars



SIMZ-YM1-311
Fig. 88R Nevius
upper molars,
right



SIMZ-YM1-312
Fig. 85A Tomas
lower canines,
premolars and molars



SIMZ-YM1-313
Fig. 88L Nevius
upper molars,
left

EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма) • Щипці екстракційні (американська форма)
- Daviers - forme américaine • Pinzere per estrazione - forma americana
- Férceps de extracción - forma americana • استخراج الملقط (النمط الأمريكي)



SIMZ-YM1-314
Fig. 99A
Kells
 upper incisors,
 canines, premolars



SIMZ-YM1-315
Fig. 104
 upper third molars



SIMZ-YM1-316
Fig. 99C
Kells
 upper incisors,
 canines, premolars



SIMZ-YM1-317
Fig. 101
Hull
 Deciduous teeth




SIMZ-YM1-318
Fig. 103
 lower incisors,
 premolars

EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма) • Щипці екстракційні (американська форма)
- Daviers - forme américaine • Pinziper per estrazione - forma americana
- Fôrceps de extracción - forma americana • استخراج المنطق (النمط الأمريكي)



SIMZ-YM1-319
Fig. 150
Cryer
upper incisors,
premolars,
roots



SIMZ-YM1-320
Fig. 150A
Cryer
upper incisors,
premolars,
roots



SIMZ-YM1-321
Fig. 150S
Cryer
upper teeth,
for children



SIMZ-YM1-322
Fig. 150X
Upper premolars,
incisors and roots




SIMZ-YM1-326
Fig. 151
Cryer
lower teeth,
for children



SIMZ-YM1-323
Fig. 150XAS
Upper premolars,
incisors and
roots -
also for children



SIMZ-YM1-324
Fig. 151A
Cryer
lower incisors,
premolars,
roots



SIMZ-YM1-325
Fig. 151 S
Cryer
lower incisors,
premolars,
roots

EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Щипці екстракційні (американська форма)
- Daviers - forme américaine
- Pinzere per estrazione - forma americana
- Férceps de extracción - forma americana
- استخراج المنطق (النمط الأمريكي)



SIMZ -YM1-327
Fig. 203
 lower incisors
 and roots



SIMZ -YM1-328
Fig. 210S
 upper third
 molars



SIMZ -YM1-330
Fig. 210H
 upper third
 molars



SIMZ -YM1-331
Fig. 213
 upper centrals
 and canines




SIMZ -YM1-329
Fig. 201
 lower incisors,
 canines,
 premolars



SIMZ -YM1-332
Fig. 217
 lower molars

EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Щипці екстракційні (американська форма)
- Pinzere per estrazione - forma americana
- Férceps de extracción - forma americana
- سكراج الملقط (النمط الأمريكي)




SIMZ -YM1-333
Fig. 222
 lower third molars



SIMZ -YM1-335
Fig. 287 Swing
 lower molars



SIMZ -YM1-336
Fig. 288
 Lower canines
 and premolars




SIMZ -YM1-334
Fig. 286
 upper incisors
 and roots

EXTRACTING FORCEPS (American Pattern)


- Щипцы экстракционные (американская форма)
- Щипці екстракційні (американська форма)
- Pinzere per estrazione - forma americana
- Férceps de extracción - forma americana
- استخراج الملقط (النمط الأمريكي)



SIMZ-YM1-337
Cohen
universal forceps



SIMZ-YM1-338
Fig. 10 S Nevius
upper canines



SIMZ-YM1-339
Fig. 1 A Henahan
upper incisors



SIMZ-YM1-340
Fig. 10 H Nevius
upper canines



SIMZ-YM1-341
Fig. 4 Standard
upper third molars

EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Щипці екстракційні (американська форма)
- Daviers - forme américaine
- Pinzere per estrazione - forma americana
- Fórceps de extracción - forma americana
- استخراج المنقط (النمط الأمريكي)



SIMZ-YM1-342
Fig. FH Woodward
Lower molars



SIMZ-YM1-343
Fig. FS Woodward
Lower molars



SIMZ-YM1-344
Fig. 89
Upper molars,
right



SIMZ-YM1-345
Fig. 90
Upper molars,
left

EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Щипці екстракційні (американська форма)
- Daviers - forme americane
- Pinzoper per estrazione - forma americana
- Fórceps de extracción - forma americana
- استخراج المنطق (النمط الأمريكي)



SIMZ-YM1-346
Fig. 23 S
Lower molars




SIMZ-YM1-347
Fig. 17 SK
Lower molars



SIMZ-YM1-348
Fig. 150 SK
Upper anteriors
and roots



SIMZ-YM1-349
Fig. 151 SK
lower anteriors
and roots



SIMZ-YM1-350
Fig. 23 SK
Lower molars

EXTRACTING FORCEPS (Anatomical Handle)


- Щипцы экстракционные (с эргономичной формой ручки)
- Давиеры (avec manche d'ajustage)
- Zahnzangen (mit Paßformgriff)
- Щипці екстракційні (з ергономічною формою ручки)
- Forceps de Extracción (con mango adaptado a la mano)
- استخراج الملقط (مع تركيب مقابض)




SIMZ -YM1-351
Fig. 1
upper laterals
and canines



SIMZ -YM1-352
Fig. 2
upper laterals
and canines



SIMZ -YM1-353
Fig. 3
upper incisors
and canines



SIMZ -YM1-354
Fig. 4
lower incisors
and canines



SIMZ -YM1-355
Fig. 7
upper premolars



SIMZ -YM1-356
Fig. 13
lower premolars



SIMZ -YM1-357
Fig. 13 S
lower premolars,
for children

EXTRACTING FORCEPS (Anatomical Handle)

- Щипцы экстракционные (с эргономичной формой ручки)
- Daviers (avec manche d'ajustage)
- Zahnzangen (mit Paßformgriff)

- Щипці екстракційні (з ергономічною формою ручки)
- Forceps de Extracción (con mango adaptado a la mano)
- استخراج المائلط (مع تركيب مقابض)



SIMZ-YM1-358
Fig. 17
upper molars,
right



SIMZ-YM1-359
Fig. 18
upper molars,
left



SIMZ-YM1-360
Fig. 18 A
upper molars



SIMZ-YM1-361
Fig. 18 B
upper molars



SIMZ-YM1-362
Fig. 19
upper third molars



SIMZ-YM1-363
lower molars



SIMZ-YM1-364
Fig. 22 S
lower molars,
for children

EXTRACTING FORCEPS (Anatomical Handle)

- Щипцы экстракционные (с эргономичной формой ручки)
- Щипці екстракційні (з ергономічною формою ручки)
- Forceps de Extracción (con mango adaptado a la mano)
- استفرآج الملظط (مع تركيب مقابض)
- Daviers (avec manche d'ajustage)
- Zahnzangen (mit Paßformgriff)




SIMZ-YM1-365
Fig. 29
upper roots



SIMZ-YM1-366
Fig. 29 S
upper roots,
for children



SIMZ-YM1-367
Fig. 30
upper roots



SIMZ-YM1-368
Fig. 30 S
upper roots,
for children



SIMZ-YM1-369
Fig. 33
lower roots



SIMZ-YM1-370
Fig. 33 A
lower roots



SIMZ-YM1-371
Fig. 33 B
lower roots

EXTRACTING FORCEPS (Anatomical Handle)

- Щипцы экстракционные (с эргономичной формой ручки)
- Щипці екстракційні (з ергономічною формою ручки)
- Daviers (avec manche d'ajustage)
- Forceps de Extracción (con mango adaptado a la mano)
- Zahnzangen (mit Paßformgriff)
- استخراج الملقط (مع تركيب مقابض)



SIMZ -YM1-372
Fig. 33 L
lower roots



SIMZ -YM1-373
Fig. 33 S
lower roots,
for children



SIMZ -YM1-374
Fig. 37
upper incisors,
for children



SIMZ -YM1-375
Fig. 39
upper molars,
for children



SIMZ -YM1-376
Fig. 39 A
upper molars,
for children



SIMZ -YM1-377
Fig. 39 R
upper molars right,
for children



SIMZ -YM1-378
Fig. 39 L
upper molars left,
for children

EXTRACTING FORCEPS (Anatomical Handle)

- Щипцы экстракционные (с эргономичной формой ручки)
- Щипці екстракційні (з ергономічною формою ручки)
- Daviers (avec manche d'ajustage)
- Forceps de Extracción (con mango adaptado a la mano)
- Zahnzangen (mit Paßformgriff)
- استخراج الملقط (مع تركيب مقابض)



SIMZ-YM1-379
Fig. 41
upper roots



SIMZ-YM1-380
Fig. 44
upper roots



SIMZ-YM1-381
Fig. 45
lower roots



SIMZ-YM1-382
Fig. 48
lower roots



SIMZ-YM1-383
Fig. 51
upper roots



SIMZ-YM1-384
Fig. 51 A
upper roots



SIMZ-YM1-385
Fig. 51 S
upper roots,

EXTRACTING FORCEPS (Anatomical Handle)

- Щипцы экстракционные (с эргономичной формой ручки)
- Daviers (avec manche d'ajustage)
- Zahnzangen (mit Paßformgriff)

- Щипці екстракційні (з ергономічною формою ручки)
- Forceps de Extracción (con mango adaptado a la mano)
- استخراج الملقط (مع تركيب مطابق)



SIMZ -YM1-386
Fig. 51 L
upper roots



SIMZ -YM1-387
Fig. 52
upper roots



SIMZ -YM1-388
Fig. 54
for separating
upper roots



SIMZ -YM1-389
Fig. 65 R
upper molars,
right



SIMZ -YM1-390
Fig. 65 L
upper molars, left



SIMZ -YM1-391
Fig. 67
upper third molars



SIMZ -YM1-392
Fig. 67 A
upper third molars

EXTRACTING FORCEPS (Anatomical Handle)

- Щипцы экстракционные (с эргономичной формой ручки)
- Щипці екстракційні (з ергономічною формою ручки)
- Daviers (avec manche d'ajustage)
- Forcets de Extracción (con mango adaptado a la mano)
- Zahnzangen (mit Paßformgriff)
- استخراج المنقط (مع تركيب مقابض)



SIMZ-YM1-393
Fig. 67 1/2 R
upper third molars,
right



SIMZ-YM1-394
Fig. 67 1/2 L
lower third molars,
left.



SIMZ-YM1-395
Fig. 73
lower molars



SIMZ-YM1-396
Fig. 74
lower roots



SIMZ-YM1-397
Fig. 74 M
lower roots



SIMZ-YM1-398
Fig. 74 N
lower roots



SIMZ-YM1-399
Fig. 74 XN
lower roots and crowded
incisors

EXTRACTING FORCEPS (Anatomical Handle)

- Щипцы экстракционные (с эргономичной формой ручки)
- Щипці екстракційні (з ергономічною формою ручки)
- Forceps de Extracción (con mango adaptado a la mano)
- استخراج المفلط (مع تركيب مقابض)
- Daviers (avec manche d'ajustage)
- Zahnzangen (mit Paßformgriff)



SIMZ-YM1-400
Fig. 75
lower premolars



SIMZ-YM1-401
Fig. 76
upper roots



SIMZ-YM1-402
Fig. 76N
upper roots



SIMZ-YM1-403
Fig. 81
lower premolars
SIMZ and roots



SIMZ-YM1-404
Fig. 79A
lower third molars



SIMZ-YM1-405
Fig. 79
lower third molars

EXTRACTING FORCEPS (Anatomical Handle)

- Щипцы экстракционные (с эргономичной формой ручки)
- Щипці екстракційні (з ергономічною формою ручки)
- Daviers (avec manche d'ajustage)
- Forceps de Extracción (con mango adaptado a la mano)
- Zahnzangen (mit Paßformgriff)
- استخراج المثلث (مع تركيب مقابض)



SIMZ-YM1-406
Fig. 86 A
lower molars



SIMZ-YM1-407
Fig. 86 B
lower molars



SIMZ-YM1-408
Fig. 86 C
lower molars



SIMZ-YM1-409
Fig. 99
lower molars



SIMZ-YM1-410
Fig. 99 1/4
lower molars



SIMZ-YM1-411
Fig. 136
upper premolars
and roots



SIMZ-YM1-412
Fig. 143
lower molars

EXTRACTING FORCEPS (American Pattern Hole Handle)

- Щипцы экстракционные (американская форма)
- Щипці екстракційні (американська форма)
- Daviers - forme américaine
- Pinzaper per estrazione - forma americana
- Fórceps de extracción - forma americana
- استخراج الملقط (النمط الأمريكي)



SIMZ-YM1-413
Fig. 1
Upper Anterior



SIMZ-YM1-414
Lower Molar
Fig. 23 Cow Horn



SIMZ-YM1-415
Fig. 88L
Upper Molar Left



SIMZ-YM1-416
Fig. 88R
Upper Molar Right



SIMZ-YM1-417
Fig. 150
Upper Universal

EXTRACTING FORCEPS (American Pattern Hole Handle)

- Щипцы экстракционные (американская форма)
- Щипці екстракційні (американська форма)
- Daviers - forme américaine
- Pinzere per estrazione - forma americana
- Fórceps de extracción - forma americana
- استخراج الملقط (النمط الأمريكي)



SIMZ-YM1-418
Fig. 150S
Upper Universal



SIMZ-YM1-419
Fig. 151
Lower Universal



SIMZ-YM1-420
Fig. 151S
Lower Universal



SIMZ-YM1-421
Fig. 106
Lower Incisors



SIMZ-YM1-422
Fig. 174
Extra Narrow Lower Anterior
Apical Retention Forceps

EXTRACTING FORCEPS (Anatomical Handle)

- Щипцы экстракционные (с эргономичной формой ручки)
- Щипці екстракційні (з ергономічною формою ручки)
- Forceps de Extracción (con mango adaptado a la mano)
- استخراچ المثلث (مع تركيب مقابض)
- Daviers (avec manche d'ajustage)
- Forceps de Extracción (con mango adaptado a la mano)



SIMZ-YM1-358
Fig. 17
upper molars,
right



SIMZ-YM1-359
Fig. 18
upper molars,
left



SIMZ-YM1-360
Fig. 18 A
upper molars




SIMZ-YM1-361
Fig. 18 B
upper molars



SIMZ-YM1-362
Fig. 19
upper third molars



SIMZ-YM1-363
lower molars



SIMZ-YM1-364
Fig. 22 S
lower molars,
for children

EXTRACTING FORCEPS (Anatomical Handle)

- Щипцы экстракционные (с эргономичной формой ручки)
- Щипці екстракційні (з ергономічною формою ручки)
- Daviers (avec manche d'ajustage)
- Forceps de Extracción (con mango adaptado a la mano)
- Zahnzangen (mit Paßformgriff)
- استخراج المثلث (مع تركيب مقابض)




SIMZ-YM1-365
Fig. 29
upper roots



SIMZ-YM1-366
Fig. 29 S
upper roots,
for children



SIMZ-YM1-367
Fig. 30
upper roots



SIMZ-YM1-368
Fig. 30 S
upper roots,
for children



SIMZ-YM1-369
Fig. 33
lower roots



SIMZ-YM1-370
Fig. 33 A
lower roots



SIMZ-YM1-371
Fig. 33 B
lower roots

FUNCTION OF ROOT ELEVATORS IN SURGERY

To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.


HANDLES FOR ROOT ELEVATORS POINTS A, B, C, D, E, F



ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elevadores para raices
- Wurzelheber
- Elevator corneveiy ekstraktsionnyy
- المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.



- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Elevator corneviy ekstraktsionniy
- Elevadores para raices
- المعاصر الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Atraumatic Extraction Elevators & Socket Scoops
Socket Scoops For Cleaning Out Extraction Sites



SIMZ-YM2-018
 Fig. 5B SA
 Schmeckebeier
 Apexo
 Concave tip



SIMZ-YM2-019
 Fig. 4 SA
 Schmeckebeier
 Apexo
 Flat tip



SIMZ-YM2-020
 Fig. 5 SA
 Schmeckebeier
 Apexo
 Flat tip



SIMZ-YM2-021
 Misch Socket Scoop
 4mm Angled Left



SIMZ-YM2-022
 Misch Socket Scoop
 4mm Angled Right

SIMZ-YM2-017
 Fig. 4B SA
 Schmeckebeier Apexo
 Concave tip



SIMZ-YM2-023
 Potts Elevator-
 Left



SIMZ-YM2-024
 Potts Elevator-
 Right



SIMZ-YM2-025
 Potts "T-Handle"
 Elevator-Left



SIMZ-YM2-026
 Potts "T-Handle"
 Elevator-Right

ROOT ELEVATOR

- Элеватор корневой экстракционный
 - Elévateurs à racines
 - Wurzelheber
- Elevator corneviy ekstraktsionniy
 - Elevadores para raices
 - المعاصر الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Apexo Root Elevators



Apexo Root Elevator



FUNCTION: To loosen tooth or root from bony socket prior to placement of the extraction forceps.



Luxator Root Elevators



Luxator Root Elevator




- Элеватор корневой экстракционный
- Elevateurs à racines
- Wurzelheber

- Elevator cornevier ekstraktsionnyy
- Elevadores para raices
- المصاعد الجذر


FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Bein Root Elevators




Bein Root Elevator

3.0mm	4.0mm	3.0mm	4.0mm	3.2mm	3.8mm
SIMZ-YM2-041 Bein Fig. 1B	SIMZ-YM2-042 Bein Fig. 2B	SIMZ-YM2-043 Bein Fig. 3B	SIMZ-YM2-044 Bein Fig. 4B	SIMZ-YM2-045 Bein Fig. 5B	SIMZ-YM2-046 Bein Fig. 6B



Seldin Root Elevators



Seldin Root Elevator


4 mm	5 mm	10 mm	10 mm	13 mm	13 mm	3.5mm
SIMZ-YM2-047 Seldin 34S	SIMZ-YM2-048 Seldin 34	SIMZ-YM2-049 Seldin 4L (Left)	SIMZ-YM2-050 Seldin 4R (Right)	SIMZ-YM2-051 Seldin 1L (Left)	SIMZ-YM2-052 Seldin 1R (Right)	SIMZ-YM2-053 Seldin

ROOT ELEVATOR


- Elevator корневой экстракционный
- Elevator корневий экстракційний
- Elévateurs à racines
- Elevadores para raíces
- Wurzelheber
- المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.


Flohr Root Elevators




SIMZ-YM2-054
Fig. 1F
Str.
Flohr Root Elevators




SIMZ-YM2-055
Fig. 2F
Left
Flohr




SIMZ-YM2-056
Fig. 3F
Right
Flohr




SIMZ-YM2-059
Fig. 4 SA
Schmeckebier
Apexo



SIMZ-YM2-060
Fig. 5 SA
Schmeckebier
Apexo




SIMZ-YM2-057
Fig. 4B SA
Schmeckebier
Apexo




SIMZ-YM2-058
Fig. 5B SA
Schmeckebier
Apexo


Coupland Root Elevators



SIMZ-YM2-061
Coupland No.1
Straight




SIMZ-YM2-062
Coupland No.2




SIMZ-YM2-063
Coupland No.3

Apical Root Elevators




SIMZ-YM2-064
Apical 301


Apical Root Elevator



SIMZ-YM2-065
Apical 302
(Left)




SIMZ-YM2-066
Apical 303
(Right)




SIMZ-YM2-067
Apical 301


Apical Root Elevator




SIMZ-YM2-068
Apical 302




SIMZ-YM2-069
Apical 303



SIMZ-YM2-070
Apical 302A



SIMZ-YM2-071
Apical 303A



SIMZ-YM2-072
Fig. 304

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Elevator corneveiy ekstraktsionnyy
- Elevadores para raices
- المعصاع الجذر

Laxutor Root Elevators

Long Type



Laxutor Root Elevator



SIMZ-YM2-073
Laxutor
1.5mm
Str.



SIMZ-YM2-074
Laxutor
2mm
Str.



SIMZ-YM2-075
Laxutor
3mm
Str.



SIMZ-YM2-076
Laxutor
5mm
Str.



SIMZ-YM2-077
Laxutor
3mm
Cvd.

FUNCTION:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



SIMZ-YM2-078
Laxutor
5mm
Cvd.



SIMZ-YM2-079
Laxutor
3mm
Upper
Angle



SIMZ-YM2-080
Laxutor
3mm
Double
Angle



SIMZ-YM2-081
Laxutor
3mm
Str.




SIMZ-YM2-082
Laxutor
5mm
Str.

ROOT ELEVATOR


- Элеватор корневой экстракционный
- Elevateur à racines
- Wurzelheber
- Elevator corneviu extracționii
- Elevadores para raices
- المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Cryer Root Elevators




Miller Root Elevators



- Элеватор корневой экстракционный
 - Elévateurs à racines
 - Wurzelheber
- Elevator corneviu extracționiu
 - Elevadores para raices
 - المعصاع الجذر


FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Lindo Levian Root Elevators



Lindo Levian Root Elevator

Friedmann Roots Elevator



Friedmann Roots Elevator

Heidbrink Root Elevators



Heidbrink Root Elevator


ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Elevator corneviy ekstraktsionnyy
- Elevadores para raices
- المصاعد الجذر

Function of Root Elevators: To loosen tooth or root from bony socket prior to placement of the extraction forceps.




SOS-YM2-107
Fig. 1
Apical 301



SIMZ-YM2-108 Fig. 2 Apical 302
SIMZ-YM2-109 Fig. 3 Apical 303
SIMZ-YM2-110 Fig. 4 Apical 304
SIMZ-YM2-111 Fig. 5 12 A
SIMZ-YM2-112 Fig. 6 34 S Seldin
SIMZ-YM2-113 Fig. 7 34 Seldin
SIMZ-YM2-114 Fig. 8

Function of Root Elevators: To loosen tooth or root from bony socket prior to placement of the extraction forceps.



SIMZ-YM2-115 Fig. 9 Flohr
SIMZ-YM2-116 Fig. 10 Flohr
SIMZ-YM2-117 Fig. 11 Flohr
SIMZ-YM2-118 Fig. 12 3 Seldin
SIMZ-YM2-119 710/14 3 A Cryer
SIMZ-YM2-120 710/13
SIMZ-YM2-121 Fig. 15 2 Seldin

- Элеватор корневой экстракционный
- Elevator cornevier ekstraktsionnyy
- Elévateurs à racines
- Elevadores para raíces
- Wurzelheber
- المعصاع الجذر

Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



SIMZ-YM2-122
Fig. 301
Apex
SIMZ

For displacement of deeply seated roots



SIMZ-YM2-123
Fig. 303
Apexo



SIMZ-YM2-124
Fig. 301W
Seldin



SIMZ-YM2-125
Fig. 304W
Seldin



SIMZ-YM2-126
Fig. 1
Heidbrink
Straight
Upper and lower anterior teeth and roots



SIMZ-YM2-127
Fig. 2
Heidbrink
Straight
Upper and lower anterior teeth and roots

HEIDBRINK



upper root elevator



upper root elevator



SIMZ-YM2-128
Fig. 80



SIMZ-YM2-129
Fig. 81



SIMZ-YM2-130
Fig. 81A



SIMZ-YM2-131
Fig. 81E



SIMZ-YM2-132
Fig. 40
Heidbrink
Straight



SIMZ-YM2-133
Fig. 41
Heidbrink
Straight

Upper anterior roots and small teeth For displacement of anterior roots

Upper and lower anterior teeth and roots

ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Elevator corneviu extracționii
- Elevadores para raíces
- المصاعد الجذر



SIMZ-YM2-134	SIMZ-YM2-135	SIMZ-YM2-136	SIMZ-YM2-137	SIMZ-YM2-138	SIMZ-YM2-139	SIMZ-YM2-140	SIMZ-YM2-141	SIMZ-YM2-142	SIMZ-YM2-143
Luxator 4mm straight	Luxator 3.5mm straight	Apexo 3mm Straight	Apexo 2.5mm Straight	Apexo 2mm Straight	Luxator 4mm curved	Luxator 3.5mm curved	Luxator 3mm curved	Luxator 2.5mm curved	Luxator 2mm curved

Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



SIMZ -YM2-144
Straight
3mm



SIMZ-YM2-145	SIMZ-YM2-146	SIMZ-YM2-147	SIMZ-YM2-148	SIMZ-YM2-149	SIMZ-YM2-150	SIMZ-YM2-151
Curved 3mm	Straight 5mm	Curved 5mm	Straight 2mm	Curved 2mm	Straight 4mm	Curved 4mm

LUXATING ELEVATORS



Upper root elevator with helicoidal tip



Upper root elevator with helicoidal tip

- Элеватор корневой экстракционный
- Elevator cornevei ekstraktsionnyi
- Elévateurs à racines
- Elevadores para raíces
- Wurzelheber
- المصاعد الجذر

Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



SIMZ-YM2-152
Fig. 190

SIMZ-YM2-153
Fig. 191



SIMZ-YM2-154
Fig. 27
27 Cryer



SIMZ-YM2-155
Fig. 28
28 Cryer



SIMZ-YM2-156
Fig. 27 S
27 S Cryer



SIMZ-YM2-157
Fig. 28 S
28 S Cryer



SIMZ-YM2-158
710/29
(344)



SIMZ-YM2-159
710/30
(345)



SIMZ-YM2-160
Fig. 1
Hylin 4,5mm



SIMZ-YM2-161
Fig. 2
Hylin
3mm



SIMZ-YM2-162
Fig. 3
Hylin
3mm



SIMZ-YM2-163
Fig. 4
Hylin
3mm



SIMZ-YM2-164
Hylin 153mm



SIMZ-YM2-165
Read



SIMZ-YM2-166
Read



SIMZ-YM2-168
Coleman



SIMZ-YM2-167
London Hospital



SIMZ-YM2-169



Hylin lower roots elevator

ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Elevator corneviy ekstraktsionniy
- Elevadores para raices
- المصاعد الجذر

Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



Cryer root or fragment elevator

CRYER



ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elevateurs à racines
- Wurzelheber

- Elevator corneveiy ekstraktsionnyy
- Elevadores para raices
- المصاعد الجذر



Function of Root Elevators:
To loosen tooth or root from bony socket prior to placement of the extraction forceps.



Seldin root or fragment elevator

SELDIN



SIMZ-YM2-186
Fig. 34S
Thinner tip and more narrow than the 34 shank



SIMZ-YM2-187 Fig. 1L
SIMZ-YM2-188 Fig. 1R
SIMZ-YM2-189 Fig. 4L
SIMZ-YM2-190 Fig. 4R

Upper and lower molar roots



SIMZ-YM2-191 Fig. 51 Curtis
SIMZ-YM2-192 Fig. 52 Curtis
SIMZ-YM2-193 Fig. 11 Stout 3rd Molars
SIMZ-YM2-194 Fig. 11A Stout 3rd Molars



SIMZ-YM2-195 Fig. 34 Upper and lower anterior roots
SIMZ-YM2-196 Fig. 2S Lower 3rd Molars
SIMZ-YM2-197 Fig. 3S Upper and lower anterior roots



SIMZ-YM2-198 Fig. 8 Crane Pick
SIMZ-YM2-199 Fig. 44C Bradley For displacement of roots and 3rd Molars
SIMZ-YM2-200 A Cogswell
SIMZ-YM2-201 B Cogswell

ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elevateurs à racines
- Wurzelheber
- Elevator corneviy ekstraktsionnyy
- Elevadores para raices
- المصاعد الجذر

Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



SIMZ-YM2-202
Fig. 4
Schmeckebeier
Apexo

Schmeckebeier and Miller designs used for displacement of broken or deep-seated root tips.



SIMZ-YM2-203
Fig. 5
Schmeckebeier
Apexo



SIMZ-YM2-204
Fig. 71
Miller
Apexo



SIMZ-YM2-205
Fig. 72
Miller
Apexo



SIMZ-YM2-206
Fig. 6
Potts



SIMZ-YM2-207
Fig. 7
Potts

For impacted upper 3rd Molars without bony obstruction.



SIMZ-YM2-208
Fig. 73
Miller
Apexo



SIMZ-YM2-209
Fig. 74
Miller
Apexo



SIMZ-YM2-210
Fig. 19



SIMZ-YM2-211
Fig. 20

Shorter leverage for lower 1st and 2nd Molar roots. Also useful for premolar and molar deciduous teeth.



lower roots elevator

- Elevator корневой экстракционный
- Elévateurs à racines
- Wurzelheber

- Elevator coronevii ekstraktsionny
- Elevadores para raíces
- المصاعد الجذر

Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



SIMZ-YM2-212
Fig. 38
Serrated
For removal of root fragments

SIMZ-YM2-213
Fig. 39
Serrated



SIMZ-YM2-214
Fig. 190

SIMZ-YM2-215
Fig. 191



SIMZ-YM2-216
Fig. 15
Woodward
Lower Molar roots

SIMZ-YM2-217
Fig. 18
Woodward
Lower Molar roots

ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elevadores para raices
- Wurzelheber
- Elevator corneviu extracționii
- Elevationes para raices
- المضاع الجذر

Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



SIMZ-YM2-218
Fig. 23



SIMZ-YM2-219
Fig. H5
Hourigan
Serrated



SIMZ-YM2-220
Fig. H6
Hourigan
Serrated



SIMZ-YM2-221
Fig. 1C
Coupland
Gouge



SIMZ-YM2-222
Fig. 2C
Coupland
Gouge
Lower Molars



SIMZ-YM2-223
Fig. 3C
Coupland
Gouge
Lower 3rd Molars



SIMZ-YM2-224
Fig. 69W
Woodward Hand
Pressure Chisel



SIMZ-YM2-225
Fig. 31F
Friedman
Craig
Gouge Serrated



SIMZ-YM2-226
Fig. 32F
Friedman
Craig
Gouge Serrated



SIMZ-YM2-227
Fig. 11M
MacMillan
Gouge



SIMZ-YM2-228
Fig. 12M
MacMillan
Gouge



SIMZ-YM2-229
Fig. 13M
MacMillan
Gouge

For Upper and Lower Cuspids

ENGLISH PATTERNS



Upper root elevator with helicoidal tip



SIMZ-YM2-230
Upper and Lower
3rd Molars



SIMZ-YM2-231
Fig. A1C
Coupland
Gouge



SIMZ-YM2-232
Fig. A2C
Coupland
Gouge



SIMZ-YM2-233
Fig. A3C
Coupland
Gouge



SIMZ-YM2-234
Fig. A30
Cryer



SIMZ-YM2-235
Fig. A31
Cryer

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber

- Elevator coronevii ekstraktsionnyy
- Elevadores para raíces
- المعاضد الجذر

ERGO DYNAMIC ELEVATORS

Ergo Dynamic Elevators offer an alternative to standard size elevator handles.



SIMZ-YM2-236
Fig. 301
Apexo



SIMZ-YM2-237
Fig. 301
Apexo
Ergo Dynamic



SIMZ-YM2-238
Fig. 77R Elevator
Serrated
Ergo Dynamic



SIMZ-YM2-239
Fig. B Cogswell
Elevator
Ergo Dynamic




SIMZ-YM2-240
Fig. 34S
Seldin Elevator
Straight
Ergo Dynamic



SIMZ-YM2-241
Luxating Elevator
Straight
3mm
Ergo Dynamic



SIMZ-YM2-242
Luxating Elevator
Curved
3mm
Ergo Dynamic



SIMZ-YM2-243
Luxating Elevator
Straight
5mm
Ergo Dynamic



SIMZ-YM2-244
Luxating Elevator
Curved
5mm
Ergo Dynamic

ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elevator cornevil ekstraktsionnyy
- Elévateurs à racines
- Elevadores para raíces
- Wurzelheber
- المعاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Laxutor Root Elevators Set of 9 Pcs



SIMZ-YM2-245
2mm Str.



SIMZ-YM2-246
3mm Str.



SIMZ-YM2-247
5mm Str.



SIMZ-YM2-248
3mm Cvd.



SIMZ-YM2-249
5mm Cvd.



SIMZ-YM2-250
3mm Angle



SIMZ-YM2-251
3mm Angle



SIMZ-YM2-252
3mm Str.



SIMZ-YM2-253
5mm Str.



- Элеватор корневой экстракционный
 - Elévateurs à racines
 - Wurzelheber
- Elevator cornevei ekstraktsionny
 - Elevadores para raices
 - المعاصر الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Mesial & Distal Root Elevators



Root Elevators



Apexo Root Elevators



ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elevateurs à racines
- Wurzelheber
- Elevator corneviy ekstraktsionnyy
- Elevadores para raices
- المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.


Wisdom Roots Elevator



SIMZ -YM2-265
Wisdom
for 18 and 38

SIMZ -YM2-266
Wisdom
for 28 and 48

Roots Elevator



SIMZ -YM2-268
Fig. 61

SIMZ -YM2-269
Fig. 62

SIMZ -YM2-270
Fig. 63

SIMZ -YM2-267
Fig. 60



Luxier Roots Elevator



SIMZ -YM2-271
Luxier
3mm, Straight

SIMZ -YM2-272
Luxier
3mm, Curved

SIMZ -YM2-273
Luxier
5mm, Straight

SIMZ -YM2-274
Luxier
3mm, Angular

- Elevator корневой экстракционный
- Elévateurs à racines
- Wurzelheber


- Elevator корневой экстракційний
- Elevadores para raices
- المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Hylin Roots Elevator



Molars Roots Elevator



Schlemmer Roots Elevator



ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elevateurs à racines
- Wurzelheber
- Elevator corneviy ekstraktsionnyy
- Elevadores para raices
- المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.


Medan Roots Elevator
















Medan Roots Elevator

							
SIMZ-YM2-284 Bein 2.5mm	SIMZ-YM2-285 Bein 3mm Flat	SIMZ-YM2-286 Bein 3mm Round	SIMZ-YM2-287 Bein 3.5mm	SIMZ-YM2-288 Flohr 3.5mm (Left)	SIMZ-YM2-289 Flohr 3.5mm (Right)	SIMZ-YM2-290 Cryer 10mm (Left)	SIMZ-YM2-291 Cryer 10mm (Right)
							
	SIMZ-YM2-292 Seldin 4L (Left)	SIMZ-YM2-293 Seldin 4R (Right)		SIMZ-YM2-294 Seldin 1L (Left)	SIMZ-YM2-295 Seldin 1R (Right)		
							
							SIMZ-YM2-296 Seldin

Roots Elevator



							
SIMZ-YM2-297 Bein 2.5mm	SIMZ-YM2-298 Bein 3mm Flat	SIMZ-YM2-299 Bein 3mm Round	SIMZ-YM2-300 Bein 3.5mm	SIMZ-YM2-301 Flohr 3.5mm (Left)	SIMZ-YM2-302 Flohr 3.5mm (Right)	SIMZ-YM2-303 Cryer 10mm (Left)	SIMZ-YM2-304 Cryer 10mm (Right)
							
	SIMZ-YM2-305 Seldin 4L (Left)	SIMZ-YM2-306 Seldin 4R (Right)		SIMZ-YM2-307 Seldin 1L (Left)	SIMZ-YM2-308 Seldin 1R (Right)		
							
							SIMZ-YM2-309 Seldin

ROOT ELEVATOR



- Элеватор корневой экстракционный
- Elevator cornevier ekstrakcyjny
- Elévateurs à racines
- Elevadores para raices
- Wurzelheber
- المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.



Root Elevators



SIMZ-YM2-310 Bein 2mm round
 SIMZ-YM2-311 Bein 3mm round
 SIMZ-YM2-312 Bein 4mm round
 SIMZ-YM2-313 Bein 2mm Flat
 SIMZ-YM2-314 Bein 3mm Flat
 SIMZ-YM2-315 Bein 4mm Flat
 SIMZ-YM2-316 White 4mm
 SIMZ-YM2-317 White 4mm



SIMZ-YM2-318 Berten 5mm Lower roots
 SIMZ-YM2-319 3mm
 SIMZ-YM2-320 3mm
 SIMZ-YM2-321 3mm
 SIMZ-YM2-322 3mm



SIMZ-YM2-323 Flohr 3mm



SIMZ-YM2-324 Flohr 3mm
 SIMZ-YM2-325 Flohr 3mm Flat
 SIMZ-YM2-326 Flohr 3mm Flat
 SIMZ-YM2-327 Cryer 10mm
 SIMZ-YM2-328 Cryer 10mm

FUNCTION: To loosen tooth or root from bony socket prior to placement of the extraction forceps.





ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Elevator corneviu extracționii
- Elevadores para raíces
- المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Kopp Roots Elevator



ROOT ELEVATOR



- Элеватор корневой экстракционный
- Elevator corneveiy ekstraktsionnyiy
- Elévateurs à racines
- Elevadores para raíces
- Wurzelheber
- المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Kopp Roots Elevator



SIMZ-YM2-346
Kopp
Right

SIMZ-YM2-347
Kopp
Left

Heber Roots Elevator



SIMZ-YM2-348
Heber
Bayonet

SIMZ-YM2-349
Heber
Bayonet


Bein Roots Elevator



SIMZ-YM2-350
Bein
3mm

SIMZ-YM2-351
Bein
4mm

Bein Roots Elevator



SIMZ-YM2-352 Bein 2mm, curved serrated	SIMZ-YM2-356 Bein 2mm, curved serrated	SIMZ-YM2-360 Bein 2mm, curved serrated
--	--	--

SIMZ-YM2-353 2mm straight serrated	SIMZ-YM2-354 3mm straight	SIMZ-YM2-355 4mm Curved
SIMZ-YM2-357 2mm straight serrated	SIMZ-YM2-358 3mm straight	SIMZ-YM2-359 4mm Curved
SIMZ-YM2-361 2mm straight serrated	SIMZ-YM2-362 3mm straight	SIMZ-YM2-363 4mm Curved

ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Слеватор корневий екстракційний
- Elevadores para raíces
- المصاعد الجذر

Atraumatic Extraction Bayonet Root Elevators



SIMZ-YM2-364
3mm Luxating
Elevator Curved



SIMZ-YM2-365
3mm Luxating
Elevator Straight



SIMZ-YM2-366
5mm Luxating
Elevator Curved



SIMZ-YM2-367
5mm Luxating
Elevator Straight



SIMZ-YM2-368
2mm Grooved
Luxating Elevator



SIMZ-YM2-369
3mm Grooved
Luxating Elevator



SIMZ-YM2-370
4mm Grooved
Luxating Elevator



SIMZ-YM2-371
5mm Grooved
Luxating Elevator

- Элеватор корневой экстракционный
- Élévateurs à racines
- Wurzelheber

- Єлеватор кореневий екстракційний
- Elevadores para raices
- المعاصر الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Luxating Elevator With Periotome Tip



SIMZ-YM2-372
3mm Luxating
Elevator / Periotome

SIMZ-YM2-373
4mm Luxating
Elevator / Periotome

SIMZ-YM2-374
5mm Luxating
Elevator / Periotome

Atraumatic Extraction Bayonet Root Elevators



SIMZ-YM2-375
Bayonet Elevator
Straight

SIMZ-YM2-376
Bayonet Elevator
Offset "In" Small

SIMZ-YM2-377
Bayonet Elevator
Offset "Out" Small

SIMZ-YM2-378
Bernard Bayonet
Elevator Offset Left

SIMZ-YM2-379
Bernard Bayonet
Elevator Offset Right

ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elevadores para raíces
- Wurzelheber
- Elevator хореновий экстракційний
- Elevadores para raíces
- المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.


W. Barry Roots Elevator



Winter Roots Elevator




Pott Roots Elevator



- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Elevator cornevil extracționii
- Elevadores para raíces
- المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Heidbrink Roots Elevator



SIMZ-YM2-395 SIMZ-YM2-396 SIMZ-YM2-397 SIMZ-YM2-398
Fig. 1

Heidebrink Root tip Picks




SIMZ-YM2-399 Fig. 2 Left
SIMZ-YM2-400 Fig. 3 Right

Syndesmotomes Roots Elevator




SIMZ-YM2-401 Straight
SIMZ-YM2-402 Left
SIMZ-YM2-403 Right

Periotome Roots Elevator



SIMZ-YM2-404 Periotome 2.1mm Straight
SIMZ-YM2-405 Periotome 2.6mm Straight
SIMZ-YM2-406 Periotome 2.1mm Angular
SIMZ-YM2-407 Periotome 2mm Curved

Syndesmotomes




SIMZ-YM2-408 Periotome 2.0mm
SIMZ-YM2-409 Periotome 2.5mm
SIMZ-YM2-410 Periotome 3.0mm
SIMZ-YM2-411 Periotome 4.0mm

ROOT ELEVATOR

- Элеватор корневой экстракционный
 - Elévateurs à racines
 - Wurzelheber
- Єлеватор корневий екстракційний
 - Elevadores para raíces
 - المعصاع الجذر

FUNCTION: To loosen small root fragments from bony socket



- Элеватор корневой экстракционный
- Élévateurs à racines
- Wurzelheber

- Elevator cornevil ekstraktsionniy
- Elevadores para raices
- المصاعد الجذر

Function of Root Elevators: To loosen tooth or root from bony socket prior to placement of the extraction forceps.



Heidbrink universal mini-elevator for fragments

HOWARD

Slightly heavier design than Heidbrink picks.

HEIDBRINK


APICAL



ROOT ELEVATOR

- Elevator корневой экстракционный
- Elevator coronevii ekstraktsionnyy
- Elévateurs à racines
- Elevadores para raíces
- Wurzelheber
- المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.



- Элеватор корневой экстракционный
- Élévateurs à racines
- Wurzelheber

- Elevator corneveiey extraicéionnyy
- Elevadores para raices
- المصاعد الجذر



Lecluse elevator for lower third molar dislocation



SIMZ-YM2-449
Fig. 1
Chompret



SIMZ-YM2-450
Fig. 2
Chompret



SIMZ-YM2-451
Fig. 3
Chompret



SIMZ-YM2-452
Fig. 4
Chompret



Chompret syndesmotome, sickle pattern

CROSS-BAR HANDLE



SIMZ-YM2-453
Fig. 11LX

SIMZ-YM2-454
Fig. 11RX



SIMZ-YM2-455
Fig. 12LX



SIMZ-YM2-456
Fig. 12RX



SIMZ-YM2-457
Fig. 14LX



SIMZ-YM2-458
Fig. 14RX



SIMZ-YM2-459
Fig. 1LX



SIMZ-YM2-460
Fig. 1RX

ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Elevator corneviu extracționii
- Elevadores para raíces
- المصاعد الجذر

Function of Root Elevators: To loosen tooth or root from bony socket prior to placement of the extraction forceps.



SIMZ-YM2-461
775/3
Winter
11 L



SIMZ-YM2-462
775/4
Winter
11 R



SIMZ-YM2-463
775/5
Winter
12 L



SIMZ-YM2-464
775/6
Winter
12 R



SIMZ-YM2-465
775/7
Winter
13 L



SIMZ-YM2-466
775/8
Winter
13 R



SIMZ-YM2-467
winter
Left

SIMZ-YM2-468
winter
right



SIMZ-YM2-469
775/1
Winter
1 L



SIMZ-YM2-470
775/2
Winter
1 R



SIMZ-YM2-471
Winter
1L




SIMZ-YM2-472
Winter
1R

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber

- Elevator corneevy ekstraktsionnyy
- Elevadores para raíces
- المعصاع الجذر

LUXACUT



Luxacut Luxation Instrument

Over recent years, there has been an increased demand for reduced trauma during extraction procedures. The use of elevators has become widespread for the loosening of teeth before the application of forceps. Luxation instruments tips are more finely tapered than elevators tips. Their characteristics are a combination of a Periostome and an elevator. The sharp tip enables more efficient cutting of the periodontal ligament, whilst the tapering shaft assists with controlled expansion of the bony socket. The extraction operation may be performed with the minimum amount of tissue damage. Healing is therefor quicker and the operation is less distressing for the patient and surgical team alike.

The Luxacut Kit comprises of seven Luxation instruments, a straight, curved, inverted curved & contra blade version in Three different widths. Also available in dual edge blades

- SIMZ-YM2-473-L2S Straight Blade- Green-2 mm- Apical and/or Interproximal
- SIMZ-YM2-475-L3C Curved Blade- Black- 3 mm General use/Starter
- SIMZ-YM2-474-L3S Straight Blade- Grey- 3 mm Interproximal
- SIMZ-YM2-477-L3CA Contra Angle- Blue- 3 mm Lingual and/or distal-Molars
- SIMZ-YM2-476-L3IC Inverted Curved- Light Blue- 3 mm- General, Lingual and/or Distal
- SIMZ-YM2-479-L5C Curved Blade- Brown- 5 mm- Large molar roots- General use
- SIMZ-YM2-478-L5S Straight Blade- Yellow- 5 mm- Interproximal
- SIMZ-YM2-480-DE3 dual edge/straight grey 3 mm/1.5 mm
- SIMZ-YM2-481-DE5 dual edge/straight brown 5 mm Deep fractured/Decayed Roots

Luxacut Kit



SIMZ-YM2-482
LWC7



SIMZ-YM2-483
LGWC7

ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Elevator corneviy ekstraktsionnyy
- Elevadores para raices
- المعصاة الجذر

-YM2-484

LUXATING ELEVATORS SET OF 10 PCS



SIMZ-YM2-495

LUXATING ELEVATORS SET OF 10 PCS - BLUE TITANIUM COATED TIPS



- Элеватор корневой экстракционный
- Elevateurs à racines
- Wurzelheber
- Elevator corneviy ekstraktsionniy
- Elevadores para raices
- المصاعد الجذر

SIMZ-YM2-506

LUXATING ELEVATORS SET OF 10 PCS - GOLD COATED TIPS



SIMZ-YM2-517

LUXATING ELEVATORS SET OF 08 PCS - BLUE TITANIUM COATED TIPS



ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Elevator corneviy ekstraktsionnyy
- Elevadores para raices
- المصاعد الجذر

SIMZ-YM2-526

PERIOTOME LUXATING ELEVATORS SET OF 04 PCS - GOLD COATED TIPS



SIMZ-YM2-527
PERIOTOME GOLD TIP - L3S



SIMZ-YM2-528
PERIOTOME GOLD TIP - De3



SIMZ-YM2-529
PERIOTOME GOLD TIP - L2S



SIMZ-YM2-530
PERIOTOME GOLD TIP - L3CA



SIMZ-YM2-531 (Double Curved 1.5mm)
SIMZ-YM2-532 (Double Curved 2mm)
SIMZ-YM2-533 (Double Curved 3mm)
SIMZ-YM2-534 (Double Curved 4mm)
SIMZ-YM2-535 (Double Curved 5mm)



SIMZ-YM2-536
(Inverted Curved 3mm)




SIMZ-YM2-537
(Dual Edge 3mm / 1.5mm)



SIMZ-YM2-538
(Dual Edge 5mm / 3mm)

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Elevator corneviu extracționii
- Elevadores para raíces
- المصاعد الجذر

PDL ELEVATORS



ROOT ELEVATOR

- Elevator корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Elevator corneveiy ekstraktsionnyiy
- Elevadores para raices
- المصاعد الجذر



Area: distal




Area: mesial



SIMZ-YM2-541 PDL Elevators Set of 10 Pcs



			SIMZ-YM2-542 Luxa-Tool, straight, 3.0 mm,
			SIMZ-YM2-543 Luxa-Tool, curved, 3.0 mm,
			SIMZ-YM2-544 Proximator, curved, 2.5 mm, distal elevator,
			SIMZ-YM2-545 Proximator, curved, 2.5 mm, mesial elevator,
			SIMZ-YM2-546 Proximator, curved, 4.5 mm, distal elevator,
			SIMZ-YM2-547 Proximator, curved, 4.5 mm, mesial elevator,
			SIMZ-YM2-548 Syndesmotome, curved, 4.0 mm, for distal access
			SIMZ-YM2-549 Syndesmotome, curved, 4.0 mm, for mesial access
			SIMZ-YM2-550 Syndesmotome, curved left, 4.0 mm, for lingual/buccal luxation,
			SIMZ-YM2-551 Syndesmotome, curved right, 4.0 mm, for lingual/buccal luxation,

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber

- Elevator corneviy ekstraktsionnyy
- Elevadores para raíces
- المعصاع الجذر

PDL ELEVATORS

Delicate Arrow tip



4,0mm straight
4,0mm curved

Luxivators straight
atraumatic luxation of the tooth



2,5mm straight
3,5mm straight
4,5mm straight

Luxivators curved
atraumatic luxation of the tooth



2,5mm curved
3,5mm curved
4,5mm curved



SIMZ-YM2-552
PDL Elevators
Set of 8 Pcs

Dental Hand Instruments

The types of hand instruments can be categorized based on their function:

- Examining the mouth and teeth
- Scaling or professional cleaning of teeth
- Cutting teeth and removing cavities
- Placing and condensing filling materials
- Carving and finishing dental fillings
- Miscellaneous types
- Instruments used for examining the mouth and teeth

Mouth mirrors are used to examine the condition of the teeth and mouth, as well as to retract the tongue and cheeks upon examination.

Probes— There are many probes available but the main types are the straight probe, Briault probe and periodontal probe.

The straight probe is sharp and it is used to examine cavities in the dentine region as well as to check the margins of restorations or fillings.

The Briault probe is also sharp and it is used to detect the junction between the enamel and dentine. Briault can also be used to detect tartar hidden in the gum pockets.

A periodontal probe can be used to measure the depth of the gum pockets and for other measurements. Instruments for cutting and scraping tissues and deposits

Scalers remove tartar (calculus) on the tooth surface as well as those hidden in the gum pockets. They can also be used to remove other surface deposits and temporary crowns.

Excavators are used to remove softened dentine and temporary fillings. The back of the blade can also be used for placing linings between filling materials. Sometimes excavators can be used for carving metal or amalgam fillings.

Chisels are used to remove unsupported enamel during cavity preparation. Unsupported enamel can be harmful for it can compromise the longevity of the filling if not removed.

Hatchets and hoes have similar purposes to chisels. Instruments used for placing and condensing filling materials. Plastic instruments convey and shape filling materials and are usually flat blades that do not use heavy pressure.

Some instruments have rounded ends for shaping and burnishing the surfaces of fillings. Plastic instruments are usually made of stainless steel for most materials though there are also thin Teflon-coated or titanium nitride instruments which are non-stick for composite fillings. **Condensers or pluggers** are used for compressing and forming filling materials. Heavy pressure is applied on use therefore these instruments are usually used for metal or amalgam fillings. Instruments used for carving and finishing fillings. Carving and finishing instruments have sharp or semi-sharp blades of various shapes. They are used to carve materials by cutting and scraping.



Symptoms of Gum and Dental Problems



VALUED BUYERS

Compare our scalers and curettes with the best in the industry and you will immediately recognize we are equal in every way but one... price. Note the deep precision knurling for better tactile feel. The same high-grade stainless steel tips, the same light weight (under 20 grams) for less hand fatigue, the same professional craftsmanship, and surpassing the highest industry standards. Of course, any manufacturer can make these claims, but we guarantee and certify ours. We regularly manufacturing below articles in showing handles & pattern A to L for better grip & perfection. We are also fully equipped & capable to develop these instruments as per your own specific requirements & demands in materials, sizes, shapes and patterns.



A-Hollow Handle Ø 10 mm



B-Hollow Handle Ø 10 mm



C-Hollow Handle Ø 10 mm



D-Hollow Handle Ø 10 mm



E-Hollow Handle Ø 10 mm



F-Hollow Handle Ø 10 mm



G-Hollow Handle Ø 10 mm



H-Hollow Handle Ø 6.4 mm



I-Hollow Handle Ø 8 mm



J-Hollow Handle Ø 10 mm



K-Soild Handle Ø 6 mm



L-Soild Handle Ø 6 mm



M-Hollow Handle Ø 10 mm



N-Hollow Handle Ø 10 mm



O-Hollow Handle Ø 10 mm



P-Hollow Handle Ø 10 mm



Q-Hollow Handle Ø 10 mm



R-Hollow Handle Ø 10 mm



S-Hollow Handle Ø 10 mm



T-Hollow Handle Ø 10 mm



U-Hollow Handle Ø 10 mm



V-Hollow Handle Ø 10 mm



W-Hollow Handle Ø 10 mm

MOUTH MIRRORS

- зеркала
- Espejos de la boca
- Mundspiegel

- дзеркала
- Miroirs à la bouche
- مرايا الفم



SIMZ-YM3-001
Set of 12 Pcs

Plain



SIMZ-YM3-002 Fig. 3 Ø20mm
SIMZ-YM3-003 Fig. 4 Ø22mm
SIMZ-YM3-004 Fig. 5 Ø24mm
SIMZ-YM3-005 Fig. 6 Ø26mm



SIMZ-YM3-006
Set of 12 Pcs

Cone-Socket Plain



SIMZ-YM3-007 Fig. 3 Ø20mm
SIMZ-YM3-008 Fig. 4 Ø22mm
SIMZ-YM3-009 Fig. 5 Ø24mm



SIMZ-YM3-010
Set of 12 Pcs

Magnifying



SIMZ-YM3-011 Fig. 3 Ø20 mm
SIMZ-YM3-012 Fig. 4 Ø22 mm
SIMZ-YM3-013 Fig. 5 Ø24 mm



SIMZ-YM3-014
Set of 12 Pcs

Cone-Socket Magnifying



SIMZ-YM3-015 Fig. 3 Ø20mm
SIMZ-YM3-016 Fig. 4 Ø22mm
SIMZ-YM3-017 Fig. 5 Ø24mm

MOUTH MIRRORS

- зеркала
- Espejos de la boca
- Mundspiegel
- дзеркала
- Miroirs à la bouche
- مرايا الفم



SIMZ-YM3-018
Set of 12 Pcs



SIMZ-YM3-022
Set of 12 Pcs

Rhodium Front Surface



SIMZ-YM3-019 Fig. 3 Ø20mm
SIMZ-YM3-020 Fig. 4 Ø22mm
SIMZ-YM3-021 Fig. 5 Ø24mm

Cone-Socket Rhodium



SIMZ-YM3-023 Fig. 3 Ø20mm
SIMZ-YM3-024 Fig. 4 Ø22mm
SIMZ-YM3-025 Fig. 5 Ø24mm



SIMZ-YM3-026
Set of 12 Pcs



SIMZ-YM3-031
Set of 12 Pcs

Rhodium Front Surface



SIMZ-YM3-027 Fig. 3 Ø20mm
SIMZ-YM3-028 Fig. 4 Ø22mm
SIMZ-YM3-029 Fig. 5 Ø24mm
SIMZ-YM3-030 Fig. 6 Ø26mm

Cone-Socket Front Surface



SIMZ-YM3-032 Fig. 3 Ø20mm
SIMZ-YM3-033 Fig. 4 Ø22mm
SIMZ-YM3-034 Fig. 5 Ø24mm

HANDLES FOR MOUTH MIRRORS

- Ручки для зеркал
- Manches pour miroirs à bouche
- Mundspiegelgriffe

- Ручки для дзеркал
- Mangos para espejos de boca
- مقابض لمرآة الفم



SIMZ-YM3-035
Hollow Mirror
Handle Round
Ø 10 mm



SIMZ-YM3-036
Hollow Mirror
Handle
Hexagonal
Ø 8 mm



SIMZ-YM3-037
Hollow
Ø 8 mm



SIMZ-YM3-038
Hollow Mirror
Handle
Hexagonal
Ø 10 mm



SIMZ-YM3-039
Hollow Mirror
Handle
Ø 8 mm



SIMZ-YM3-040
Hollow Mirror
Handle
Ø 10 mm



SIMZ-YM3-041
Hollow Mirror
Handle
Ø 8 mm



SIMZ-YM3-042
Mirror Handle
Autoclavable
Material
Available in
Various Colours
Ø 9.5 mm



SIMZ-YM3-043
Mirror Handle
Autoclavable
Material
Available in
Various Colours
Ø 9.5 mm



SIMZ-YM3-044
Mirror Handle
Silicon Grip
Autoclavable
Available in
Various Colours



SIMZ-YM3-045
Mirror Handle
Solid
Ø 6 mm



SIMZ-YM3-046
Mirror Handle
Solid
Ø 6 mm



SIMZ-YM3-047
Mirror Handle
Solid Octagonal
Ø 6 mm



SIMZ-YM3-048
Mirror Handle
Solid
Stem with
Rule 35mm



SIMZ-YM3-049
Cone-Socket
Solid Octagonal
Ø 6 mm



SIMZ-YM3-050
Cone-socket
Solid
Ø 6 mm



SIMZ-YM3-051
Cone-socket
Hexagonal
Handle
Ø 6 mm

Cone-Socket Mirror Handle

TWEEZERS

- Пинцеты
- Précelles - Bruxelles
- Zahnpinzetten
- Пинцети
- Pinzas para algodón
- ملاقط



SIMZ-YM3-052
London College
150mm



SIMZ-YM3-053
Meriam
160mm



SIMZ-YM3-054
Semken Taylor
125mm
Straight



SIMZ-YM3-055
Semken Taylor
125mm
Curved



SIMZ-YM3-056
Semken Taylor
125mm
1x2 Straight



SIMZ-YM3-057
Semken Taylor
125mm
1x2 Curved



SIMZ-YM3-058
Flagg
150mm



SIMZ-YM3-059
Meriam
160mm



SIMZ-YM3-060
Amsterdam
155mm



SIMZ-YM3-061
Perry
130mm
Curved



SIMZ-YM3-062
Self-Locking
150mm
Angled



SIMZ-YM3-063
Endo-Locking
Grooved
150mm

- Пинцеты
- Prêcelles - Brucelles
- Zahnpinzetten
- Пинцети
- Pinzas para algodón
- ملقح

Zahnpinzetten



SIMZ-YM3-064 90mm
SIMZ-YM3-065 120mm
SIMZ-YM3-066 150mm
SIMZ-YM3-067 180mm
Taylor



SIMZ-YM3-068 90mm
SIMZ-YM3-069 120mm
SIMZ-YM3-070 150mm
SIMZ-YM3-071 180mm
Taylor

1x2
ø mm 0,9



SIMZ-YM3-072 0,6mm
SIMZ-YM3-073 0,9mm
Gerald
175mm

1x2
ø mm 0,9
ø mm 0,6



SIMZ-YM3-074 0,6mm
SIMZ-YM3-075 0,9mm
Gerald
175mm

1x2
ø mm 0,9
ø mm 0,6



SIMZ-YM3-076
Gerald
175mm



SIMZ-YM3-077
Gerald
175mm



SIMZ-YM3-078
150mm
College



SIMZ-YM3-079
150mm
Meriam



SIMZ-YM3-080
De Backy
Atraumatic
150mm
Straight



SIMZ-YM3-081
De Backy
Atraumatic
150mm
Angled



SIMZ-YM3-082
De Backy
Atraumatic
175mm



SIMZ-YM3-083
Gilles
150mm 1x2

TWEEZERS

- Пинцеты
- Précelles - Brucelles
- Zahnpinzetten
- Пинцети
- Pinzas para algodón
- منجلج



SIMZ-YM3-084
Meriam Tweezers
smooth 16cm



SIMZ-YM3-085
London College
smooth 16cm



Adson With Teeth
SIMZ-YM3-086 12cm
SIMZ-YM3-087 15cm



SIMZ-YM3-088
Perry
10.5 cm



SIMZ-YM3-089
Perry
10.5 cm



SIMZ-YM3-090
Fig. 1
Perry
13 cm



Potts-Smith Dressing Forceps
SIMZ-YM3-091 18 cm
SIMZ-YM3-092 20 cm
SIMZ-YM3-093 25 cm
straight serrated



Potts-Smith Dressing Forceps
SIMZ-YM3-094 18 cm
SIMZ-YM3-095 20 cm
SIMZ-YM3-096 25 cm
curved serrated



Potts-Smith Tissue Forceps 1x2 teeth
SIMZ-YM3-097 18 cm
SIMZ-YM3-098 20 cm
SIMZ-YM3-099 25 cm
straight



Potts-Smith Tissue Forceps 1x2 teeth
SIMZ-YM3-100 18 cm
SIMZ-YM3-101 20 cm
SIMZ-YM3-102 25 cm
curved



De-Bakey With Tungsten Carbide (TC)
SIMZ-YM3-103 17.5cm



Cushing Taylor 18cm
SIMZ-YM3-104



TWEEZERS

- Пинцеты
- Pinças para algodón
- Pinçettes - Brucelles
- Zahnpinzetten
- ملاقف



SIMZ-YM3-105
Adson-TC
120mm



SIMZ-YM3-106
Adson-TC
120mm 1x2



SIMZ-YM3-107
Adson-TC
150mm 1x2



SIMZ-YM3-108
Adson
120mm



SIMZ-YM3-109
Adson
120mm 1x2



SIMZ-YM3-110
Adson
150mm 1x2



SIMZ-YM3-111
Adson Brown
120mm



SIMZ-YM3-112
Crane Kaplan
150mm



SIMZ-YM3-113 115mm
SIMZ-YM3-114 140mm
Dressing Forceps



SIMZ-YM3-115 115mm
SIMZ-YM3-116 140mm
Tissue Forceps



SIMZ-YM3-117
Hemingway
150mm



SIMZ-YM3-118
Crane Kaplan
150mm



SIMZ-YM3-119
Crane Kaplan
150mm



SIMZ-YM3-120
Gerald
175mm

TWEEZERS

- Пинцеты
- Précelles - Brucelles
- Zahnpinzetten
- Пинцети
- Pinzas para algodón
- ملقاط



Adson-Brown



Adson-Brown



Brown




Ewald / Hudson



Adson /
Micro-Adson



Adson
T.C



TWEEZERS

- Пинцеты
- Pinças para algodón
- Pinzette
- Prêcelles - Brucelles
- Zahnpinzetten
- ملاقط



SIMZ-YM3-134
Adson
1 mm
12 cm

SIMZ-YM3-135
Adson T.C
1,7 mm
12 cm

SIMZ-YM3-136
Adson
1,7 mm
15 cm



SIMZ-YM3-137
Micro-Adson
very fine
0,8 mm
15 cm



SIMZ-YM3-138
Adson
1x2, 1,5 mm
12 cm



SIMZ-YM3-139
Micro-Adson
1x2, 1 mm
15 cm
very fine



SIMZ-YM3-140
Micro-Adson
1x2, 1 mm
15 cm
curved
very fine



SIMZ-YM3-141
Cushing
1,2 mm
18 cm



SIMZ-YM3-142 13 cm
1x2, 1,5 mm
SOS-YM3-143 14 cm
1x2, 1,7 mm



SIMZ-YM3-144
Taylor
1x2, 1,5 mm
17,5 cm



SIMZ-YM3-145 13 cm
1x2, 2,3 mm
SOS-YM3-146 14 cm
1x2, 2,3 mm



SIMZ-YM3-147
Micro-Cushing
1x2, 0,9 mm
17,5 cm
very fine



SIMZ-YM3-148
Suture Forceps
15 cm

TWEEZERS

- Пинцеты
- Précelles - Brucelles
- Zahnpinzetten
- Пинцети
- Pinzas para algodón
- ملقف



SIMZ-YM3-149
15 cm
Fig. 1
Thin Pattern



SIMZ-YM3-150
15 cm
Fig. 2
Wide Pattern



SIMZ-YM3-151
16 cm



SIMZ-YM3-152
15 cm
Fig. 1



SIMZ-YM3-153
15 cm
Fig. 2



SIMZ-YM3-154
11.5 cm
Fig. 1

Micro Surgery Instrument



SIMZ-YM3-155
18 cm
Fig. 1



SIMZ-YM3-156
18 cm
Fig. 2



SIMZ-YM3-157
18 cm
Fig. 3



SIMZ-YM3-158
18 cm
Fig. 4

SOLDERING PLIERS - TWEEZERS & SPLINTER FORCEPS

- Пинцеты для удаления осколков
- Pincettes pour la vidalennia oskolokiv
- Pincas à souder - Pincas et pincas à échardes
- Alicates Soldadura - Pinzas y pinzas de Splinter
- Löt-Zange - Pinzette & SplinterZange
- لحام كعاشة - ملاقط منشفة الملقط



SIMZ-YM3-159
Soldering Pliers
160mm



SIMZ-YM3-160
Soldering Pliers
160mm



SIMZ-YM3-161
Posigrip Articulating
Paper Forceps
150mm



SIMZ-YM3-162
Standard Articulating
Paper Forceps, Economy Grade
150mm



SIMZ-YM3-163
Miller
150mm
Curved



SIMZ-YM3-164
125mm



SIMZ-YM3-165
125mm



SIMZ-YM3-166
115mm



SIMZ-YM3-167
125mm



SIMZ-YM3-168
125mm



SIMZ-YM3-169
Peet
115mm

EXAMINATION SETS

- Комплект для обследования
- Экзаменаційний набір
- Ensemble d'examen
- Conjunto de examen
- Prüfungssatz
- مجموعة الفحص

SIMZ-YM3-170

First Aid Examination kit Set of 3 Pieces



SIMZ-YM3-171

First Aid Examination Set of 4 Pieces



SIMZ-YM3-172

First Aid Examination Set of 3 Pieces




SIMZ-YM3-173

First Aid Examination Set of 3 Pieces



SIMZ-YM3-174

First Aid Examination Set of 3 Pieces



SIMZ-YM3-175

First Aid Examination Set of 3 Pieces



EXAMINATION SETS

- Комплект для обследования
- Экзаменационный набор
- Ensemble d'examen
- Conjunto de examen
- Prüfungssatz
- مجموعة الفحص



SIMZ-YM3-176
First Aid Examination kit Set of 3 Pieces




SIMZ-YM3-177
First Aid Examination Set of 4 Pieces




SIMZ-YM3-178
First Aid Examination Set of 3 Pieces



SIMZ-YM3-179
First Aid Examination Set of 3 Pieces



SIMZ-YM3-180
First Aid Examination Set of 3 Pieces




SIMZ-YM3-181
First Aid Examination Set of 3 Pieces




EXPLORERS

- Зонды
- Sondes
- Zahnsonden
- Зонди
- Sondas
- المستكشفين


Single Ended Probe Solid Handle Ø 6 mm, Length 150mm




Single Ended Probe Hollow Handle Ø 8 mm Length 150mm



**EXPLORERS INSTRUMENTS AVAILABLE IN
Hollow Handle Ø 8 mm, Ø 10mm & Solid Handle Ø 6mm
Length 170MM**




Available in Hollow Handle
 Ø 8 mm, Ø 10mm & Solid Handle Ø 6mm



EXPLORERS


- Зонды
- Sondes
- Zahnsonden
- Зонди
- Sondas
- المستكشفين



Pointed Probes Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm




Pointed Probes Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm



- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Инструменти пародонтологічні
- instrumentos para periodoncia
- الصغوك اللثة


SIMZ-YM3-248
 Periodontal Set of 10 Pieces Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm8



SIMZ-YM3-249
 Mouth Mirror

SIMZ-YM3-250 EXP,23/17-A SIMZ-YM3-251 Exavators SIMZ-YM3-252 Filling Instruments SIMZ-YM3-253 EXP,23/CP-18 SIMZ-YM3-254 EXP,23/CP-QW SIMZ-YM3-255 EXP,23/CP-11 SIMZ-YM3-256 EXP,23/CP-10 SIMZ-YM3-257 EXP,23/CP-12 SIMZ-YM3-258 EXP,23/CP-11.5-B

SIMZ-YM3-259
William Probe Set Double Ended With Solid Handle Ø 5.5 mm, Length 170mm



SIMZ-YM3-260
 EXP, 23/CO,QW

SIMZ-YM3-261
 EXP, 23/CP,10

SIMZ-YM3-262
 EXP, 23/CP,18

SIMZ-YM3-263
 EXP, 23/CP,12

SIMZ-YM3-264
 EXP, 23/CP,2


SIMZ-YM3-265
 EXP, 23/UNC,15

SIMZ-YM3-266
 EXP, 23/CP,11

SIMZ-YM3-267
 EXP, 23/CP, 12 ,S

Explorers Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm

PLUGGERS




- Зонды пародонтологические
- Sondes paradontales
- sonde paradontali
- зонди пародонтологічні
- Sondas paradontales
- المعايير التتة

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm


SPREADERS

HEAT CARRIERS




Endodontics Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm

ANTERIOR SINGLE ENDED



SIMZ-YM3-294
8A
29mm

LUKS



SIMZ-YM3-303
L1
18mm

SIMZ-YM3-304
L2
18mm

SIMZ-YM3-305
L3
18mm


SIMZ-YM3-306
L4
18mm

PERIODONTAL PROBES

- Зонды пародонтологические
- Sondas paradontales
- sonde parodontali
- зонди пародонтологічні
- Sondas paradontales
- المعاسير اللثة

Endodontics Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm
POSTERIOR SINGLE ENDED

RETRO FILLING PLUGGERS




ROOT CANAL SPREADER & AMALGAM INSTRUMENTS

- Инструменты для расширения корневых каналов и амальгамы
- Épandeur de canal radiculaire et instruments d'amalgame
- Wurzelkanalverteiler & Amalgaminstrumente

- Инструменти для розповсюдження корневих каналів та амальгами
- Esparcidor de conducto radicular e instrumentos de amalgama
- آلة نشر قناة الجذر وانوات الملعق

Silicone Handles



SIMZ-YM3-319
Explorer 1

SIMZ-YM3-320
Explorer 2

SIMZ-YM3-321
Explorer 3

SIMZ-YM3-322
Root Canal Spreader 1

SIMZ-YM3-323
Root Canal Spreader 2

SIMZ-YM3-324
Root Canal Spreader 3

SIMZ-YM3-325
Root Canal Spreader 4



SIMZ-YM3-326
Amalgam
Instruments 1



SIMZ-YM3-327
Amalgam
Instruments 2



SIMZ-YM3-328
Amalgam
Instruments 3



SIMZ-YM3-329
Amalgam
Instruments 4




SIMZ-YM3-330
Amalgam
Instruments 5



SIMZ-YM3-331
Amalgam
Instruments 6

PROBES



SIMZ-YM3-332
Probe

SIMZ-YM3-333
Probe

SIMZ-YM3-334
Probe

SIMZ-YM3-335
Probe

SIMZ-YM3-336
Probe


SIMZ-YM3-337
Probe

SIMZ-YM3-338
Probe

SIMZ-YM3-339
Probe

PERIODONTAL INSTRUMENTS

- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Інструменти пародонтологічні
- instrumentos para periodoncia
- المسكوك للثة



SIMZ-YM3-340
Excavator

SIMZ-YM3-341
Excavator

SIMZ-YM3-342
Excavator

SIMZ-YM3-343
Filling
Instruments

SIMZ-YM3-344
Filling
Instruments


SIMZ-YM3-345
Filling
Instruments

SIMZ-YM3-346
Filling
Instruments

SIMZ-YM3-347
Filling
Instruments

SIMZ-YM3-348
Cement Spatula

MODELLING INSTRUMENTS



SIMZ-YM3-349
Modelling
Instruments

SIMZ-YM3-350
Modelling
Instruments

SIMZ-YM3-351
Modelling
Instruments

SIMZ-YM3-352
Modelling
Instruments

SIMZ-YM3-353
Modelling
Instruments


SIMZ-YM3-354
Modelling
Instruments

SIMZ-YM3-355
Modelling
Instruments

SIMZ-YM3-356
Modelling
Instruments

- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Инструменти пародонтологічні
- instrumentos para periodoncia
- الصلوك اللثة

Periodontal Instruments Set




SIMZ-YM3-357
Set of 10 Pieces.

CURETTES

- Кюрети
- Curetas
- Kuretten
- Кюрети
- Curette
- كرتين

HOLLOW HANDLE PLASMA BLUE COATED Ø 10 MM, LENGTH 170MM



SIMZ-YM3-358
GRACEY
1/2

SIMZ-YM3-359
GRACEY
3/4

SIMZ-YM3-360
GRACEY
5/6

SIMZ-YM3-361
GRACEY
7/8

SIMZ-YM3-362
GRACEY
9/10

SIMZ-YM3-363
GRACEY
11/12


SIMZ-YM3-364
GRACEY
13/14

SIMZ-YM3-365

Gracey Curettes Set of 7 Pcs

- Gracey Curettes 1/2
- Gracey Curettes 3/4
- Gracey Curettes 5/6
- Gracey Curettes 7/8
- Gracey Curettes 9/10
- Gracey Curettes 11/12
- Gracey Curettes 13/14

HOLLOW HANDLE Ø 10 MM, LENGTH 170MM WITH SILICON RING



SIMZ-YM3-366
GRACEY
1/2

SIMZ-YM3-367
GRACEY
3/4

SIMZ-YM3-368
GRACEY
5/6

SIMZ-YM3-369
GRACEY
7/8

SIMZ-YM3-370
GRACEY
9/10

SIMZ-YM3-371
GRACEY
11/12

SIMZ-YM3-372
GRACEY
13/14


SIMZ-YM3-373

Gracey Curettes Set of 7 Pcs

- Gracey Curettes 1/2
- Gracey Curettes 3/4
- Gracey Curettes 5/6
- Gracey Curettes 7/8
- Gracey Curettes 9/10
- Gracey Curettes 11/12
- Gracey Curettes 13/14

CURETTES

- Кюрети
 - Curetas
 - Kuretten
- Кюрети
 - Curette
 - كرتيس




SIMZ-YM3-381
 Gracey Curette Set of 7 Pcs
 Gracey Curette 1-2
 Gracey Curette 3-4
 Gracey Curette 5-6
 Gracey Curette 7-8
 Gracey Curette 9-10
 Gracey Curette 11-12
 Gracey Curette 13-14

SIMZ-YM3-389
 Gracey Curette
 (Black & Gold Coating)
 Set of 7 Pcs
 Gracey Curette 1-2
 Gracey Curette 3-4
 Gracey Curette 5-6
 Gracey Curette 7-8
 Gracey Curette 9-10
 Gracey Curette 11-12
 Gracey Curette 13-14

CURETTES

- Кюрети
- Curetas
- Küretten

- Кюрети
- Curette
- کرس




Please Ask for Special Packing of Gracey Curettes Above

CURETTES

- Кюрети
- Curetas
- Küretten
- Кюрети
- Curette
- کورتس

Gracey Curettes Hollow Handle Ø 10 mm, Length 170mm



SIMZ-YM3-406
1/2

SIMZ-YM3-407
3/4

SIMZ-YM3-408
5/6

SIMZ-YM3-409
7/8

SIMZ-YM3-410
9/10

SIMZ-YM3-411
11/12

SIMZ-YM3-412
13/14

SIMZ-YM3-413
15/16

SIMZ-YM3-414
17/18


CURETTES

- Кюрети
- Curetas
- Küretten

- Кюрети
- Curette
- كرش



Gracey Curettes Solid Handle Ø 6 mm, Length 170mm



SIMZ-YM3-415
1/2

SIMZ-YM3-416
3/4

SIMZ-YM3-417
5/6

SIMZ-YM3-418
7/8

SIMZ-YM3-419
9/10

SIMZ-YM3-420
11/12

SIMZ-YM3-421
13/14

SIMZ-YM3-422
15/16


SIMZ-YM3-423
17/18

CURETTES

- Кюрети
- Curetas
- Kuretten
- Кюрети
- Curette
- كرتس

GRACEY MINI FIVE
HOLLOW HANDLE Ø 10 MM, LENGTH 170MM

GRACEY AFTER FIVE
HOLLOW HANDLE Ø 10 MM, LENGTH 170MM



SIMZ-YM3-424	SIMZ-YM3-425	SIMZ-YM3-426	SIMZ-YM3-427	SIMZ-YM3-428	SIMZ-YM3-429	SIMZ-YM3-430	SIMZ-YM3-431	SIMZ-YM3-432	SIMZ-YM3-433	SIMZ-YM3-434	SIMZ-YM3-435
Gracey 1/2	Gracey 3/4	Gracey 5/6	Gracey 7/8	Gracey 11/12	Gracey 13/14	Gracey 1/2 AF	Gracey 3/4 AF	Gracey 5/6 AF	Gracey 7/8 AF	Gracey 11/12 AF	Gracey 13/14 AF

GRACEY CURETTES

MICRO MINI




SIMZ-YM3-436
Micro Mini
Five 1/2

SIMZ-YM3-437
Micro Mini
Five 7/8

SIMZ-YM3-438
Micro Mini
Five 11/12

SIMZ-YM3-439
Micro Mini
Five 13/14

GRACEY PROPHY CURETTES



SIMZ-YM3-440
P1/2

SIMZ-YM3-441
P3/4

SIMZ-YM3-442
P5/6

SIMZ-YM3-443
P7/8

CURETTES

- Кюрети
- Curetas
- K retten
- Кюрети
- Curette
- كرتس



GRACEY 1/2



GRACEY 7/8



GRACEY 11/12



GRACEY 13/14



GRACEY 15/16



GRACEY 17/18



SIMZ-YM3-444
1/2 (shown)

SIMZ-YM3-446
3/4 (shown)

SIMZ-YM3-448
5/6 (shown)

SIMZ-YM3-451
7/8 (shown)

SIMZ-YM3-454
9/10 (shown)

SIMZ-YM3-456
11/12 (shown)

SIMZ-YM3-459
13/14 (shown)

SIMZ-YM3-462
15/16 (shown)

SIMZ-YM3-464
17/18 (shown)

SIMZ-YM3-445
1/2 Rigid

SIMZ-YM3-447
3/4 Rigid

SIMZ-YM3-449
5/6 Rigid

SIMZ-YM3-452
7/8 Rigid

SIMZ-YM3-455
9/10 Rigid

SIMZ-YM3-457
11/12 Rigid

SIMZ-YM3-460
13/14 Rigid

SIMZ-YM3-463
15/16 Rigid

SIMZ-YM3-465
17/18 Rigid

SIMZ-YM3-450
5/6 Extra Rigid

SIMZ-YM3-453
7/8 Extra Rigid

SIMZ-YM3-458
11/12 Extra Rigid


SIMZ-YM3-461
13/14 Extra Rigid

- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Инструменти пародонтологічні
- instrumentos para periodoncía
- الصكوك اللثة

PERIODONTAL INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm

SCALER CURETTE



SIMZ-YM3-466
11/12
Scaler

SIMZ-YM3-467
11A/12A
Scaler

SIMZ-YM3-468
11/12S
Modied
McCall

SIMZ-YM3-469
13S/14S
Pointed
Curette

SIMZ-YM3-470
MTC 13/14
Pointed
Curette

SIMZ-YM3-471
17S/18S
Curette


SIMZ-YM3-472
19/20
Curette

- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Інструменти пародонтологічні
- instrumentos para periodoncia
- الصوكك التتة

PERIODONTAL INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm


SCALER CURETTE



- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Инструменти пародонтологічні
- instrumentos para periodoncia
- الصكوك للتة

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm

**SIMZ-YM3-480
Columbia Set**



Columbia 2L/2R

Columbia 4L/4R


Columbia 13/14

SIMZ-YM3-481

SIMZ-YM3-482

SIMZ-YM3-483

**SIMZ-YM3-484
McCall Set**



McCalls 14S/18S

McCalls 17/18

McCalls 13S/14S

McCalls 13/14

SIMZ-YM3-485

SIMZ-YM3-486


SIMZ-YM3-487

SIMZ-YM3-488

- Кюрети
- Curetas
- Kuretten
- Кюрети
- Curette
- كورتس

SCALER / CURETTE SETS

Mc CALL SET



SIMZ-YM3-489
11/12
Scaler

SIMZ-YM3-490
11A/12A
Scaler

SIMZ-YM3-491
11/12S
Modified
McCall

SIMZ-YM3-492
13S/14S
Pointed
Curette

SIMZ-YM3-493
MTC 13/14
Pointed
Curette

SIMZ-YM3-494
17/18
Curette

SIMZ-YM3-495
17S/18PT
Curette


SIMZ-YM3-496
19/20
Curette

SIMZ-YM3-497
Goldman-Fox

CURETTES

- Кюрети
- Curetas
- Küretten

- Кюрети
- Curette
- کرش



SIMZ-YM3-498
204

SIMZ-YM3-499
204
IUFW

SIMZ-YM3-500
204S

SIMZ-YM3-501
204SD

SIMZ-YM3-502
1
Springfield
Anterior

SIMZ-YM3-503
2/3
Springfield


SIMZ-YM3-504
2/3
Springfield

SIMZ-YM3-505
J1/2
Johnson

SIMZ-YM3-506
Sickle-Scaler
Posterior

Sickle Scalers Instruments

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm




- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Инструменти пародонтологічні
- instrumentos para periodoncía
- المسكوك للتانة

SICKLE SCALERS INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm

SICKLE SCALERS



SIMZ-YM3-514
H6/H7
Hygienist

SIMZ-YM3-515
6/7
Goldman

SIMZ-YM3-516
U15
Towner

SIMZ-YM3-517
H5/U15
Hygienist-Towner

SIMZ-YM3-518
U15/30
Towner-Jacquette

SIMZ-YM3-519
H15/33†
Hygienist-Jacquette

SIMZ-YM3-520
U15/33
Towner-Jacquette

- Кюрети
- Curetas
- Küretten
- Кюрети
- Curette
- کرتس

Large Prophy Kit Root Planing Kit



5 - Ø 24 mm
Box of 12

SIMZ-YM3-521
Solid - Cone Socket
125 mm



SIMZ-YM3-522
2045D
175 mm



SIMZ-YM3-523
H6/H7
170 mm



SIMZ-YM3-524
COLUMBIA
4RL
175 mm



SIMZ-YM3-525
GOLDMAN-FOX
3
170 mm



SIMZ-YM3-526
16/17
185 mm



SIMZ-YM3-527
11/12 ODU
165 mm



SIMZ-YM3-528
MERIAM
160 mm


- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Инструменти пародонтологічні
- instrumentos para periodoncia
- المسكوك للتة

PERIODONTAL INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm

JACQUETTE SCALERS

CRANE-KAPLAN




- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Инструменти пародонтологічні
- instrumentos para periodoncia
- الصكوك اللثة

PERIODONTAL INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm


JACQUETTE SCALERS

CRANE-KAPLAN



SCALERS

JACQUETTE SCALERS



SIMZYM3-541 Nevi 1 Anterior	SIMZ-YM3-542 Nevi 1/H5 Anterior	SIMZ-YM3-543 Nevi 2 Posterior	SIMZ-YM3-544 Nevi 3 Posterior	SIMZ-YM3-545 Nevi 4 Posterior	SIMZ-YM3-546 30/33 Jacquette	SIMZ-YM3-547 31/32 Jacquette	SIMZ-YM3-548 34/35 Jacquette	SIMZ-YM3-549 H3/H4 Jacquette	SIMZ-YM3-550 14/15 Posterior
--	--	--	--	--	---	---	---	---	---

- Кюрети
- Curetas
- Kuretten
- Кюрети
- Curette
- کرتس

SICKLE SCALERS

POSTERIOR



SIMZ-YM3-551
135

SIMZ-YM3-552
REM 3/4
Remington



SIMZ-YM3-553
129

SIMZ-YM3-554
T2/T3
Taylor

ANTERIOR



SIMZ-YM3-555
137

SIMZ-YM3-556
0/00
Morse

CURETTES

- Кюрети
- Curetas
- Küretten

- Кюрети
- Curette
- كرش



SIMZ-YM3-557
137
Mini Five



SIMZ-YM3-558
OD1/33
Odont/
Jacquette



SIMZ-YM3-559
128
Nebraska



SIMZ-YM3-560
USC-128



SIMZ-YM3-561
WH-2
Whiteside




SIMZ-YM3-562
YG-15
Younger-
Good

CURETTES

- Кюрети
 - Curetas
 - Küretten
- Кюрети
 - Curette
 - كرتيس


BALL SCALERS

CRANE-KAPLAN



SIMZ-YM3-563 CI-1
 SIMZ-YM3-564 CI-2/3
 SIMZ-YM3-565 152 Cumine
 SIMZ-YM3-566 B-1
 SIMZ-YM3-567 B-4
 SIMZ-YM3-568 B-2/3
 SIMZ-YM3-569 B-5/6
 SIMZ-YM3-570 CK-6
 SIMZ-YM3-571 CK-6D
 SIMZ-YM3-572 333/334 Mini Kaplan

SCALER / CURETTE SETS
GOLDMAN-FOX SET



SIMZ-YM3-573
GF-1
Universal Scaler

SIMZ-YM3-574
GF-2
Curette

SIMZ-YM3-575
GF-3
Curette

SIMZ-YM3-576
GF-4
Curette

SIMZ-YM3-577
GF-5 Hoe
Buccal/Lingual

SIMZ-YM3-578
GF-6
Hoe

SIMZ-YM3-579
GF-21

SIMZ-YM3-580
GF-21S/33
Combination

CURETTES


- Кюрети
 - Curetas
 - Küretten
- Кюрети
 - Curette
 - كرتس

SCALER / CURETTE SETS

GOLDMAN-FOX SET

CATTONI SET

COLUMBIA



SIMZ-YM3-581
103/106

SIMZ-YM3-582
104/105
Hoe/Chisel

SIMZ-YM3-583
112/113

SIMZ-YM3-584
107/108

SIMZ-YM3-585
107/108S
Modified

SIMZ-YM3-586
2R/2L
Columbia

SIMZ-YM3-587
4R/4L
Columbia


SIMZ-YM3-588
13/14
Columbia
(shown)

SIMZ-YM3-589
13/14
Columbia
Rigid

SIMZ-YM3-590
13/14
Columbia
Curved Tip

UNIVERSAL CURETTES

BARNHART



SIMZ-YM3-591
IU13/14 Indiana
University
Anterior

SIMZ-YM3-592
IU17/18 Indiana
University
Posterior

SIMZ-YM3-593
1/2 Bunting
Curette-Chisel
combination

SIMZ-YM3-594
5/6 Bunting

SIMZ-YM3-595
YG7/8
Younger-Good

SIMZ-YM3-596
1/2 Barnhart
(shown)

SIMZ-YM3-597
1/2 Barnhart
Rigid

SIMZ-YM3-598
5/6 Barnhart
(shown)

SIMZ-YM3-599
5/6 Barnhart
Rigid

CURETTES

- Кюрети
- Curetas
- Küretten

- Кюрети
- Curette
- كرتس



SIMZ-YM3-600
1/2



SIMZ-YM3-601
3/4



SIMZ-YM3-602
5/6



SIMZ-YM3-603
7/8



SIMZ-YM3-604
9/10



SIMZ-YM3-605
11/12



SIMZ-YM3-606
13/14

Please Ask for Special Packing of Gracey Curettes Above

SIMZ-YM3-607

**Gracey Curettes Set of 7 Pcs
with Silicon Grip**

Gracey Curettes 1/2

Gracey Curettes 3/4

Gracey Curettes 5/6

Gracey Curettes 7/8

Gracey Curettes 9/10

Gracey Curettes 11/12

Gracey Curettes 13/14

CURETTES

- Кюрети
- Curetas
- Kuretten

- Кюрети
- Curette
- كرش



SIMZ-YM3-608
Filling Instruments 1



SIMZ-YM3-609
Filling Instruments 2



SIMZ-YM3-610
Filling Instruments 3



SIMZ-YM3-611
Filling Instruments 4



SIMZ-YM3-612
Filling Instruments 5



SIMZ-YM3-613
Lucas Curettes 1



SIMZ-YM3-614
Lucas Curettes 2



SIMZ-YM3-615
Lucas Curettes 3




SIMZ-YM3-616
Lucas Curettes 4




SIMZ-YM3-617
Lucas Curettes 5

SCALERS Single Ended Solid Handle Ø 6mm



Tartar Remover Instruments
 Solid Handle Ø 6mm




- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Инструменти пародонтологічні
- instrumentos para periodoncia
- الصكوك اللثة

Periodontal Files

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm

Sugarman File Set

Excellent to remove and to crush strong calculus deposits



SIMZ-YM3-633
9/10
Schluger
Curved File

SIMZ-YM3-634
11/12
Buck
Straight File

SIMZ-YM3-635
1S/2S
Sugarman File

SIMZ-YM3-636
3S/4S
Sugarman File

SIMZ-YM3-637
3-7
Hirschfeld

SIMZ-YM3-638
5-11
Hirschfeld


SIMZ-YM3-639
9-10
Hirschfeld

Periodontal Files

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm

Sugarman File Set

Excellent to remove and to crush strong calculus deposits



SIMZ-YM3-640
9/10
Schluger
Curved File

SIMZ-YM3-641
11/12
Buck
Straight File

SIMZ-YM3-642
1S/2S
Sugarman File

SIMZ-YM3-643
3S/4S
Sugarman File


SIMZ-YM3-644
3-7
Hirschfeld

SIMZ-YM3-645
5-11
Hirschfeld

SIMZ-YM3-646
9-10
Hirschfeld

PERIODONTAL KNIVES

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm



SIMZ-YM3-647
1
Reddy Thin

SIMZ-YM3-648
2
Reddy Thin

SIMZ-YM3-649
3/4
Buck

SIMZ-YM3-650
5/6
Buck

SIMZ-YM3-651
1/2
Solt

SIMZ-YM3-652
1/2
Sanders


SIMZ-YM3-653
USC Towner
19/20

SIMZ-YM3-654
CK-3
Crane Kaplan

SIMZ-YM3-655
Dr. Allen
End-Cutting
Intrasulcular Knife

PERIODONTAL KNIVES

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm



SIMZ-YM3-656
15K/16K
Kirkland

SIMZ-YM3-657
1/2
Orban

SIMZ-YM3-658
1/2 Orban
Knife

SIMZ-YM3-659
7
Goldman-Fox

SIMZ-YM3-660
8
Goldman-Fox


SIMZ-YM3-661
9
Goldman-Fox

SIMZ-YM3-662
11
Goldman-Fox

- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Інструменти пародонтологічні
- instrumentos para periodoncia
- لصكوك اللثة

PERIODONTAL CHISELS

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm



- Инструменты для костной ткани
- Limes a os
- Knochenfeilen

- Инструменти для кісткової тканини
- Limas para huesos
- ملفات العظام



SIMZ-YM3-671
Fig. 1
PoloKoff



SIMZ-YM3-672
Fig. 33



SIMZ-YM3-673
Fig. 45
Miller



SIMZ-YM3-674
Fig. 64
Miller



SIMZ-YM3-675
Fig. 10



SIMZ-YM3-676
Fig. 1
Wahl



SIMZ-YM3-677
Fig. 2
Wahl



SIMZ-YM3-678
Fig. 12
Howard



SIMZ-YM3-679
Fig. 12A




SIMZ-YM3-680
Fig. 12CA

- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Инструменти пародонтологічні
- instrumentos para periodoncia
- الصكوك اللثة

PERIODONTAL CHISELS

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm



SIMZ-YM3-681
1/2
Kramer-Nevins

SIMZ-YM3-682
36/37
Rhodes
Back-Action

SIMZ-YM3-683
44/45
Remington
Hoe

SIMZ-YM3-684
TG

SIMZ-YM3-685
TG-0
Modified TG with
Ochsendein cutouts


SIMZ-YM3-686
3/4
Solt

SIMZ-YM3-687
13K/TG

SIMZ-YM3-688
13K/13KL

BONE CURETTES

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm



SIMZ-YM3-689
Hemingway

SIMZ-YM3-690
Hemingway

SIMZ-YM3-691
Hemingway

SIMZ-YM3-692
Martin

SIMZ-YM3-693
Martin


SIMZ-YM3-694
Martin

CURETTES

- Кюрети
- Curetas
- Kuretten
- Кюрети
- Curette
- كرتس

LUCAS SPOON SHAPE

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm



SIMZ-YM3-695
84

SIMZ-YM3-696
85

SIMZ-YM3-697
86

SIMZ-YM3-698
87

SIMZ-YM3-699
88

SIMZ-YM3-700
85/86
Lucas
Surgical Curette


SIMZ-YM3-701
86/87
Lucas
Surgical Curette

SIMZ-YM3-702
Lucas

SIMZ-YM3-703
Lucas

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصمغ


Amalgam and Plastic Filling Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm



FILLING INSTRUMENTS

- Инструменты для реставрации
- Instrumenti pentru restaurații
- Instruments pour plombage
- Instrumentos para obturaciones
- Füllungsinstrumente (mit Kontrawinkel)
- ملء الصكوك

Amalgam and Plastic Filling Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm




FILLING INSTRUMENTS

• Инструменты для реставрации
 • Instruments pour plombage
 • Füllungsinstrumente (mit Kontrawinkel)

• Инструменти для реставрації
 • Instrumentos para obturaciones
 • ملء المسوك




Amalgam and Plastic Filling Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm



FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

Amalgam Instruments and Burnishers Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm




FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصنوك




Amalgam Instruments and Burnishers Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm



FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

Plastic Filling Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm




FILLING INSTRUMENTS

• Инструменты для реставрации
 • Instruments pour plombage
 • Füllungsinstrumente (mit Kontrawinkel)

• Инструменти для реставрації
 • Instrumentos para obturaciones
 • ملء الصلوك



Plastic Filling Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm




FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك


CARVERS INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm



CARVERS INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm



FILLING INSTRUMENTS

- Инструменты для реставрации
- Инструменти для реставрації
- Instruments pour plombage
- Instrumentos para obturaciones
- Füllungsinstrumente (mit Kontrawinkel)
- ملء الصكوك

CARVERS INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm

CLEOIDS

1,7 mm	2 mm	2,5 mm	2,1 mm	2,7 mm	2 mm
0	1 89/92	2 4/5	3 W-A-C-D	4 3-6	5 110
1,7 mm	2,4 mm	2 mm	2,6 mm	3 mm	2 mm
SIMZ-YM3-791	SIMZ-YM3-792	SIMZ-YM3-793	SIMZ-YM3-794	SIMZ-YM3-795	SIMZ-YM3-796

Nystrom

Nystrom 1	Nystrom 2	Nystrom 3
SIMZ-YM3-797	SIMZYM3-798	SIMZ-YM3-799


CHISELS

2	4	6
Wedelstaedt 1/2	Wedelstaedt 3/4	Wedelstaedt 5/6
SIMZ-YM3-800	SIMZYM3-801	SIMZYM3-802


ENAMEL INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm

MARGINTRIMMERS




HATCHES



FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm
CHISELS




FILLING INSTRUMENTS

• Инструменты для реставрации
 • Instruments pour plombage
 • Füllungsinstrumente (mit Kontrawinkel)

• Инструменти для реставрації
 • Instrumentos para obturaciones
 • ملء الصنوك

















Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm















EXCAVATORS


- Экскаваторы
- Excavateurs
- Excavatoren
- Excavatori
- Excavadores
- حفارات

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm

Oval Spoon			Spoon			
1.2 mm	1.5 mm	2 mm	1.5 mm	1.2 mm	1.5 mm	2 mm
						
Oval spoon 17W	Oval spoon 18W	Oval spoon 19W	14 Spoon	17 Spoon	18 Spoon	19 Spoon
						
SIMZ-YM3-834	SIMZ-YM3-835	SIMZ-YM3-836	SIMZ-YM3-837	SIMZ-YM3-838	SIMZ-YM3-839	SIMZ-YM3-840


DENTINS REMOVALS			Spoon		
1,25 mm	1,5 mm	2 mm	1,15 mm	1,5 mm	1 mm
					
31L Endodontic excavators	32L Endodontic excavators	33L Endodontic excavators	38/39 Spoon	63/64	65/66
					
SIMZ-YM3-841	SIMZ-YM3-842	SIMZ-YM3-843	SIMZ-YM3-844	SIMZ-YM3-845	SIMZ-YM3-846

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm



Spoon Excavators Round Shape

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm



- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء المسكوك

Composite Gold Coated Instruments



SIMZ-YM3-869
Markley 2



SIMZ-YM3-870
Westcott N.21



SIMZ-YM3-871
NYU/B2



SIMZ-YM3-872
Filling instruments



SIMZ-YM3-873
Ball Filling Instruments



SIMZ-YM3-874
Filling Instruments



SIMZ-YM3-875
Filling Instruments



SIMZ-YM3-876
Filling Instruments



SIMZ-YM3-877
Plastic Filling Instruments



SIMZ-YM3-878
Heidemann



SIMZ-YM3-879
Plastic Filling Instruments

FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

Composite Gold Tip Instruments



SIMZ-YM3-880
Markley 2



SIMZ-YM3-881
Westcott N.21



SIMZ-YM3-882
NYU/B2



SIMZ-YM3-883
Filling instruments



SIMZ-YM3-884
Ball Filling Instruments



SIMZ-YM3-885
Filling Instruments



SIMZ-YM3-886
Filling Instruments



SIMZ-YM3-887
Filling Instruments



SIMZ-YM3-888
Plastic Filling Instruments



SIMZ-YM3-889
Heidemann



SIMZ-YM3-890
Plastic Filling Instruments

- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Инструменти пародонтологічні
- instrumentos para periodoncia
- المسكوكات اللثة

Heidemann Spatula

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



Excavators

Excavators Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



SIMZ-YM3-900



SIMZ-YM3-901
Fig. 3



SIMZ-YM3-902
Fig. 4



SIMZ-YM3-903
Exc # 14



SIMZ-YM3-904
Exc # 17



SIMZ-YM3-905
Exc # 18



SIMZ-YM3-906
Exc # 19




SIMZ-YM3-907
Exc # 22

- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Инструменти пародонтологічні
- instrumentos para periodoncia
- الصكوك اللثة

Composite Instruments

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Инструменти пародонтологічні
- instrumentos para periodoncia
- الصكوك اللثة

Crane Kaplan

Available in Hollow Handle Ø 8 mm



SIMZ-YM3-919



SIMZ-YM3-920
Jaquette



SIMZ-YM3-921
Columbia 13/14



SIMZ-YM3-922
Kirkland 15/16



SIMZ-YM3-923
Goldman Fox



SIMZ-YM3-924
Gracey 11/12



SIMZ-YM3-925
Gracey 13/14



SIMZ-YM3-926
Columbia 2R/2L



SIMZ-YM3-927
Orban Knife 1/2

Lucas Curettes


Available in Hollow Handle Ø 8 mm




SIMZ-YM3-928
Ø 2 mm



SIMZ-YM3-929



SIMZ-YM3-930
Ø 2.5 mm



SIMZ-YM3-931
Ø 3 mm


- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Инструменти пародонтологічні
- Intrumentos para periodoncia
- الصكوك اللثة

Sinus Lift Instruments Blue Plasma Coating

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm



SIMZ-YM3-932
Fig. 1



SIMZ-YM3-933
Fig. 2



SIMZ-YM3-934
Fig. 3



SIMZ-YM3-935
Fig. 4



SIMZ-YM3-936
Fig. 5



SIMZ-YM3-937
Fig. 6



SIMZ-YM3-938
Fig. 7



SIMZ-YM3-939
Fig. 8

Gingival Retractors

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm



SIMZ-YM3-940




SIMZ-YM3-941



SIMZ-YM3-942

- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Інструменти пародонтологічні
- instrumentos para periodoncia
- المستوكه اللثة


Blue Plasma Coating & Operative instruments
Available in Hollow Handle Ø 8 mm



- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء المسكوك

P.K. THOMAS WAXING INSTRUMENTS

P.K. THOMAS WAX CARVER SET OF 5 PCS
Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm



SIMZ-YM3-956 PKT-1 SIMZ-YM3-957 PKT-2 SIMZ-YM3-958 PKT-3 SIMZ-YM3-959 PKT-4 SIMZ-YM3-960 PKT-5

P.K. THOMAS WAX CARVER SET OF 5 PCS
Aluminum Anodizing Colours Handle Ø 8 mm



SIMZ-YM3-962 PKT-1 SIMZ-YM3-963 PKT-2 SIMZ-YM3-964 PKT-3 SIMZ-YM3-965 PKT-4 SIMZ-YM3-966 PKT-5


SIMZ-YM3-967
P.k. Thomas Wax Carvers Set of 5 Pcs
with Sterilized Cassette



SIMZ-YM3-961
P.k. Thomas Wax Carvers Set of 5 Pcs

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

COMPOSITE INSTRUMENTS TITANIUM BLUE ENDODONTIC



SIMZ-YM3-974
Composite Instruments Blue
Titanium with Cassettes

CEMENT SPATULAS



- Шпатели для цемента
- Spatules à ciment
- Zementspatel
- Шпателі для цементу
- Espátulas para cemento
- ملاعق الأسمنت

Available in Hollow Handle Ø 8 mm, Ø 10 mm & Solid Handle Ø 6 mm



SIMZ-YM3-975
Fig. 1



SIMZ-YM3-976
Fig. 2



SIMZ-YM3-977
Fig. 3



SIMZ-YM3-978
Fig. 00



SIMZ-YM3-979
Fig. 0



SIMZ-YM3-980
Fig. 1



SIMZ-YM3-981
Fig. 2



SIMZ-YM3-982
Fig. 3



SIMZ-YM3-983
Fig. 4



SIMZ-YM3-984
Weston




SIMZ-YM3-985



SIMZ-YM3-986

- Шпатели для цемента
- Spatules à ciment
- Zementspatel
- Шпателі для цементу
- Espátulas para cemento
- ملاط الاسمنت

Wax and Modelling Carvers



SIMZ-YM3-987
Le Cron

SIMZ-YM3-988
Zahle

SIMZ-YM3-989
Zahle

SIMZ-YM3-990
Zahle

SIMZ-YM3-991
Gritman

SIMZ-YM3-992
Beale

SIMZ-YM3-993
Roach

SIMZ-YM3-994
Vehe

CARVERS INSTRUMENTS


- Инструменты для реставрации
- Instruments a modeler
- Modellierinstrumente

- Инструменти для реставрації
- Instrumentos para modelar
- أدوات نطقع



SIMZ-YM3-995

Graded Handle Hylin Carver



SIMZ-YM3-996
Hylin

SIMZ-YM3-997
Hylin-Le Cron

SIMZ-YM3-998
Hylin-Ceramica

SIMZ-YM3-999
Hylin

SIMZ-YM3-1000

SIMZ-YM3-1001
Le Cron

SIMZ-YM3-1002
Le Cron

SIMZ-YM3-1003
Le Cron

WAX KNIVES

- Ножи для воска
- Couteaux a cire
- Wachsmesser
- Ножі для воску
- Cuchillos para cera
- مسكين شمع

WAX KNIVES



SIMZ-YM3-1004
Fig. 1
Fahnenstock
125mm

SIMZ-YM3-1005
Fig. 2
Fahnenstock
170mm

SIMZ-YM3-1006
Fig. 1
Lessmann
125mm

SIMZ-YM3-1007
Fig. 2
Lessmann
170mm

SIMZ-YM3-1008
Fig. 1
with small spoon
125mm

SIMZ-YM3-1009
Fig. 2
with large spoon
170mm

- Ножи для воска
- Couteaux a cire
- Wachsmesser

- Ножі для воску
- Cuchillos para cera
- السكاكين الشمع

PLASTER KNIVES



SIMZ-YM3-1010
160mm

SIMZ-YM3-1011
160mm

SIMZ-YM3-1012
Buffalo

SIMZ-YM3-1013
Buffalo

SIMZ-YM3-1014
Plaster Knife

SIMZ-YM3-1015
170mm

WAX KNIVES

- Ножи для воска
- Couteaux a cire
- Wachsmesser
- Ножі для воску
- Cuchillos para cera
- السكاكين الشمع

WAX KNIVES



SIMZ-YM3-1016
Sachs



SIMZ-YM3-1017
130mm
Ermert



SIMZ-YM3-1018
Hook Knives



SIMZ-YM3-1019
125mm
Antkoviak



SIMZ-YM3-1020



Special Inox
SIMZ-YM3-1021 170mm
SIMZ-YM3-1022 200mm

- Шпатель для замешивания
- Spatules a alginate et a platre
- Alginat- und Gipsspatel
- Шпатель для замішування
- Espatulas para alginato y para yeso
- ملاعق الجبس الجوانات

WAX SPATULAS



SIMZ-YM3-1023
Fig. 1
18.5 cm

SIMZ-YM3-1024
Fig. 2
18 cm

SIMZ-YM3-1025
Wooden Handle
195mm curved

SIMZ-YM3-1026
Wooden Handle
200mm curved

SIMZ-YM3-1027
Wooden Handle
225mm curved

SIMZ-YM3-1028
Wooden Handle
195mm

SIMZ-YM3-1029
Wooden Handle
215mm

SIMZ-YM3-1030
Wooden Handle
205mm



Denture Flask Lower, Brass
SIMZ-YM3-1031
Close View



Denture Flask Upper, Brass
SIMZ-YM3-1032
Close View



Round Dental Flask, Brass
SIMZ-YM3-1033
Close View



Open View



Small Flask with Compressor, Brass
SIMZ-YM3-1034



Compressor Single Flask
SIMZ-YM3-1035



Compressor Double Flask
SIMZ-YM3-1036



Compressor Double Flask
SIMZ-YM3-1037



Open View



Dental Flask
4 Parts, Brass
SIMZ-YM3-1038
Close View



Open View



Dental Single
Flask, Brass
SIMZ-YM3-1039
Close View



Open View



Dental Flask,
4 Parts Silver
SIMZ-YM3-1040
Close View



Dyerline Flask,
Silver
SIMZ-YM3-1041
Close View



Open View



Round Dental
Flask, Silver
SIMZ-YM3-1042
Close View



Open View

ARTICULATORS

- Артикуляторы
- Artocilar
- Artikulations
- Артикулятори
- Articular
- ارنیو لائوس



SIMZ-YM3-1043
Articulator



SIMZ-YM3-1044
Articulator



SIMZ-YM3-1045
Articulator
S



SIMZ-YM3-1047
Articulator



SIMZ-YM3-1046
Articulator
M



SIMZ-YM3-1048
Articulator
L

ARTICULATORS

- Артикуляторы
- Artocilar
- Artikulations

- Артикулятори
- Articular
- آرٹیکولائٹورس



SIMZ-YM3-1049
Articulator



SIMZ-YM3-1050
Articulator



SIMZ-YM3-1051
Articulator



SIMZ-YM3-1052
Articulator



SIMZ-YM3-1053
Articulator



SIMZ-YM3-1054
Articulator

SIMZ-YM3-1055

Technical Datas:

Intercondylar distance	120 mm.
condylar inclination	fixed at 30°
Bennet angle	fixed at 7°
Incisal table inclination	10°
Distance between mounting plates	110 mm.
Weight	700 gr.
Opening	180°

features:

Lightweight anatomic articulator made of aluminum die-casting. It can be used following the geometrical laws of medium values (Bonwill triangle). Fixing according to centric. Automatic return to the hinge axis with springs. Perfect stability also in reversed position. Easy access and visibility to the lingual surfaces. interchangeable self centring aluminum mounting plates.



SIMZ-YM3-1056

Technical Datas:

Intercondylar distance	120 mm.
condylar inclination	fixed at 30°
Bennet angle	fixed at 7°
Incisal table inclination	flat
Distance between mounting plates	110 mm.
Weight	700 gr.
Opening	180°

Special Features:

Adjustable back by registration screw to work without incisal rod by keeping constant trajectories and inclinations during side and back/forward movements. Pure translation movement (key and instructions included). Anti-shock hinge rod in case articulator falls down. Bigger size lock screws for mounting plates.



Increased height and width

- (1) More space when using removable stumps systems.
- (2) More space for installation of mobile prosthesis with high chewing or big base.
- (3) More space when using magnetic plates



SIMZ-YM3-1057

Technical Datas:

Intercondylar distance	120 mm.
condylar inclination	fixed at 30°
Bennet angle	fixed at 7°
Incisal table inclination	10°
Distance between mounting plates	85 mm.
Weight	400 gr.
Opening (item 5000)	125°
Opening (item 5010)	180°



SIMZ-YM3-1058



Technical Datas:


Intercondylar distance	110 mm.
condylar inclination	fixed at 30°
Bennet angle	fixed at 7°
Incisal table inclination	flat
Distance between mounting plates	85 mm.
Weight	450 gr.
Opening	180°



SIMZ-YM3-1059

Anvil
Dental Laboratory
Instrument.





Maxiofacial Surgery Instruments

A Bard-Parker handle with a No. 15 blade is the usual scalpel used to make incisions on the gums where a flap can be raised to gain access to tooth.

Retractors

Tongue, cheek and flap retractors are needed to provide adequate access to the tooth and are legion in number.

Needle Holders and Sutures

Sutures are used to replace and secure the gum flap in its former position, to limit the size of the wound and therefore restrict subsequent contamination with food debris, and to assist in controlling bleeding. The usual suture size for intra-oral work is 4/0 and the material may be absorbable or non-absorbable Scissors
Dissecting scissors and suture-cutting scissors are usually used.

Artery Forceps

These are designed to hold soft tissue without damaging it. College tweezers are used to lift up sutures prior to removing them.

Periosteal Elevators

A number are available and are similarly used like a normal elevator but periosteal elevators have two heads – one at each end. One end is sharp-tipped to expand the tooth socket while the other has a blunt end to remove granulation tissues or raise a gum flap.

Chisels and Mallet

Chisels are more appropriate in the operating theatre and are particularly useful for disto-angular third molars and upper wisdom teeth extraction.

Bone Curettes

The Mitchells trimmer is probably the most valuable instrument in the category and is used to remove granulation tissue in the tooth socket.

Bone Ronguer and Bone File


Used to smoothen any sharp bone edges after tooth extraction is done.



SCALPEL HANDLES

- Ручки для скальпелей
 - Manches de bistouris
 - Skalpellgriffe
- Ручки для скальпелів
 - Mangos de bisturi
 - مشرط مقابض

Scalpel Handles



SIMZ-YM4-001
Fig. 3
125mm



SIMZ-YM4-002
Fig. 3
125mm



SIMZ-YM4-003
Fig. 7
165mm



SIMZ-YM4-004
Fig. 5
145mm



SIMZ-YM4-005
Fig. 5A
145mm



SIMZ-YM4-006
6 Way Adjustable
Angled
Scalpel Handle



SIMZ-YM4-007
Fig. 3L
210mm



SIMZ-YM4-008
Fig. 4L
218mm



SIMZ-YM4-009
Fig. 4
135mm

- Ножи для скальпеля и скальпеля
- Scalpel Blade & Scalpel Knives
- Scalpel Blade & Scalpel Messer

- Ножи скальпельного леза та скальпеля
- Cuchillos de Escalpele y Escalpele
- مشرط شفرة و مشرط المسكين

STERIL SCALPEL BLADES



SIMZ-YM4-010
10



SIMZ-YM4-011
11



SIMZ-YM4-012
20



SIMZ-YM4-013
21



SIMZ-YM4-014
22



SIMZ-YM4-015
12



SIMZ-YM4-016
15



SIMZ-YM4-017
23



SIMZ-YM4-018
24



SIMZ-YM4-019
165mm



Kit of scalpel blade
remover and orientable
blade scalpel handle

SCALPEL KNIVES



DIEFFENBACH
125mm



SIMZ-YM4-020
Fig. 1
125mm



SIMZ-YM4-021
Fig. 2
130mm



SIMZ-YM4-022
Fig. 4
135mm



SIMZ-YM4-023
Fig. 5
140mm



SIMZ-YM4-024
Fig. 7
125mm



SIMZ-YM4-025
Fig. 8
130mm



SIMZ-YM4-026
Fig. 9
135mm



SIMZ-YM4-027
Fig. 10
140mm

SCALPEL KNIVES

- Скальпельные ножи
- Couteaux à scalpel
- Skalpellmesser
- Ножи скальпеля
- Cuchillos Scalpel
- سكاكين مشرط



Dieffenbach
160mm



SIMZ-YM4-028
Fig. 1
160mm



SIMZ-YM4-029
Fig. 2
160mm



SIMZ-YM4-030
Fig. 3
160mm



SIMZ-YM4-031
Fig. 4
160mm



SIMZ-YM4-032
Fig. 5
160mm



SIMZ-YM4-033
Fig. 6
160mm



SIMZ-YM4-034
Fig. 7
160mm



SIMZ-YM4-035
Fig. 8
160mm



SIMZ-YM4-036
Fig. 9
160mm



SIMZ-YM4-037
Fig. 10
160mm



SIMZ-YM4-038
Fig. 15
160mm



SIMZ-YM4-039
Fig. 16
160mm

RETRACTORS - MOUTH GAGS

- Роторасширители
- Ecarteurs - Ouvre-bouches
- Vestibulumabhalter / Mundspanner
- Роторозширювачі
- Separadores - Abre-bocas
- الكاشطات، الكشحات الفم



SIMZ-YM4-040 105mm
SIMZ-YM4-041 135mm
Weitlaner
Fig. 1



SIMZ-YM4-042 105mm
SIMZ-YM4-043 135mm
Weitlaner
Fig. 2



SIMZ-YM4-044
Stainless Steel
Double
Ended Cheek
Retractors



SIMZ-YM4-045
Stainless Steel
Single
Ended Cheek
Retractors



SIMZ-YM4-046 110mm
SIMZ-YM4-047 130mm
SIMZ-YM4-048 150mm
Weitlander



SIMZ-YM4-049 110mm
SIMZ-YM4-050 130mm
SIMZ-YM4-051 150mm
Jennings



SIMZ-YM4-052
Denhard
120mm



SIMZ-YM4-053
Roser-Konig
160mm



SIMZ-YM4-054
Roser-Konig
190mm

RETRACTORS - MOUTH GAGS

- Роторасширители
- Ecarteurs - Ouvre-bouches
- Vestibulumabhalter / Mundspanner
- Роторозширювачі
- Separadores - Abre-bocas
- الكاشطات - الكاشطات الفم



SIMZ-YM4-055
HMG
Heister
130mm



SIMZ-YM4-056
Mouth Gag Molt
110mm for (Children)



SIMZ-YM4-057
Mouth Gag Molt
140mm for (Adults)



SIMZ-YM4-058
140mm

SIMZ-YM4-059
150mm

Cheek Retractor



SIMZ-YM4-060
Wangen Holder
150mm, with or without
lip protection



SIMZ-YM4-061
Wangen Holder
150mm, with
or without lip protection



SIMZ-YM4-062
Wangen Holder
Wangen holder with dual
function, 170mm



SIMZ-YM4-063
Retractor with
Ergo Grip
240mm
ergonomic handle



SIMZ-YM4-064
Cheek Retractor
Retractor Ivory
Width adjustable
& flexible, 200mm

- Ранорасширители
 - Ecarteurs
 - Wundhaken
- Ранорозширювачі
 - Separadores
 - الكاشحات

RETRACTORS



SIMZ-YM4-065
Columbia
14.5 cm



SIMZ-YM4-066
Sternberg
14 cm
Double Ended



SIMZ-YM4-067
Sternberg
16 cm



SIMZ-YM4-068
Sternberg
14 cm
Single Ended



SIMZ-YM4-069
Oringer Medium
11x8 cm



SIMZ-YM4-070
Brunings
Cheek Retractor
190mm



SIMZ-YM4-071
Minnesota
Small
140mm



SIMZ-YM4-072
Minnesota
Large
160mm



SIMZ-YM4-073
Cheek Retractor
200mm
Easy to Bend
by Hand



SIMZ-YM4-074
Cheek Retractor
Retractor and Tongue
Pushers Bishop with
Comfortable Margin,
160mm



SIMZ-YM4-075
Cheek Retractor
Retractor Combination
250mm

RETRACTORS

- Ранорасширители
 - Ecarteurs
 - Wundhaken
- Ранорозширювачі
 - Separadores
 - رانوراسشیرکن



SIMZ-YM4-076
Cheek Retractor
Long



SIMZ-YM4-077
Cheek Retractor
Long



SIMZ-YM4-078
Cheek Retractor
Long



SIMZ-YM4-079
Unique Retractor
Length 175mm
angled 45°



SIMZ-YM4-080
Unique Retractor
Length 175mm
angled 90°



SIMZ-YM4-081
Retractor
160mm



SIMZ-YM4-082
Retractor
190mm



SIMZ-YM4-083
Mini-Langenbeck
Retractor
160mm, Hooks 22 mm



SIMZ-YM4-084
Cheek Retractor
Retractors Combination
210mm, 2 Pcs



SIMZ-YM4-085
Mini Retractor
About 80mm



SIMZ-YM4-086
Senn Miller Retractor
16cm Hooks 22 mm



SIMZ-YM4-087
75x75
Austin Tissue
Retractors




SIMZ-YM4-088
75x75
Austin Tissue
Retractors

130mm
Available with
round tip on
request

80mm
Available with
round tip on
request

RETRACTORS

- Ранорасширители
 - Ecarteurs
 - Wundhaken
- Ранорозширювачі
 - Separadores
 - فككيات



SIMZ-YM4-098 30x11 mm
SIMZ-YM4-099 30x16 mm
 Langenbeck Retractor
 210mm, hook



Middeldorpf Retractor
SIMZ-YM4-100 17 cm
SIMZ-YM4-101 22 cm



SIMZ-YM4-102
 Farabeuf
 (120 X 12mm)
 Set of 2 Pcs



SIMZ-YM4-103
 Seldin Retractor



SIMZ-YM4-104
 Feldman
 Retractor
 140mm

GUM SCISSORS

- Ножницы десневые
- Ciseaux a gengives
- Zahnfleischscheren
- Ножиці для ясен
- Tijeras para encias
- مقص صمغ



SIMZ-YM4-105
Gum Scissors
120mm
Straight



SIMZ-YM4-106
Gum Scissors
120mm
Curved



SIMZ-YM4-107
Micro Iris
90mm
Straight



SIMZ-YM4-108
Micro Iris
90mm
Curved



SIMZ-YM4-109
Iris
115mm
Curved



SIMZ-YM4-110
Iris
110mm
Straight



SIMZ-YM4-111
Iris
110mm
Curved



SIMZ-YM4-112
Angular, Saw edge
115mm



SIMZ-YM4-113
Angular
115mm



Metzzenbaum
SIMZ-YM4-114 140mm Straight
SIMZ-YM4-115 200mm Straight
SIMZ-YM4-116 140mm Curved
SIMZ-YM4-117 200mm Curved

GUM SCISSORS

- Ножницы десневые
- Ciseaux à gencives
- Zahnfleischscheren
- Ножиці для ясен
- Tijeras para encías
- مقصس مساج



SIMZ-YM4-118
Goldman-Fox
130mm
Saw edge,
straight



SIMZ-YM4-119
Goldman-Fox
130mm
Saw edge,
curved



SIMZ-YM4-120
Goldman-Fox
130mm
Saw edge,
double curved



SIMZ-YM4-121
La Grange
115mm
Saw edge,
double curved



SIMZ-YM4-122
Neumann
120mm
Saw edge,
double curved



SIMZ-YM4-123
Kilner (Ragnell)
120mm
Straight



SIMZ-YM4-124
Kilner (Ragnell)
120mm
Curved



SIMZ-YM4-125
Joseph
140mm
Straight

GUM SCISSORS

- Ножницы десневые
- Ciseaux a gencives
- Zahnfleischscheren
- Ножиці для ясен
- Tijeras para encias
- مقص سمع



SIMZ-YM4-126
Joseph
140mm
Curved



SIMZ-YM4-127
Kelly
160mm
Straight



SIMZ-YM4-128
Kelly
160mm
Curved



SIMZ-YM4-129
Locklin
160mm
Saw edge,
angular



Locklin
SIMZ-YM4-130 120mm
SIMZ-YM4-131 160mm
Saw edge,
angular



SIMZ-YM4-132
Dean
170mm
Saw edge,
angular



Spencer
SIMZ-YM4-133 90mm
SIMZ-YM4-134 110mm



SIMZ-YM4-135
Spencer
130mm



SIMZ-YM4-136
Spencer
115mm



SIMZ-YM4-137
Straight
120mm

GUM SCISSORS

- Ножницы десневые
- Ciseaux a gengivas
- Zahnfleischscheren
- Ножницы для ясен
- Tijeras para encias
- مقص سمع



SIMZ-YM4-138
Straight
115mm



SIMZ-YM4-139
Curved
115mm



Curved
SIMZ-YM4-143 140mm
SIMZ-YM4-144 160mm
SIMZ-YM4-145 180mm

Standard Operating Scissors
Straight
SIMZ-YM4-140 140mm
SIMZ-YM4-141 160mm
SIMZ-YM4-142 180mm



Standard Operating Scissors
Sharp/Sharp Straight
SIMZ-YM4-146 140 mm
SIMZ-YM4-147 160 mm
SIMZ-YM4-148 180 mm



Standard Operating Scissors
Sharp/Sharp Curved
SIMZ-YM4-149 140 mm
SIMZ-YM4-150 160 mm
SIMZ-YM4-151 180 mm



Mayo
SIMZ-YM4-152 140mm
SIMZ-YM4-153 170mm
SIMZ-YM4-154 200mm
Straight
SIMZ-YM4-155 140mm
SIMZ-YM4-156 170mm
SIMZ-YM4-157 200mm
Curved

GUM SCISSORS

- Ножницы десневые
- Ciseaux a gengivas
- Zahnfleischscheren
- Ножницы для ясен
- Tijeras para encias
- مقص صمغ



SIMZ-YM4-158
Noyes
110mm, Straight



SIMZ-YM4-159
Noyes
110mm
Curved



SIMZ-YM4-160
Noyes
110mm
Angular



Castroviejo
SIMZ-YM4-161 150mm
SIMZ-YM4-162 180mm
Curved




SIMZ-YM4-163
Westcott Scissors
110mm



SIMZ-YM4-164
155mm

CROWN SCISSORS

- Ножницы для короны
- Ciseaux de couronne
- Krone Schere
- Корональні ножиці
- Tijeras de la corona
- مقص التاج



SIMZ-YM4-165
Universal
120mm
Saw edge



SIMZ-YM4-166
Beebee
110mm
Saw edge, straight



SIMZ-YM4-167
Beebee
110mm
Saw edge, curved



SIMZ-YM4-168
Beebee
120mm
Saw edge, straight



SIMZ-YM4-169
Beebee
120mm
Saw edge, curved



Lister Bandage Scissors
SIMZ-YM4-170 90mm
SIMZ-YM4-171 110mm
SIMZ-YM4-172 140mm
SIMZ-YM4-173 180mm
SIMZ-YM4-174 240mm



SIMZ-YM4-175
Dental scissors
large 175mm



SIMZ-YM4-176
180mm
universal
autoclavable up to
143 °C

GUM SCISSORS

- Ножницы десневые
- Ciseaux a gengives
- Zahnfleischscheren
- Ножиці для ясен
- Tijeras para encias
- مقص صمغ

T.C GOLD



SIMZ-YM4-177
Iris
110mm
Straight



SIMZ-YM4-178
Iris
110mm
Curved



SIMZ-YM4-179
Goldman-Fox
130mm
Saw edge, straight



SIMZ-YM4-180
Goldman-Fox
130mm
Saw edge, Curved



SIMZ-YM4-181
Iris Double Bend
11cm



SIMZ-YM4-182
La Grange
115mm
Saw edge,
double Curved



Mayo
SIMZ-YM4-183 140mm
SIMZ-YM4-184 170mm
SIMZ-YM4-185 230mm
Straight



Mayo
SIMZ-YM4-186 140mm
SIMZ-YM4-187 170mm
SIMZ-YM4-188 230mm
Curved

GUM SCISSORS

- Ножицы десневые
- Ciseaux a gengivas
- Zahnfleischscheren
- Ножиці для ясен
- Tijeras para encias
- مقص مساج



T.C GOLD



SIMZ-YM4-189
Dean Scissors TC
17cm curve



Metzenbaum
SIMZ-YM4-190 145mm
SIMZ-YM4-191 180mm
SIMZ-YM4-192 200mm
Straight



Metzenbaum
SIMZ-YM4-193 145mm
SIMZ-YM4-194 180mm
SIMZ-YM4-195 200mm
Curved



SIMZ-YM4-196
Kelly
160mm
Straight



SIMZ-YM4-197
Kelly
160mm
Curved



SIMZ-YM4-198
Universal
120mm



SIMZ-YM4-199
Beebee
110mm
Saw edge,
straight



SIMZ-YM4-200
Beebee
110mm
Saw edge,
curved



SIMZ-YM4-201
Beebee
120mm
Saw edge,
straight



SIMZ-YM4-202
Beebee
120mm
Saw edge,
curved

GUM SCISSORS

- Ножницы десневые
- Ciseaux a gencives
- Zahnfleischscheren
- Ножницы для ясен
- Tijeras para encias
- مقص صمغ



SIMZ-YM4-203
Goldman Fox Scissors
12.5cm Curved



SIMZ-YM4-204
Iris Scissors
11.5cm Curved
Serrated For Slip-Free Cutting



SIMZ-YM4-205
LaGrange Scissors
11.5cm Double Curved
Serrated For Slip-Free Cutting

Double
Curved



SIMZ-YM4-206
Castroviejo Micro Scissors
12cm Curved



SIMZ-YM4-207
Super Cut Goldman Fox
13.2cm Curved
Serrated For Slip-Free Cutting



SIMZ-YM4-208
Super Cut Iris
11cm Straight
Serrated For Slip-Free Cutting



SIMZ-YM4-209
Super Cut Iris
11cm Curved
Serrated For Slip-Free Cutting



SIMZ-YM4-210
Super Cut LaGrange
11.5cm
Double Curved
Serrated For Slip-Free Cutting

Double
Curved

- Ножницы десневые
- Ciseaux a gengivas
- Zahnfleischscheren
- Ножки для ясен
- Tijeras para encias
- مقص سمع

T.C. / CERAMIC COATED SCISSORS



Mayo Dissecting Scissors



SIMZ-YM4-211
14.5cm
SIMZ-YM4-212
17cm
SIMZ-YM4-213
23cm
B/B, Str.



SIMZ-YM4-214
14.5cm
SIMZ-YM4-215
17cm
SIMZ-YM4-216
23cm
B/B, Cvd.



Metzenbaum Standard Scissors



SIMZ-YM4-217
14.5cm
SIMZ-YM4-218
18cm
SIMZ-YM4-219
20cm
SIMZ-YM4-220
23cm
SIMZ-YM4-221
26cm
SIMZ-YM4-222
28cm
SIMZ-YM4-223
30cm
B/B, Str.



SIMZ-YM4-224
14.5cm
SIMZ-YM4-225
18cm
SIMZ-YM4-226
20cm
SIMZ-YM4-227
23cm
SIMZ-YM4-228
26cm
SIMZ-YM4-229
28cm
SIMZ-YM4-230
30cm
B/B, Cvd.



Locklin Dissecting Scissors, S/S

SIMZ-YM4-231
16cm
1 Blade
Serrated



Locklin Dissecting Scissors, S/S

SIMZ-YM4-232
16cm
1 Blade
Serrated



Dean Dissecting Scissors, S/S

SIMZ-YM4-233
17cm
1 Blade
Serrated

HAEMOSTATIC FORCEPS

- Зажимы артериальные
- Pincès hemostatiques
- Arterienklemmen
- Затискачі артеріальні
- Pinças hemostáticas
- ملقط هایموستاتیک



Hartmann (Baby Mosquito)
SIMZ-YM4-234 100mm straight
SIMZ-YM4-235 100mm curved



Halstead Mosquito
SIMZ-YM4-236 125mm straight
SIMZ-YM4-237 125mm curved



Halstead Mosquito 1x2
SIMZ-YM4-238 125mm straight
SIMZ-YM4-239 125mm Curved



Kelly
SIMZ-YM4-240 140mm straight
SIMZ-YM4-241 140mm curved



Rochester-Pean
SIMZ-YM4-242 140mm
SIMZ-YM4-243 160mm
SIMZ-YM4-244 180mm
SIMZ-YM4-245 200mm
SIMZ-YM4-246 140mm
SIMZ-YM4-247 160mm
SIMZ-YM4-248 180mm
SIMZ-YM4-249 200mm
straight curved



Rochester-Kocher 1x2
SIMZ-YM4-250 140mm
SIMZ-YM4-251 160mm
SIMZ-YM4-252 180mm
SIMZ-YM4-253 200mm
SIMZYM4-254 140mm
SIMZYM4-255 160mm
SIMZYM4-256 180mm
SIMZYM4-257 200mm
straight curved




Kelly-Rankin
SIMZ-YM4-258 160mm straight
SIMZ-YM4-259 160mm curved



SIMZ-YM4-260
Heiss
200mm curved

HAEMOSTATIC FORCEPS & TOWEL CLAMPS

- Зажимы артериальные, зажимы для салфеток
- Затискачі артеріальні, затискачі для серветок
- Pincas hemostatiques, Pincas a compreses
- Pinzas hemostaticas, Pincas para compresas
- Arterienklemmen, Tuschklemmen
- الملقط هايموستاتيک المشابك منسوجة



SIMZ-YM4-261
Backhaus
Not-perforating and
not-damaging 140mm



SIMZ-YM4-262
Backhaus
130mm



SIMZ-YM4-263
Backhaus
85mm



SIMZ-YM4-264
Jones
90mm



SIMZ-YM4-265
Gross Maier
Dressing Forceps
200mm straight



SIMZ-YM4-266
curved



SIMZ-YM4-267
Foerster
Dressing Forceps
200mm



SIMZ-YM4-268
Young
Tongue Forceps
170mm

Spare rubber




SIMZ-YM4-269
Collins
Tongue Forceps
170mm

STERILIZING & TISSUE CLAMPS

- Стерилизационные и тканевые зажимы
- Стерилизационні та тканинні затискачі
- Stérilisation et pinces pour tissus
- Esterilização e grampos de tecido
- تعقيم & المشابك الأنسجة



Micro Masquito
Fig. 1
10 cm




Micro Masquito 1x2
Fig. 1
10 cm



Cheate
SIMZ-YM4-274
270mm



Rogge
SIMZ-YM4-275 200mm
SIMZ-YM4-276 300mm



Lorna
SIMZ-YM4-277 100mm
-YM4-278 130mm



SIMZ-YM4-279
Allis-Baby 3:4
150mm



SIMZ-YM4-280
Allis 4:5
200mm



SIMZ-YM4-281
Allis 6:7
200mm

NEEDLE HOLDERS

- Иглодержатели
- Porte-aiguilles
- Nadelhalter
- Голкотримачі
- Porta-agujas
- اصحاب ابرة



T.C GOLD



NEEDLE HOLDERS

- Иглодержатели
- Porte-aiguilles
- Nadelhalter
- Голкотримачи
- Porta-agujas
- اصحاب ابرة


T.C GOLD



Mathieu-Kocher
Needle Holder




SIMZ-YM4-295 140mm SIMZ-YM4-296 170mm



SIMZ-YM4-297 140mm SIMZ-YM4-298 170mm



Mathieu
Needle Holder



SIMZ-YM4-299 140mm SIMZ-YM4-300 170mm



SIMZ-YM4-301
Olsen-Hegar Needle
Holder With Scissors
140mm



SIMZ-YM4-302
Gillies Needle
Holder With Scissors
160mm



Castroviejo
SIMZ-YM4-303 140mm straight
SIMZ-YM4-304 180mm straight
SIMZ-YM4-305 140mm curved
SIMZ-YM4-306 180mm curved



Castroviejo
SIMZ-YM4-307 140mm straight
SIMZ-YM4-308 180mm straight
SIMZ-YM4-309 140mm curved
SIMZ-YM4-310 180mm curved

NEEDLE HOLDERS

- Игподержатели
- Porte-aiguilles
- Nadelhalter
- Голкотримачи
- Porta-agujas
- اصحاب ابرة



SIMZ-YM4-311
Lichtenberg-
Ryder T.C
14 cm



SIMZ-YM4-312
Lichtenberg-
Ryder T.C
14 cm



SIMZ-YM4-313
Lichtenberg-
Ryder T.C
17 cm



SIMZ-YM4-314
Lichtenberg-
Ryder T.C
17 cm



SIMZ-YM4-315
Barraquer
18 cm



SIMZ-YM4-316
Barraquer
18 cm



SIMZ-YM4-317
Barraquer T.C
18 cm



SIMZ-YM4-318
Barraquer T.C
18 cm



SIMZ-YM4-319
Tönnis
18 cm



SIMZ-YM4-320
Castroviejo
13 cm
Straight



SIMZ-YM4-321
Kalt Needle Holders T.C
14 cm Straight



SIMZ-YM4-322
Peet Splinter Forceps
11.5 cm

SIMZ-YM4-323
La Grange
Curved




SIMZ-YM4-324
All-Purpose Pliers

NEEDLE HOLDERS

- Иглодержатели
- Голкотримачи
- Porte-aiguilles
- Porta-agujas
- Nadelhalter
- اصحاب ابرة



Castroviejo



SIMZ-YM4-325
130mm

SIMZ-YM4-326
150mm



SIMZ-YM4-327
180mm



SIMZ-YM4-328
Halsey
130mm



SIMZ-YM4-329
Baumgartner
130mm



SIMZ-YM4-330
Derf
120mm



SIMZ-YM4-331
Delicate
150mm



SIMZ-YM4-332
Adson
170mm



Mayo-Hegar
SIMZ-YM4-333 130mm
SIMZ-YM4-334 160mm
SIMZ-YM4-335 180mm
SIMZ-YM4-336 200mm



SIMZ-YM4-337
Crile-Wood
150mm

NEEDLE HOLDERS

- Игподдержатели
- Porte-aiguilles
- Nadelhalter
- Голкотримачи
- Porta-agujas
- اصحاب ابرة



SIMZ-YM4-338
Olsen-Hegar Needle
Holder With Scissors
140mm



SIMZ-YM4-339
Gillies Needle
Holder With Scissors
160mm



Mathieu



SIMZ-YM4-340
140mm



SIMZ-YM4-341
170mm



SIMZ-YM4-342
200mm



Mathieu-Kocher



SIMZ-YM4-343
140mm



SIMZ-YM4-344
170mm





SIMZ-YM4-345
Boynton
120mm



RUBBER DAM CLAMPS FORCEPS

- Инструменты для системы коффердам
- Pincas porte-crampons a digue
- Kofferdam Klammerzangen
- Инструменты для системи коффердам
- Pinzas de engrapar para dique de caucho
- السد المطاط الملقط المشابيك



RUBBER DAM CLAMPS KITS

- Резиновая дамба зажимы комплекты
- Гумова дамба затискачи комплекти
- Trousses de caoutchouc
- Conjuntos de abrazaderas de retención de caucho
- Gummidammklammern
- المطاط المد المشترك مجموعات



Rubber Dam Kits




SIMZ-YM4-357

On request the Content of the rubber dam kit can be changed with other instruments.

Content:

- 01 Anisworth
- 02 Brewer
- 03 Clamps Stand
- 04 Rubber Dam Frame L
- 05 Rubber Dam Frame M
- 06 Rubber Dam Box Lower Case
- 07 Rubber Dam Box Upper Case

- 08 Rubber Dam Clamp Fig. 212
- 09 Rubber Dam Clamp Fig. 211
- 10 Rubber Dam Clamp Fig. 8
- 11 Rubber Dam Clamp Fig. 8A
- 12 Rubber Dam Clamp Fig. 14
- 13 Rubber Dam Clamp Fig. 14A
- 14 Rubber Dam Clamp Fig. 0
- 15 Rubber Dam Clamp Fig. 00
- 16 Rubber Dam Clamp Fig. 2
- 17 Rubber Dam Clamp Fig. 2A
- 18 Rubber Dam Clamp Fig. 7




SIMZ-YM4-358

On request the Content of the rubber dam kit can be changed with other instruments.

Content:

- 01 Ivory
- 02 Ivory light weight
- 03 Clamps Stand
- 04 Rubber Dam Frame M
- 05 Rubber Dam Frame L
- 06 Rubber Dam Box Lower Case
- 07 Rubber Dam Box Upper Case

- 08 Rubber Dam Clamp Fig. 209
- 09 Rubber Dam Clamp Fig. 212
- 10 Rubber Dam Clamp Fig. 211
- 11 Rubber Dam Clamp Fig. 8
- 12 Rubber Dam Clamp Fig. 8A
- 13 Rubber Dam Clamp Fig. 14
- 14 Rubber Dam Clamp Fig. 14A
- 15 Rubber Dam Clamp Fig. 0
- 16 Rubber Dam Clamp Fig. 00
- 17 Rubber Dam Clamp Fig. 2
- 18 Rubber Dam Clamp Fig. 2A
- 19 Rubber Dam Clamp Fig. 7



SIMZ-YM4-359
Rubber Dum Template

RUBBER DAM CLAMPS KITS

- Резиновая дамба зажимы комплекты
- Trousses de caoutchouc
- Gummidammklammern
- Гумова дамба затискачі комплекти
- Conjuntos de abrazaderas de retención de caucho
- المطاط الدم المشبك مجموعات



SIMZ-YM4-360

Stainless Steel Gold Plated Rubber Dam Clamp Tray with 11 Different Clamps
Fig. 212, 211, 8, 8A, 14, 14A, 0, 00, 2, 2A, 7



SIMZ-YM4-361

Stainless Steel Rubber Dam Clamp Tray with 12 Different Clamps
Fig. 209, 212, 211, 8, 8A, 14, 14A, 0, 00, 2, 2A, 7



SIMZ-YM4-362

Stainless Steel Rubber Dam Clamp Tray with 8 Different Clamps
Fig. 1, 26, 29, 51, 60, 202, 205, 210




SIMZ-YM4-363

Stainless Steel Rubber Dam Clamp Tray
(without clamps)

RUBBER DAM CLAMPS

- Кламмера для коффердама
- Grapas para dique de caucho
- المشابك سد المطاط
- Crampons pour digue
- Kofferdam Klammern



CROWN REMOVER

- Инструменты для удаления коронок
- Arrache-couronnes
- Teleskopkronenzangen
- Инструменты для видалення коронок
- Arranca-corona
- مزيلات التاج

Complete With 3 Attachment



SIMZ-YM4-422
Wynman Crown Gripper
No Lock T.C 16 cm



SIMZ-YM4-423
Crown Remover



A B



C



SIMZ-YM4-424
Crown Spreader, Crown Remover
Special Instrument for crown spreading after separating, simple and safe operation.



SIMZ-YM4-425
Crown Remover



SIMZ-YM4-426
Fig. 7



SIMZ-YM4-427
Fig. 8



SIMZ-YM4-428
Fig. 9



SIMZ-YM4-429
Fig. 10

For Molars and Premolars



SIMZ-YM4-430
Fig. 11



SIMZ-YM4-431
Fig. 12



SIMZ-YM4-432
Fig. 13
Toffemire

CROWN REMOVER

- Инструменты для удаления коронок
- Arrache-couronnes
- Teleskopkronenzangen

- Инструменты для выдалення коронок
- Arranca-corona
- مزيلات التاج



SIMZ-YM4-433
Crown Spreaders



SIMZ-YM4-434
Crown Spreaders
with Spring



SIMZ-YM4-435
Crown Spreader



SIMZ-YM4-436
Fig. 1




SIMZ-YM4-437
Fig. 2



SIMZ-YM4-438
Crown Spreaders
Fig. 3




SIMZ-YM4-439
Crown Remover



SIMZ-YM4-440
Crown Remover



SIMZ-YM4-441
Crown Remover



CROWN REMOVER

- Инструменты для удаления коронок
- Arrache-couronnes
- Teleskopkronenzangen
- Инструменти для видалення коронок
- Arranca-corona
- مزبلات التاج



SIMZ-YM4-442
Crown Remover
Turkish Style



SIMZ-YM4-443
Crown Remover
Turkish Style



SIMZ-YM4-444
Crown Remover
Turkish Style



SIMZ-YM4-445
Crown Remover
Turkish Style

MISCELLANEOUS PRODUCTS

- Прочие товары
- Divers produits et services
- Sonstige Produkte

- Інші товари
- Varios Productos
- المنتجات المتنوعة



SIMZ-YM4-446
X Ray Film Hangers



SIMZ-YM4-447
X Ray Film Hangers



SIMZ-YM4-448
X Ray Film Hangers



SIMZ-YM4-449
X Ray Film Hangers



SIMZ-YM4-450
X Ray Film Hangers



SIMZ-YM4-451
X Ray Film Hangers



SIMZ-YM4-452
X Ray Film Hangers



SIMZ-YM4-453
X Ray Film Hanger



SIMZ-YM4-454
X Ray Film Hanger



SIMZ-YM4-455
X Ray Film Hanger



SIMZ-YM4-456
Napkin holder with metal chain



SIMZ-YM4-457
Napkin Holder with
Plastic Chain



Mouth Props
SIMZ-YM4-458
Small



Mouth Props
SIMZ-YM4-459
Medium



Mouth Props
SIMZ-YM4-460
Large

- Матрицедержатели
- Porte Matrices
- Matrizenspanner
- Матрицеутримувачи
- Porta Matrices
- مصفوفة الخدم



SIMZ-YM4-461
Fig. 5
Ivory 8
5 mm



SIMZ-YM4-462
Fig. 6
Ivory 8
6 mm



SIMZ-YM4-463
Fig. 7
Ivory 8
7 mm



SIMZ-YM4-464
Ivory 8 N



SIMZ-YM4-465
Didner 8 P



SIMZ-YM4-466
Fig. 1
Tofflemire-Universal



SIMZ-YM4-467
Fig. 2
Tofflemire-Junier



SIMZ-YM4-468
Fig. 3
Tofflemire-Senior

Amalgam Well



SIMZ-YM4-469
Stainless Steel
Amalgam Well
Medium



SIMZ-YM4-470
Stainless Steel
Amalgam Well
large

MATRIX RETAINERS

- Матрицедержатели
- Porte Matrices
- Matrizenspanner
- Матрицеутримвачи
- Porta Matrices
- مصالفة الخدم



SIMZ-YM4-471
Tofflemire-Junior
45 mm



SIMZ-YM4-472
Tofflemire-Universal
60 mm



SIMZ-YM4-473
Matrix Ivory 8
6 mm



SIMZ-YM4-474
Didner 8 P
70 mm



SIMZ-YM4-475
Nyström
Left



SIMZ-YM4-476
Nyström
Right



SIMZ-YM4-477
Narrow



SIMZ-YM4-478
Wide



SIMZ-YM4-479
Ivory
with shields



SIMZ-YM4-480
Ivory
without shields



SIMZ-YM4-481
Ivory



SIMZ-YM4-482
Elliot

AMALGAM INSTRUMENTS

- Инструменты для амальгамы
- Instruments pour amalgame
- Amalgamisinstrumente
- Инструменты для амальгами
- Instrumentos para amalgama
- الصوكات الملقم

Ø 2mm



SIMZ-YM4-483
Fig. 10
165mm
English pattern

Ø 2mm



SIMZ-YM4-484
Fig. 1
185mm
German pattern

Ø 2mm



SIMZ-YM4-485
Fig. 1
150mm
Swiss pattern

Ø 2mm



SIMZ-YM4-486
Fig. 1
150mm

Ø 2mm



SIMZ-YM4-487
Fig. 2
150mm, 45°

Ø 2mm



SIMZ-YM4-488
Fig. 3
140mm, 90°



SIMZ-YM4-489
Chip Blower



SIMZ-YM4-490



SIMZ-YM4-491



SIMZ-YM4-492



SIMZ-YM4-493
Water Syringe

- Инструменты для амальгамы
- Instruments pour amalgame
- Amalgamisinstrumente
- Инструменти для амальгами
- Instrumentos para amalgama
- الصكوك المثلج

AMALGAM CARRIERS

TIP SIZES FOR DISTAL OR SINGLE-END

1.5 mm



SIMZ-YM4-494
Mini/
Single-End

2.0 mm



SIMZ-YM4-495
Regular
Single-End

2.8 mm



SIMZ-YM4-496
Large
Single-End

1.5 mm



SIMZ-YM4-497
Mini
Distal

2.0 mm



SIMZ-YM4-498
Regular
Distal

TIP SIZES FOR DOUBLE-END

1.5 mm



SIMZ-YM4-499
Mini/
Regular

2.8 mm



SIMZ-YM4-500
Large/
Regular

3.2 mm



SIMZ-YM4-501
Jumbo/
Regular



SIMZ-YM4-502
Amalgam
Well



Single-End

TIP SIZES FOR DOUBLE-END


SIMZ-YM4-503 - 1.5mm-2.0mm **SIMZ-YM4-506** - 1.0mm-2mm
SIMZ-YM4-504 - 2.0mm-2.8mm **SIMZ-YM4-507** - 2.0mm-3mm
SIMZ-YM4-505 - 2.0mm-3.2mm **SIMZ-YM4-508** - 2.0mm-4mm

SIMZ-YM4-509 - 1.5mm-2.0mm **SIMZ-YM4-512** - 1.0mm-2mm
SIMZ-YM4-510 - 2.0mm-2.8mm **SIMZ-YM4-513** - 2.0mm-3mm
SIMZ-YM4-511 - 2.0mm-3.2mm **SIMZ-YM4-514** - 2.0mm-4mm

SIMZ-YM4-515 - 1.5mm-2.0mm **SIMZ-YM4-518** - 1.0mm-2mm
SIMZ-YM4-516 - 2.0mm-2.8mm **SIMZ-YM4-519** - 2.0mm-3mm
SIMZ-YM4-517 - 2.0mm-3.2mm **SIMZ-YM4-520** - 2.0mm-4mm

SIMZ-YM4-521 - 1.5mm-2.0mm **SIMZ-YM4-524** - 1.0mm-2mm
SIMZ-YM4-522 - 2.0mm-2.8mm **SIMZ-YM4-525** - 2.0mm-3mm
SIMZ-YM4-523 - 2.0mm-3.2mm **SIMZ-YM4-526** - 2.0mm-4mm

SIMZ-YM4-527 - 1.5mm
SIMZ-YM4-528 - 2.0mm
SIMZ-YM4-529 - 2.0mm
SIMZ-YM4-530 - 3.0mm



ROOT PLANNING CHISELES & BONE FILES



SIMZ-YM4-531
Freer
160mm



SIMZ-YM4-532
Freer
160mm



SIMZ-YM4-533
Kramer-Nevins
Bone Chisel



SIMZ-YM4-534
Kramer-Nevins
Bone Chisel



SIMZ-YM4-535
Miller



SIMZ-YM4-536
Miller

BONE FILES



SIMZ-YM4-537
Miller 21



SIMZ-YM4-538
Miller 1X
cross cut



SIMZ-YM4-539
Miller 2X
cross cut



SIMZ-YM4-540
Miller 52



SIMZ-YM4-541
Wahl 1



Williger
Curette
13.5 cm



CHISELS



SIMZ-YM4-545



SIMZ-YM4-546



SIMZ-YM4-547



SIMZ-YM4-548

Bone Splitting Chisels

2.0mm



SIMZ-YM4-549

2.5mm



SIMZ-YM4-550

3.0mm



SIMZ-YM4-551

3.5mm



SIMZ-YM4-552

4.0mm



SIMZ-YM4-553

4.5mm



SIMZ-YM4-554

5.0mm



SIMZ-YM4-555

PERIODONTAL INSTRUMENTS

- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Інструменти пародонтологічні
- Instrumentos para periodoncia
- المسكوك الثلاثة




SIMZ-YM4-556 1 Ochsenbin
 SIMZ-YM4-557 2 Ochsenbin
 SIMZ-YM4-558 Wakefield 1W 170mm
 SIMZ-YM4-559 Wakefield 2W 170mm



ø 28 mm
 Teflon interchangeable inserts
 Autoclavable.

SIMZ-YM4-560 Mead 190 mm



ø 26 mm
 With interchangeable inserts
 320 gr – 11 ¼ oz

SIMZ-YM4-561 Mead 170 mm



ø 22 mm
 Lead filled
 200 gr – 7 oz


SIMZ-YM4-562 Partson 180mm




SIMZ-YM4-563 Vickers 19 cm



SIMZ-YM4-564 Herbert 19 cm



Partsch 140 mm
 SIMZ-YM4-565 3mm
 SIMZ-YM4-566 4mm
 SIMZ-YM4-567 5mm
 SIMZ-YM4-568 6mm



Partsch 170 mm
 SIMZ-YM4-569 3mm
 SIMZ-YM4-570 4mm
 SIMZ-YM4-571 5mm
 SIMZ-YM4-572 6mm



Partsch 140 mm
 SIMZ-YM4-573 3mm
 SIMZ-YM4-574 4mm
 SIMZ-YM4-575 5mm
 SIMZ-YM4-576 6mm



Partsch 170 mm
 SIMZ-YM4-577 3mm
 SIMZ-YM4-578 4mm
 SIMZ-YM4-579 5mm
 SIMZ-YM4-580 6mm

PERIOSTEEL ELEVATORS



SIMZ-YM4-581
Buser



SIMZ-YM4-582
Ohl 10



SIMZ-YM4-583
Howard 9H



SIMZ-YM4-584
Allen 9A



SIMZ-YM4-585
Prichard 3



SIMZ-YM4-586
Seldin 23



SIMZ-YM4-587



SIMZ-YM4-588
Raspator, Raspatory
Length 175mm



SIMZ-YM4-589



SIMZ-YM4-590




SIMZ-YM4-591
Raspator, Williger
Length 170mm, different widths
(4mm, 5mm or 7mm)



SIMZ-YM4-592



SIMZ-YM4-593
Trelat
140mm



- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Інструменти пародонтологічні
- Instrumentos para periodoncia
- الصكوك اللثة



SIMZ-YM4-594
Freer



SIMZ-YM4-595
Freer 15



SIMZ-YM4-596
Freer 16



SIMZ-YM4-597
Goldman-Fox
14



SIMZ-YM4-598
P 24G



SIMZ-YM4-599
Molt
175mm



SIMZ-YM4-600
Molt
175mm



SIMZ-YM4-601
Meat
175mm



SIMZ-YM4-602
Meat
175mm



SIMZ-YM4-603
Meat
175mm

- Шприцы стоматологические карпульные
- Seringues à usage dentaire
- Dentalspritzen

- Шприці стоматологічні карпульні
- Jeringas dentales
- المعالین الأسنان



Type A
SIMZ-YM4-604 1.8ml
SIMZ-YM4-605 2.2ml



Type C
SIMZ-YM4-606 1.8ml
SIMZ-YM4-607 2.2ml



Type C
SIMZ-YM4-608 1.8ml
SIMZ-YM4-609 2.2ml



Type Euro
SIMZ-YM4-610 1.8ml
SIMZ-YM4-611 2.2ml



3-Ring
SIMZ-YM4-612 1.8ml
SIMZ-YM4-613 2.2ml



2-Ring
SIMZ-YM4-614 1.8ml
SIMZ-YM4-615 2.2ml



Euro
Breech loading
SIMZ-YM4-616 1.8ml
SIMZ-YM4-617 2.2ml



Euro
Breech loading
SIMZ-YM4-618 1.8ml
SIMZ-YM4-619 2.2ml



Manual Aspirating Syringe
SIMZ-YM4-620 1.8ml
SIMZ-YM4-621 2.2ml



Quatro
SIMZ-YM4-622 1.8ml
SIMZ-YM4-623 2.2ml



Trio
SIMZ-YM4-624 1.8ml
SIMZ-YM4-625 2.2ml

DENTAL SYRINGES

- Шприцы стоматологические карпульные
- Шприці стоматологічні карпульні
- Seringues à usage dentaire
- Jeringas dentales
- المحاقن الأسنان
- Dentalspritzen



SIMZ-YM4-626 1.8ml.
SIMZ-YM4-627 2.2ml.
 Pistol-Ject
 injects 0.20ml with every click
 Lever action for high injection pressure



Autoclavable



SIMZ-YM4-636 1.8ml
SIMZ-YM4-637 2.2ml
 Fig. 1001
 Intraligamental
 Anesthesia Syringe in
 Black Colour Auto
 Cleavable



Autoclavable



Ject Ultima
SIMZ-YM4-628 1.8ml
SIMZ-YM4-629 2.2ml
 Injects 0.05 ml with every click
 Supplied with angled and straight tubes



Swedish Type
SIMZ-YM4-630 1.8ml
SIMZ-YM4-631 2.2ml



Ject angled
SIMZ-YM4-632 1.8ml
SIMZ-YM4-633 2.2ml



Ject straight
SIMZ-YM4-634 1.8ml
SIMZ-YM4-635 2.2ml



Disc System
SIMZ-YM4-638 1.8ml
SIMZ-YM4-639 2.2ml



Disc System
SIMZ-YM4-640 1.8ml
SIMZ-YM4-641 2.2ml



Disc System
SIMZ-YM4-642 1.8ml
SIMZ-YM4-643 2.2ml



Projection System
SIMZ-YM4-644 1.8ml
SIMZ-YM4-645 2.2ml



Projection System
SIMZ-YM4-646 1.8ml
SIMZ-YM4-647 2.2ml

- Шприцы стоматологические карпульные
- Шприці стоматологічні карпульні
- Seringues à usage dentaire
- Jeringas dentales
- Dentalspritzen
- المحاقن الأسنان



Stainless Steel - Dark
SIMZ-YM4-648 1.8 ml
SIMZ-YM4-649 2.2 ml



Stainless Steel- Silver
SIMZ-YM4-650 1.8 ml
SIMZ-YM4-651 2.2 ml



Stainless Steel- Orange
SIMZ-YM4-652 1.8 ml
SIMZ-YM4-653 2.2 ml



Stainless Steel- Green
SIMZ-YM4-654 1.8 ml
SIMZ-YM4-655 2.2 ml



Stainless Steel- Blue
SIMZ-YM4-656 1.8 ml
SIMZ-YM4-657 2.2 ml

BONE ROUNGER

- Костные кусачки
- Кусачки для кісток
- Hohlmeißelzangen
- Pinzas-gubias
- Pincès-gouges
- روتجورس العظام



SIMZ-YM4-658
Blumenthal
155 mm



SIMZ-YM4-659
Blumenthal 45°
150 mm



SIMZ-YM4-660
Blumenthal 90°
140 mm



SIMZ-YM4-661
Micro Luer
150 mm



SIMZ-YM4-662
Luer
150 mm



3 x 3SIMZ-YM4-663

4 x 4SIMZ-YM4-664

5 x 5SIMZ-YM4-665

Kerrison

BONE ROUNGER

- Костные кусачки
- Knochelmeißelzangen
- Pincès-gouges
- Кусачки для кісток
- Pinças-gubias
- رولنجورس العظام



SIMZ-YM4-666
Friedman
140 mm

SIMZ-YM4-667
Cleveland 4
140 mm

SIMZ-YM4-668
Micro Friedman
145 mm

SIMZ-YM4-669
Micro Friedman 90°
145 mm



SIMZ-YM4-670
Lempert Delicate
160 mm

SIMZ-YM4-671
Cleveland
155 mm

SIMZ-YM4-672
Beyer
170 mm

BONE ROUNGER

- Костные кусачки
- Hohlmeißelzangen
- Pincés-gouges
- Кусачки для кісток
- Pinzas-gubias
- روتجويرس العظام



SIMZ-YM4-673
Fig. 40
Ruskin



SIMZ-YM4-674
Fig. 41
Beyer



SIMZ-YM4-675
Fig. 50
Liston
140mm



SIMZ-YM4-676
Fig. 51
Liston
190mm



Cohen Tissue Nipper
SIMZ-YM4-677 100mm
SIMZ-YM4-678 140mm



SIMZ-YM4-679
Punch



SIMZ-YM4-680
Tomati
120mm

- Инструменты для удаления коронок
- Arrache-couronnes
- Teleskopkronenzangen

- Инструменты для удаления коронок
- Arranca-corona
- مزيلات التاج



SIMZ-YM4-681
160 mm

Crown Plier with plastic inserts Straight
For ceramical and plastic crowns and bridges, with plastic inserts (sterilization up to 180° C).



SIMZ-YM4-682
145 mm

Crown Plier with plastic inserts Angled



SIMZ-YM4-683
140 mm

Universal plier for the removal of ceramic crowns and bridges



SIMZ-YM4-684
140 mm

Universal plier for the removal of ceramic crowns and bridges



SIMZ-YM4-685
Furrer
Copper Ring
Removing Plier.
150 mm
Curved



SIMZ
SIMZ
SIMZ-YM4-686 125 mm
SIMZ-YM4-687 150 mm
Böhm
Pliers for telescope
crowns.
Diamond points
coated.



SIMZ-YM4-688
Metal tip
220 mm
Hammer With Wooden Handle



SIMZ-YM4-689
Metal tip
220 mm



SIMZ-YM4-690
1
For plaster



SIMZ-YM4-691
2
For metal



SIMZ-YM4-692
Fret saw
250 mm



SIMZ-YM4-693
Small
Saw Blades



SIMZ-YM4-694
Large
Saw Blades



SIMZ-YM4-695
Ceri Saw



SIMZ-YM4-696
Plaster Shears
200 mm



Crucible Tongs
SIMZ-YM4-697 250 mm
SIMZ-YM4-698 320 mm



SIMZ-YM4-699
Ceri Saw
180 mm




SIMZ-YM4-700
Blade
0,12x5 mm
Conf. 12 pcs



SIMZ-YM4-701
Blade
0,25x5 mm
Conf. 12 pcs



SIMZ-YM4-702
Blade
0,25x3 mm
Conf. 12 pcs



orthodontics pliers in Clinical Dentistry Orthodontics and Dentofacial Orthopaedics

Driven by the demand for tooth alignment, general practitioners are increasingly being asked to provide orthodontic treatment.


Dentists keen to participate in this area of expertise need to demonstrate the skills to identify suitable patients, the ability to provide conventional and contemporary orthodontic treatment and be aware of its limitations.

This is exactly what this programme delivers, along with the vital core competencies to excel in this unique postgraduate area.

The objective of the Orthodontics and Dentofacial Orthopaedics is to provide dental surgeons with the core knowledge and clinical skills essential for the delivery of competent, high-quality dental care.

It encompasses a blended learning approach and consists of a combination of lectures, seminars, workshops, clinical case discussions and laboratory sessions.

The flexible modular structure enables the programme to meet a wide range of student needs.



ORTHODONTIC INSTRUMENTS

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente

- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان



SIMZ-YM5-001
Height Gauge 0.18
Aluminum



SIMZ-YM5-002
Height Gauge 0.22
Aluminum



SIMZ-YM5-003
Height Gauge 0.18
Stainless Steel



SIMZ-YM5-004
Height Gauge 0.22
Stainless Steel



SIMZ-YM5-005
Moveable
Height Gauge



SIMZ-YM5-007
Endo Finger
Ruler Span
Measure Scale



SIMZ-YM5-006
Bone Bracket Gauge



SIMZ-YM5-008
Endo Finger
Ruler Span
Measure Scale



SIMZ-YM5-009
Endo Finger
Ruler Span
Measure Scale



SIMZ-YM5-010
Bone Thickness Gauge
40mm, 180mm



SIMZ-YM5-011
Osteometer
Measuring Forceps For Jaw
Ridge Width measurement, scale
on both sides, rustfree, autoclavable



SIMZ-YM5-012
Positive
Postioner



SIMZ-YM5-013
Strip Saw



SIMZ-YM5-014
Compasso
De Willis




SIMZ-YM5-015
Endo Finger Ruler
Span Measure Scale



SIMZ-YM5-016
Endo Finger Ruler
Span Measure Scale

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان

Bracket Holder Tweezers



SIMZ-YM5-017 Deluxe Bracket Holder Tweezers	SIMZ-YM5-018 Deluxe Bracket Holder Tweezers	SIMZ-YM5-019 Direct Bond Bracket Holder Tweezers	SIMZ-YM5-020 Bracket Tweezers Slot Aligner	SIMZ-YM5-021 Direct Bond Bracket Holder Tweezers, Small Tip	SIMZ-YM5-022 Exact Bracket Holder Tweezers	SIMZ-YM5-023 Bracket Placing Tweezers	SIMZ-YM5-024 Easy Access Tweezers
--	--	---	---	---	---	--	--



SIMZ-YM5-025 Bracket Placer Tweezers (Extra Fine Tip)	SIMZ-YM5-026 Deluxe Bracket Holder Tweezers 135mm	SIMZ-YM5-027 Easy Access Tweezers 136mm	SIMZ-YM5-028 Slipper Tweezers 136mm	SIMZ-YM5-029 Direct Bonding Bracket Holders 135mm	SIMZ-YM5-030 Direct Bonding Bracket Holders 136mm	SIMZ-YM5-031 Direct Bonding Bracket Holders 135mm	SIMZ-YM5-032 Deluxe Bracket Holder 140mm	SIMZ-YM5-033 Direct Bonding Bracket Holders 140mm
---	---	---	---	---	---	---	--	---

ORTHODONTIC INSTRUMENTS

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente

- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان



SIMZ-YM5-034
Elastic instruments



SIMZ-YM5-035
Twiel-on elastic module instrument




SIMZ-YM5-036
Schure Band Seater



SIMZ-YM5-037
Ortho light scaler



SIMZ-YM5-038
Ortho light scaler



SIMZ-YM5-039
Ortho instrument



SIMZ-YM5-040
Ortho instrument



SIMZ-YM5-041
Ligature instruments



SIMZ-YM5-042
Ligature instruments

ORTHODONTIC INSTRUMENTS

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان



SIMZ-YM5-043
Ortho Light double ended
ligature director



SIMZ-YM5-044
Ortho Light double ended ligature scaler director tivate coil spring
are remove modules with one instruments



SIMZ-YM5-045
Ortho band pusher scaler



SIMZ-YM5-046
Ortho ligature director single ended



SIMZ-YM5-047
Ortho Light double ended ligature scaler/director
tivate coil spring are remove modules with one instruments



SIMZ-YM5-048
Ortholite Ligature Director Scaler



SIMZ-YM5-049
Ortho light ligature director double ended



SIMZ-YM5-050
Double ended Tampers

ORTHODONTIC INSTRUMENTS

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان



SIMZ-YM5-051
Band seater
auto cleavable



SIMZ-YM5-052
Band seater
auto cleavable
Triangle Shape



SIMZ-YM5-053
Band seater
auto cleavable
Square Shape



SIMZ-YM5-054
Band seater
auto cleavable
Round Shape



SIMZ-YM5-055
Mershon Band
Pusher
140mm



SIMZ-YM5-056
Band seating &
removing
instruments



SIMZ-YM5-057
Ortho Light double
ended
ligature director
scaler



SIMZ-YM5-058
Ortho Light Band
Pusher Scaler



SIMZ-YM5-059
Band Seating &
Removing
Silicon Grip



SIMZ-YM5-060
Ligature Instrument
Silicon Grip



SIMZ-YM5-061
Elastic Instrument
Silicon Grip



SIMZ-YM5-062
Distal Bender
Silicon Grip



SIMZ-YM5-063
Distal Bender
Silicon Grip



SIMZ-YM5-064
Distal Bender
Silicon Grip

PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- מכשירי לתקويم الأسنان بروסתھيس



SIMZ-YM5-065
Aderer
Large wire bending
and loop forming pliers
130mm




SIMZ-YM5-066
Aderer Plier
For the manufacturing of appliances and
activation of wires for wires up to .030" .fine
tip and up to .036" mean tip.
Length 130mm




SIMZ-YM5-067
Angle Plier
Fig. 1
135mm



SIMZ-YM5-068
Angle Plier
Fig. 2
140mm



SIMZ-YM5-069
Angle Plier
Fig. 3
130mm



SIMZ-YM5-070
Angle Plier
Fig. 4
140mm

PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان بروتیسیس



SIMZ-YM5-071
Adams Plier
Adams pliers for bending wire arches and bends right angles for wires up to .022", .025 tip to beak. Length 130mm



SIMZ-YM5-072
Flat Nose Pliers
140mm



SIMZ-YM5-073
Johnson Contouring Pliers
140mm



SIMZ-YM5-074
Falcon Beak
Contouring Pliers.
130mm



SIMZ-YM5-075
How Plier
Universally usable tool for the multiband technology for wires up to .028". Length 145mm



SIMZ-YM5-076
How Plier curved
Universally usable tool for the multiband technology for wires up to .028". Length 145mm

PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pinces pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان برومشتیوس



SIMZ-YM5-077
De La Rosa
Fig. 1
135mm



SIMZ-YM5-078
Arch Counting Plier Grooved
for double back and triple back
bands in 0.30\" / 0.76mm & 0.36\"
/ 0.91mm lingual arch wires
135mm



SIMZ-YM5-079
Contouring Plier
Fig. 2
135mm



SIMZ-YM5-080
Nance Pliers
135mm



SIMZ-YM5-081
Tweed Loop
forming Pliers
140mm



SIMZ-YM5-082
Nance Loop
forming Pliers
140mm

PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pinces pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتفويم الأسنان بروستیس



SIMZ-YM5-083
Tweed Loop Plier
130mm



SIMZ-YM5-084
Band Removing Plier
The curved blade of the band removing plier can strip easily and gently remove directly to the tooth. Length 140mm



SIMZ-YM5-085
Crown & Shell Plier
130mm



SIMZ-YM5-086
Young Wire Bending &
Loop Forming Pliers
135mm



SIMZ-YM5-087
Ligature Cutters T.C
Hard Wire
130mm




SIMZ-YM5-088
Mini Ligature Cutters T.C
Light Tip
130mm

PLIERS FOR ORTHODONTICS & PROSTHETICS


- Щипцы для ортодонтии и протезирования
- Щипці для ортодонції та протезування
- Pincers pour l'orthodontie et la prothèse
- Alicates para ortodoncia y prótesis
- Zangen für Kieferorthopädie und Prothetik
- قماشة لتقويم الأسنان بروشيشين



SIMZ-YM5-089
Mini Ligature cutters T.C.
fine tip
130mm



SIMZ-YM5-090
Large
Fluch-cut distal and cutting Pliers T.C.
with safety hold fitted with
un-breakable, corrosion
resistant and durable
stainless steel tips
130mm



SIMZ-YM5-091
Small
Fluch-cut distal and cutting Pliers T.C.
with safety hold fitted with
un-breakable, corrosion
resistant and durable
stainless steel tips
130mm



SIMZ-YM5-092
Bracket Removing Plier
Bracket removing plier for comfortable
and quick removal of the brackets.
Length 135mm



SIMZ-YM5-093
Bracket Removing Plier curved
Bracket Removing Plier curved for
easy removal of brackets in the
posterior region. Length 145mm



SIMZ-YM5-094
Grooved Fine Plier
145mm

PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pinces pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان بروستیس



SIMZ-YM5-095
Aderer
Facebow loop forming Pliers
145mm




SIMZ-YM5-096
V Bending Pliers
135mm




SIMZ-YM5-097
Tweed / angle rectangular
arch forming & Torquing pliers
135mm



SIMZ-YM5-098
Weingart Plier medium
Weingart plier that are ideal due to
their curved tip for setting sheets and
is capable of bending the ends.
Length 145mm



SIMZ-YM5-099
Universal Plier
165mm



SIMZ-YM5-100
Jarabak Plier
145mm

PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincas para ortodoncia y prótesis
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمانة لتقويم الأسنان برومشتيس



SIMZ-YM5-101
Waldsachs Wire Bending
and Cutting Pliers
150mm



SIMZ-YM5-102
Universal Dolphin Style
Wire Bending & cutting pliers



SIMZ-YM5-103
Arrow Clasp Forming Pliers
125mm, 5"

SIMZ-YM5-104
Schwarz
130mm



SIMZ-YM5-105
Flat Nose Pliers
both sides smooth
140mm



SIMZ-YM5-106
Flat Nose Pliers
both sides serrated
140mm



SIMZ-YM5-107
Peeso
Wire Bending Pliers
135mm



SIMZ-YM5-108
Goslee
Flat Nose Pliers
140mm



SIMZ-YM5-109
Flat Nose Pliers
both sides serrated
with TC 140mm inserts

PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتفويم الأسنان بروتيسيس



SIMZ-YM5-110
Universal Pliers
for soft and hard wire
160mm



SIMZ-YM5-111
Thick grooved plier
145mm



SIMZ-YM5-112
Bird's beak claw with cutter
130mm



SIMZ-YM5-113
Elastomeric Separator
Placing Pliers 160mm



Separators placing pliers
for easy placement of radiopaque
and non-radiopaque separators



SIMZ-YM5-114
for incisors with
convex beaks
160mm




SIMZ-YM5-115
For bicuspid
and molars
with concave beaks
160mm

- Щипцы для ортодонтии и протезирования
- Pincas pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- قماشة لتقويم الأسنان برومشیس

CUTTERS



SIMZ-YM5-116
Hard Wire
Cutter 15 T.C




SIMZ-YM5-120
Mini Pin &
Ligature Cutter T.C




SIMZ-YM5-117
Distal End Cutter
Flush Cut T.C



SIMZ-YM5-121
Hard Wire Cutter
Straight T.C



SIMZ-YM5-118
Universal Cut & Hold
Distal End Cutter T.C



SIMZ-YM5-122
Micro-Mini Pin &
Ligature Cutter T.C



SIMZ-YM5-119
Pin & Ligature
Cutter T.C



SIMZ-YM5-123
Cutter T.C

PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان بروتستیس



CUTTER, UTILITY PLIERS & WIRE FORMING PLIERS




SIMZ-YM5-124
Weingart Pliers



SIMZ-YM5-129
Adhesive
Removing Pliers



SIMZ-YM5-125
Slim Weingart
Pliers



SIMZ-YM5-130
Posterior Band Removing
Pliers - Long




SIMZ-YM5-126
How Pliers
Straight



SIMZ-YM5-131
Posterior Band Removing
Pliers - Short



SIMZ-YM5-127
How Pliers Offset




SIMZ-YM5-132
Bird Beak Pliers



SIMZ-YM5-128
Band seating Pliers



Replacement tip and
hex key included.



SIMZ-YM5-133
Omega Loop
Forming Pliers

PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- مسامحة لتقويم الأسنان بروستھيس




SIMZ-YM5-134
Three Jaw Pliers



SIMZ-YM5-138
Jarabak Pliers




SIMZ-YM5-135
Light Wore Bird Beak
Pliers with Groove




SIMZ-YM5-139
Hollow Chop
Contouring Pliers





SIMZ-YM5-136
Tweed loop
Forming Pliers



SIMZ-YM5-140
Lingual Arch
Forming Pliers



SIMZ-YM5-137
Arch Bending Pliers



SIMZ-YM5-141
Rectangular Arch
Bending Pliers

- Щипцы для ортодонтии и протезирования
- Щипці для ортодонтії та протезування
- Pincers pour l'orthodontie et la prothèse
- Alicates para ortodoncia y prótesis
- Zangen für Kieferorthopädie und Prothetik
- كمانشة لتقويم الأسنان بروتستيس

L-KEY JOINT PLIERS T.C.



SIMZ-YM5-142
T.C Distal end cutter with safety hold, mini head, long handles.
fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" x .025")
- ⊕ 0.38 mm to 0.53 mm (.0155" to .0215")
- ☒ Max. 0.53 mm x 0.64 mm (.021" x .025")



SIMZ-YM5-146
T.C Soft wire and ligature cutter, 15° angled fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ☒ .25 mm to .38 mm (.010" to .015")




SIMZ-YM5-143
T.C Flush-cut distal end cutter with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" to .025")



SIMZ-YM5-147
T.C Soft wire and ligature cutter, 45° angled fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ☒ .25 mm to .38 mm (.010" to .015")



SIMZ-YM5-144
T.C Soft wire and ligature cutter, Mini fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ☒ .25 mm to .38 mm (.010" to .015")



SIMZ-YM5-148
T.C Hard wire cutter, straight Fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



SIMZ-YM5-145
T.C Soft wire and ligature cutter, standard fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ☒ .25 mm to .38 mm (.010" to .015")



SIMZ-YM5-149
T.C Hard wire cutter, 15° angled Fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



SIMZ-YM5-150
T.C Hard wire cutter, 15° angled Non-autocalavable

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincettes pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- مساملة لتقويم الأسنان بروستhetics

L-KEY JOINT PLIERS



SIMZ-YM5-151
Band removing pliers



SIMZ-YM5-152
Bracket removing pliers angled



SIMZ-YM5-153
Bracket removing pliers straight



SIMZ-YM5-154
Bracket removing pliers with plastic tip



SIMZ-YM5-155
Adhesive removing pliers



SIMZ-YM5-156
Weingart-mini
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



SIMZ-YM5-157
Weingart
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



SIMZ-YM5-158
How Straight
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



SIMZ-YM5-159
How Angled
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



SIMZ-YM5-160
Wire bending and loop forming pliers
Angled (Bird beak)

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")

PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pinces pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمانشة لتقويم الأسنان بروتیسیس

L-KEY JOINT PLIERS



SIMZ-YM5-161
Cinch Back / tie back
Distal end bending pliers

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



SIMZ-YM5-162
Detailing step pliers / Bayonet
bending pliers - 1mm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



SIMZ-YM5-163
Stop pliers V-Bebd - 1mm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.43 mm x 0.56 mm (.017" x .028")



SIMZ-YM5-164
Tweed / angel
Rectangular arch forming and
torquing pliers

● 0.30 mm to 0.71 mm (.012" to .028")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



SIMZ-YM5-165
Surgical ball hook crimping
pliers, straight

L-KEY JOINT PLIERS




SIMZ-YM5-166
Multi-purpose pliers, bird beak,
De la rosa and cutter in one

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.46 mm x 0.56 mm (.018" x .022")



SIMZ-YM5-167
Hammer head
Distal end bending pliers

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



SIMZ-YM5-168
Nance
Loop forming pliers

● 0.36 mm to 0.71 mm (.014" to .028")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



SIMZ-YM5-169
Lingual / palatal arch forming plies

● 0.76 mm (.030") / 0.91 mm (.036")



SIMZ-YM5-170
Torquing pliers



SIMZ-YM5-171
Surgical ball hook crimping
pliers, angled

PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- عسائنة لتقويم الأسنان وبروسثيسيس

L-KEY JOINT PLIERS




SIMZ-YM5-172
Aderer Clamp three jaw wire bending plier for contouring and bending wire up to .30"/.76mm.



SIMZ-YM5-178
Falcon beak contouring plier.



SIMZ-YM5-173
Bird beak plier for bending round, square and rectangle wire up to .30"/.76mm.



SIMZ-YM5-179
For placement and removal of arch wires, pins and other auxiliaries.




SIMZ-YM5-174
Bird beak plier for light wire bending round, square and rectangular up to .016"/0.41mm.



SIMZ-YM5-180
For placement and removal of arch wires, pins and other auxiliaries in angle design



SIMZ-YM5-175
Double Rounded jaw bending pliers, hooks and wire up to .020"/0.51mm.



SIMZ-YM5-181
contouring Plier.



SIMZ-YM5-182
Bending plier. In order to skirt molares bands. With grooved. Recommended wire. .016", .018" and .022" 0.41mm, 0.46mm and 0.56mm




SIMZ-YM5-176
Adams pliers for hooks and bending wire up to .030"/0.76mm.



SIMZ-YM5-183
Loop forming in four length 3,4,5 & 6 mm wire capacity 0.22", 56mm.



SIMZ-YM5-177
Johnson pliers for crown band contouring



SIMZ-YM5-184
Arch bending plier design for convenient handling of square or rectangular wire, .020" / 0.51mm.



PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان بروستھیس



PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincas pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- عسائنة لتقويم الأسنان وبروستهتیس

LAP & L-KEY JOINT PLIERS



SIMZ-YM5-193
Pin & Ligature Cutters
TC. Gold.



SIMZ-YM5-194
Mini Pin & Ligature
Cutter TC. Gold.



SIMZ-YM5-195
Standard Ligature
Cutter TC. Gold.



SIMZ-YM5-196
Distal end Cutter
TC. Gold.



SIMZ-YM5-197
Distal end Cutter
Long Handle
TC. Gold.



SIMZ-YM5-198
Hook Crimping Plier Angled
TC. Gold.



SIMZ-YM5-199
Distal end Cutter, Safety Hold Plier
(Ideal for Cutting Wire) TC. Gold.

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان بروستھتیس

L-KEY JOINT PLIERS



SIMZ-YM5-200
Direct Bracket Remover
TC. Gold.



SIMZ-YM5-201
Straight Direct Bond Bracket
Remover T.C



SIMZ-YM5-202
Pin Holding Plier



SIMZ-YM5-203
3 Prong Clasp Adjustable
Plier TC. Gold.



SIMZ-YM5-204
Utility Plier Weingart TC. Gold.



SIMZ-YM5-205
Nance Plier TC. Gold.

PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمانة لتقويم الأسنان بروستhetics



SIMZ-YM5-206
Separating Pliers For Separators
155mm



SIMZ-YM5-207
Elastic
Separating Plier



SIMZ-YM5-208 Straight
SIMZ-YM5-209 Curved
Smaha Elastomeric ligature
placing forceps
140 mm



SIMZ-YM5-210
Elastomeric ligature placing
forceps with groove
140 mm



SIMZ-YM5-211
Elastomeric ligature placing
forceps with groove
140 mm

PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Щипці для ортодонтії та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان بروسلینس
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik



SIMZ-YM5-212
Methieu Plier
With Ring



SIMZ-YM5-213
Pliers to Place Modules
140mm



SIMZ-YM5-214
Forceps For Placing Dentalastics
120mm



SIMZ-YM5-215
Mathieu
140mm




SIMZ-YM5-216
Mathieu Ligature Pliers
140mm



SIMZ-YM5-217
Ligature Pliers T.C Gold
With Tungsten Carbide Inserts
140mm



Ligature Pliers
140mm




SIMZ-YM5-218 **SIMZ-YM5-219**



SIMZ-YM5-220

PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كماناتة لتقويم الأسنان بروستھتیس



SIMZ-YM5-221
Plier Stand
Stainless Steel.



SIMZ-YM5-222
Plier Stand
Stainless Steel.



SIMZ-YM5-223
Pliers Stand
Stainless Rod
(Aluminum Anodized)
Stainless Steel




SIMZ-YM5-224
Plier Stand
Stainless Steel.



SIMZ-YM5-225
Plier Stand
Stainless Steel.



SIMZ-YM5-226
Plier Stand
Plastic



Instruments for dental implants are made of high quality stainless steel.

Our line of instruments include, but is not limited to, the following products:

Bone Grafting Instruments

- Bone Currettes (SPRATT) – designed for harvesting bone
- Bone Grinder – used in conjunction with a mallet to grind bone into appropriate particle sizes
- Bone Basin – weighted bottom design prevents tipping or spilling while mixing bone or other bone grafting materials

Bone Compactors – feature serrated tips for compacting bone

Scalpels & Biopsy Punches

- Scalpels – offers convenience of a disposable with razor sharp cutting edge for superior cutting
- Biopsy Punches – designed at a steeper angle to ensure a finely honed, sharper edge
- Bone Spreading Chisels
- Freer Chisels – double or single cut tips used for widening ridges during implant preparations

- Graduated Chisels – feature reference markings for measurement during a bone ridge augmentation and are available in straight and curved patterns

Osteotomes

- Concave Osteotomes – break the floor of the maxillary sinus through the implant cavity and are available in curved (off set) and straight (with or without stops)
- Convex Osteotomes – raise the floor of the maxillary sinus and are available in curved (off set) and straight patterns (with or without stops)
- Progressive Osteotomes – designed to prepare the implant site in the maxillary area where bone density is 3 or 4
- Bone Spreading Osteotomes – sharp or blunt configurations for progressive circumferential expansion of the implant site

Sinus Lift Instruments

- Sinus Lift Currettes and Elevators – lightweight, wide diameter stainless steel hand instruments lift the sinus membrane for bone augmentation procedures

Gauges & Calipers

- Gauges & Calipers – various measuring devices available to determine appropriate length, width, and distance for implant procedures

Collagen Products

- Collagen wound dressing product offering includes products.

Sterilization Trays

- Sterilization Trays – designed to protect and store instruments during sterilization, transportation and storage



IMPLANT GAGES



SIMZ-YM6-001
Multifunctional Meter
Length 90 mm,
different sizes,
8-10-11mm,
corrosion-resistant

SIMZ-YM6-002
Multifunctional Meter
Length 90 mm,
different sizes,
5-13-15-17 mm,
corrosion-resistant

SIMZ-YM6-003
Measuring Probe
Length 170 mm,
(6-8-10-14-16 mm),
corrosion-resistant

SIMZ-YM6-004
Implant
Depth Gauge

SIMZ-YM6-005
Implant
Depth Gauge
& Key

SIMZ-YM6-006
Implant
Depth Gauge

SIMZ-YM6-007
Implant
Depth Gauge

IMPLANT GAGES



SIMZ-YM6-008
Double Ended Osteotomy Depth Gauge

SIMZ-YM6-009
Double Ended Osteotomy Depth Gauge

SIMZ-YM6-010
Double Ended Osteotomy Depth Gauge

SIMZ-YM6-011
Depth & Gingival Gauges

Caliper / Calipers Rule



SIMZ-YM6-012
Castroviejo Calipers
20mm, 45°
90mm



SIMZ-YM6-013
Castroviejo Calipers
20mm
90mm



SIMZ-YM6-014
Castroviejo Calipers
20mm
90mm



SIMZ-YM6-015
Ridge Mapping
Double Caliper



SIMZ-YM6-016
Fig. 5
Castroviejo
18cm



SIMZ-YM6-017
Fig. 4
Castroviejo
18cm



SIMZ-YM6-018
Fig. 3
Castroviejo
18cm

- Инструменты для имплантологии
- Implantologie
- Implantologie

- Інструменти для імплантології
- Implantologia
- زراعة الأسنان



SIMZ-YM6-019
Clamp, Osteometer
Length about 155 mm,
ridge width measurement,
two-sided scale,
measuring
range 0 - 25 mm



SIMZ-YM6-020
Measuring Loop
straight length 90 mm,
measuring range
0 - 40 mm,



SIMZ-YM6-021
Measuring Loop
angled, length 90 mm,
measuring range
0 - 40 mm,



SIMZ-YM6-022
Measuring loop. Castroviejo
angled, length 175 mm,
measuring range 0 - 40 mm,
with window



SIMZ-YM6-023
Measuring Loop
angled, length
175 mm,
measuring range
0 - 40 mm,



SIMZ-YM6-024
Measuring Loop
straight length
175 mm,
measuring range
0 - 40 mm,



SIMZ-YM6-025
Fig. 1
Iwanson
spring caliper
For metal and
porcelain.



SIMZ-YM6-026
Fig. 2
Iwanson
spring caliper
For wax.



SIMZ-YM6-027
Bone Thickness Gauge
40mm, 180mm

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для имплантології
- Implantologia
- زراعة الأسنان



SIMZ-YM6-028
Boley Gauge



SIMZ-YM6-029
Vernier Gauge



SIMZ-YM6-030
Aluminium Caliper
up to 200mm



SIMZ-YM6-031
Boley Gauge



SIMZ-YM6-032
Delux Boley
Gauge



SIMZ-YM6-033
Beerendonk
130mm
Mini Gauge

- Инструменты для имплантологии
- Implantologie
- Implantologie

- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان



SIMZ-YM6-034
Micro Boley
Gauge Twisted



SIMZ-YM6-035
Micro Boley
Gauge Straight



SIMZ-YM6-036
Dentaurum Ruler




SIMZ-YM6-037
Ruler Wedge



SIMZ-YM6-038
Ruler




SIMZ-YM6-039
Ortho Ruler / Cone Socket




SIMZ-YM6-040
Dynamometer

BLACK TITANIUM COLOUR COATED INSTRUMENTS
Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



SIMZ-YM6-041	SIMZ-YM6-042	SIMZ-YM6-043	SIMZ-YM6-044	SIMZ-YM6-045	SIMZ-YM6-046	SIMZ-YM6-047	SIMZ-YM6-048	SIMZ-YM6-049	SIMZ-YM6-050
Fig. 1	Fig. 2	Fig. 3	Micro-Mini	Mini 1	Mini 4	Fig. 4	Fig. 5	Fig. 6	Mini 3
		Extra-Flex			Extra-Flex	Extra-Flex	Flexi-Thin	Flexi-Thin	Extra-Flex

BLACK TITANIUM COLOUR COATED INSTRUMENTS
 Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



SIMZ-YM6-051
 AB1
 Boghosian

SIMZ-YM6-052
 AB2
 Boghosian

SIMZ-YM6-053
 W3

SIMZ-YM6-054
 8A

SIMZ-YM6-055
 4/5
 Gregg

SIMZ-YM6-056
 Interproximal
 Carver

SIMZ-YM6-057
 Interproximal
 Carver, Long


SIMZ-YM6-058
 A6

SIMZ-YM6-059
 Contouring
 instrument

BLACK TITANIUM COLOUR COATED INSTRUMENTS
 Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

BURNISHERS

FREEDMAN BURNISHERS



SIMZ-YM6-060
21B

SIMZ-YM6-061
27/29

SIMZ-YM6-062
2
Ladmore

SIMZ-YM6-063
3
Ladmore


SIMZ-YM6-064
Freedman
"Duckhead"

SIMZ-YM6-065
Freedman
Small Contact
Forming

SIMZ-YM6-066
Freedman
Large Contact
Forming

SIMZ-YM6-067
Small/
Medium Ball
Burnisher

BLACK TITANIUM COLOUR COATED INSTRUMENTS
 Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



SIMZ-YM6-068
A
 Cosmetic
 Contouring

SIMZ-YM6-069
B
 Cosmetic
 Contouring

SIMZ-YM6-070
C
 Cosmetic
 Contouring

SIMZ-YM6-071
D
 Cosmetic
 Contouring

SIMZ-YM6-072
E
 Cosmetic
 Contouring

SIMZ-YM6-073
F
 Cosmetic
 Contouring


SIMZ-YM6-074
 Marginal Ridge
 & Embrasure
 Shaping Instrument

SIMZ-YM6-075
 Occlusal
 Anatomy
 Instrument

SIMZ-YM6-076
 Composite
 Packing
 Instrument

BLACK TITANIUM COLOUR COATED INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



SIMZ-YM6-077	SIMZ-YM6-078	SIMZ-YM6-079	SIMZ-YM6-080	SIMZ-YM6-081	SIMZ-YM6-082	SIMZ-YM6-083	SIMZ-YM6-084	SIMZ-YM6-085	SIMZ-YM6-086
Fig. 3	Mini 1	Micro-Mini	Medium	Large	OT	3	5A	Small/Medium	Medium/Large
Extra-Flex			Placing/ Condensing	Placing/ Condensing	Tanner	Hollenback		Contact Forming	Contact Forming

TITANIUM PLASMA COLOUR COATED SINUS LIFT INSTRUMENTS
 Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



SIMZ-YM6-087
Sinus Lift Elevator



SIMZ-YM6-088
Sinus Lift Elevator



SIMZ-YM6-089
Sinus Lift Elevator



SIMZ-YM6-090
Sinus Lift Elevator



SIMZ-YM6-091
Sinus Lift Elevator

TITANIUM PLASMA COLOUR COATED SINUS LIFT INSTRUMENTS
 Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



SIMZ-YM6-092
Sinus Lift Elevator



SIMZ-YM6-093
Sinus Lift Elevator




SIMZ-YM6-094
Sinus Lift Elevator



SIMZ-YM6-095
Sinus Lift Elevator

TITANIUM PLASMA COLOUR COATED IMPLANT SINUS LIFT INSTRUMENTS
 Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



SIMZ-YM6-096
Acute 6577

SIMZ-YM6-097
Sinus1

SIMZ-YM6-098
Sinus2

SIMZ-YM6-099
Sinus3


SIMZ-YM6-100
Sinus4

SIMZ-YM6-101
Sinus5

SIMZ-YM6-102
Sinus6

SIMZ-YM6-103
Sinus7

TITANIUM PLASMA COLOUR COATED SINUS LIFT INSTRUMENTS
 Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



SIMZ-YM6-104
Sinuslift 1

SIMZ-YM6-105
Sinuslift 2

SIMZ-YM6-106
Sinuslift 3

SIMZ-YM6-107
Sinuslift 4

SIMZ-YM6-108
Sinuslift 5

SIMZ-YM6-109
Sinuslift 6

SIMZ-YM6-110
Membrane
3Positioning 1

SIMZ-YM6-111
Condenser 1

TITANIUM PLASMA COLOUR COATED INSTRUMENTS
Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

Sinus Lift Instrument,
Flexible Sinus Lift instruments
reduce the risk of perforation of the
"Schneider's membrane"

Sinus Lift Elevator,
9.5mm Straight grip
Length 175 mm, Sinus lift elevators, dull,
for gentle peeling and lifting the
Schneider membrane

Sinus lift instrument
9.5mm Straight grip Length 180 mm,
Buser (B1, B2, B3) raspatories,
quality medical stainless steel,
corrosion-resistant



SIMZ-YM6-112 Fig. 1 SIMZ-YM6-113 Fig. 2 SIMZ-YM6-114 Fig. 3



SIMZ-YM6-115 Fig. 1 SIMZ-YM6-116 Fig. 2 SIMZ-YM6-117 Fig. 3




SIMZ-YM6-118 B01 SIMZ-YM6-119 B02 SIMZ-YM6-120 B03

- Инструменты для имплантологии
 - Implantologie
 - Implantologie
- Инструменти для имплантології
 - Implantologia
 - زراعة الأسنان

TITANIUM PLASMA COLOUR COATED INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

Sinus lift instrument 9.5mm Straight grip Sinus lift elevators, blunt or sharply curved (H01, H02, K03, K04, K78, K77)



SIMZ-YM6-121
Sinus lift instrument
H01

SIMZ-YM6-122
Sinus lift instrument
H02

SIMZ-YM6-123
Sinus lift instrument
K03

SIMZ-YM6-124
Sinus lift instrument
K04

SIMZ-YM6-125
Sinus lift instrument
K78

SIMZ-YM6-126
Sinus lift instrument
K77


- Инструменты для имплантологии
- Implantologie
- Implantologie

- Інструменти для імплантології
- Implantologia
- زراعة الأسنان

TITANIUM PLASMA COLOUR COATED INSTRUMENTS
Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

SINUS LIFT INSTRUMENT

Sinus Lift Elevators, 9.5mm Hollow Handle length 170mm Blunt Curved Set of 6 Pcs



SIMZ-YM6-127

SIMZ-YM6-128

SIMZ-YM6-129

SIMZ-YM6-130

SIMZ-YM6-131

SIMZ-YM6-132

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm
PERIODONTAL SURGERY SURGICAL CHISELS
KRAMER NEVINS SINUS LIFT



SIMZ-YM6-133 TG



SIMZ-YM6-134 SOLT3/4



SIMZ-YM6-135 TGO



SIMZ-YM6-136 29K



SIMZ-YM6-137 FEDI4



SIMZ-YM6-138 OCH3



SIMZ-YM6-139 OCH4



SIMZ-YM6-140 FEDI1



SIMZ-YM6-141 FEDI2



SIMZ-YM6-142 FEDI3



SIMZ-YM6-143 Gf36

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm
IMPLANT SURGICAL CHISELS



SIMZ-YM6-144 13K/TG



SIMZ-YM6-145 13K/TG Titanium Plasma Coated Point



SIMZ-YM6-146 13K/TG



SIMZ-YM6-147 13KL/R Titanium Plasma Coated Point



SIMZ-YM6-148 13KL 4/6mm



SIMZ-YM6-149 Pt2 Titanium Plasma Coated Point
A/V Regular



SIMZ-YM6-150 Pt1 Titanium Plasma Coated Point
A/V Regular



SIMZ-YM6-151 Pt3 Titanium Plasma Coated Point Angel S/E



SIMZ-YM6-152 Pt4 Titanium Plasma Coated Point STR S/E



SIMZ-YM6-153 Pt5 Titanium Plasma Coated Point TG CVD S/E

TITANIUM PLASMA COLOUR COATED INSTRUMENTS
Hollow Handle 10mm
MIXING SPATULA WITH SHOVEL



SIMZ-YM6-154
 Bone Carrier
 Titanium
 185mm



SIMZ-YM6-155
 Raspator, Prichard
 Length 180 mm,
 10 mm hollow handle,
 double-ended
 10.5 mm | 4.5 mm



SIMZ-YM6-156
 Raspator, Prichard
 Length 180 mm,
 10 mm hollow handle,
 double ended



SIMZ-YM6-157
 Raspator, Prichard
 Length 180 mm,
 10 mm hollow handle,
 double-ended
 10.5 mm | 4.5 mm



SIMZ-YM6-158
 Mixing Spatula
 with shovel
 170mm



SIMZ-YM6-159
 Mixing Spatula
 with shovel
 185mm

BONE GRAFTING INSTRUMENTS COMPACTOR /TAMPERS

Product offering includes instruments for harvesting, mixing and compacting bone Provides clinicians with the proper instrumentation for bone grafting procedures Weighted bone basin Prevents tipping and spilling of bone and other materials during mixing

Hollow Handle 9.5mm Length 170mm



SIMZ-YM6-160
Graft Packer,
Concave Tips
DE 4/6mm

SIMZ-YM6-161
Bone Spoon/
Packer DE 3mm

SIMZ-YM6-162
Graft Packer,
Serrated Tips
DE 6/8mm

SIMZ-YM6-163
Graft Packer,
Serrated Tips
DE 3/4mm

SIMZ-YM6-164
Graft Packer,
Smooth Tips
DE 5/8mm

SIMZ-YM6-165
Titanium graft packer/
placement site depth
gauge 3.3mm & 4.0mm
Diameter tips Marked
at 8, 11, 13, 15 & 18mm

SIMZ-YM6-166
Graft Material Spoon
Concave Tip For
Larger Grafts

**BONE GRAFTING INSTRUMENTS
TITANIUM PLASMA COLOUR COATED
Hollow Handle Ø 10 mm**



SIMZ-YM6-167
Blade and pusher
Length 180 mm,
for straight or angled,
blade: 15 x 5 mm,
pusher: 4 mm diameter
Straight

SIMZ-YM6-168
Blade and pusher
Length 180 mm,
for straight or angled,
blade: 15 x 5 mm,
pusher: 4 mm diameter
angular

SIMZ-YM6-169
Spoon and pusher
Length 180 mm,
for angled spoon:
10 x 5 mm, pusher:
diameter 3 mm

SIMZ-YM6-170
Tunnelierungs
Instrument
length 175 mm,

SIMZ-YM6-171
Tunnelierungs
Instrument
length 175 mm,
different designs

SIMZ-YM6-172
Tampers
Length 170 mm,
for bone material,
double-ended,
diameter 3.2 mm
(2-4-6-8-10) and diameter
2.5mm (2-4-6-8-10)

Pluggers

Hollow Handle 9.5mm Length 160 mm, with titanium lips for mixing and pouring of autologous, heterologous bone substitute material. diameter (3/4 or 6/8 mm)

BONE GRAFTING INSTRUMENTS



SIMZ-YM6-173 Ø 4mm
SIMZ-YM6-174 Ø 5mm
Compactor
150mm



SIMZ-YM6-175
Sinus Window
Bone Tampers
4 mm



SIMZ-YM6-176
Sinus Window
Bone Tampers
7 mm



SIMZ-YM6-177
6/8mm




SIMZ-YM6-178
3/4mm

- Инструменты для имплантологии
- Instrumente für Implantologie
- ɪmplantoloia
- زراعة الأسنان

SPOON AND PUSHER (TRINOVO)

Hollow Handle Ø 10 mm Length 180mm



SIMZ-YM6-179
10x5mm
Spoon and
Pusher
with safety
stop

SIMZ-YM6-180
10x5mm
Spoon and
Pusher
with safety
stop

SIMZ-YM6-181
10x5mm
Spoon and
Pusher
with safety
stop


SIMZ-YM6-182
Bone Graft/Carrier
Titanium
2-4-6-8-10mm

SIMZ-YM6-183
15x10mm
Spoon and
Pusher
with safety
stop

SIMZ-YM6-184
10x5mm
Spoon and
Pusher
with safety
stop

SIMZ-YM6-185
10x5mm
Spoon and
Pusher
with safety
stop

PERIOSTEALS ELEVATOR 10mm Hollow Handle



SIMZ-YM6-186
Raspator,
Bone Scraper
Length 185 mm,
10 mm hollow
handle,
double-ended 6
mm / 4 mm,

SIMZ-YM6-187
Raspator,
Raspatory
Length 180 mm,
10 mm hollow
handle,
double-ended
3 mm | 3 mm,

SIMZ-YM6-188
Raspator,
Raspatory
Length 180 mm,
10 mm hollow
handle,
double-ended
3 mm | 3 mm,


SIMZ-YM6-189
Raspator,
Papilla
Length 185 mm,
10 mm hollow
handle,
double-ended
3 mm | 2 mm,

SIMZ-YM6-190
Raspator,
Bone Scraper
Length 175 mm,
10 mm hollow
handle,
double-ended,

SIMZ-YM6-191
Raspator,
Micro-Freer
Length 180 mm,
10 mm hollow
double-ended
3 mm | 4 mm

SIMZ-YM6-192
Raspator,
Freer
Length 180 mm,
10 mm hollow
handle,
double-ended
5 mm | 5 mm,

PERIOSTEALS ELEVATOR 10mm Hollow Handle



SIMZ-YM6-193
Raspator,
Freer
10 mm hollow
handle,
high-quality
medical
double ended

SIMZ-YM6-194
Raspator,
Molt
Length 180 mm,
10 mm Hollow
grip,
double-ended
6 mm

SIMZ-YM6-195
Raspator,
Micro Freer
Length 185 mm,
10 mm hollow
handle,
double-ended


SIMZ-YM6-196
Raspator,
Micro Freer
Length 185 mm,
10 mm hollow
handle,
double-ended

SIMZ-YM6-197
Raspator, Micro
Raspator /
Chisel
Length 170 mm,
10 mm hollow
handle,
double-ended

SIMZ-YM6-198
Raspator, Micro
Raspator /
Chisel
Length 175 mm,
10 mm hollow
handle,
double-ended

SIMZ-YM6-199
Raspator,
Micro Raspator /
Chisel
Length 175 mm,
10 mm hollow handle,
double-ended
3.5mm | 3.5 mm

PERIOSTEALS ELEVATOR 10mm Hollow Handle



SIMZ-YM6-200
Raspatory Micro Molt
Length 175 mm,
corrosion-resistant

SIMZ-YM6-201
Raspatory Chipping
Workers:
Chisel 24G and 7.5 mm,
length about 175 mm,
corrosion-resistant

SIMZ-YM6-202
Raspatory Molt
Work ends:
5.0 mm and 4.0 mm,
length about 175 mm,
corrosion-resistant

SIMZ-YM6-203
Rasp, Molt M9
Length 175 mm,
corrosion-resistant

SIMZ-YM6-204
Raspator, Wautier
Length 175 mm,
corrosion-resistant

SIMZ-YM6-205
Raspator, Benque
Length 170 mm,
workers
5 mm | 7 mm,
corrosion-resistant

SIMZ-YM6-206
Raspator, Nordberg
Length 180 mm,
10 mm hollow handle,
double-ended
6 mm | 4 mm

Curettes
Hollow Handle 10mm Length 170mm



TISSUE MEMBRANE FORCEPS



SIMZ-YM6-216
Tissue Graft Forcep
Serrated Teeth
150mm



SIMZ-YM6-217
Self Locking
Grooved
150mm



SIMZ-YM6-218
Tissue Graft Forcep
Serrated Teeth
150mm




SIMZ-YM6-219
Titanium Metal
Locking Plier



SIMZ-YM6-220
TC Gripper Carbide
Tipped Tissue Forceps
Serrated Carbide Tips
145mm



SIMZ-YM6-221
Adson Carbide
Tipped Tissue Forceps
Serrated Carbide Tips
120mm



SIMZ-YM6-222
Chalazion Mucocele
Forceps
90mm Lenght
12mm x 23mm
Inside Ring

1x2



SIMZ-YM6-223
Micor Adson
Tissue Forcep 120mm
1x2 Teeth -
0.5mm Tips



SIMZ-YM6-224
Adson 7x7
Tissue Forceps
7x7 Teeth
120mm



SIMZ-YM6-225
Adson 7x7 Curved
Tissue Forceps
7x7 Teeth 120mm



SIMZ-YM6-226
Suture Pliers
Suture Forceps
with
1.6mm Ring Tip
160mm



SIMZ-YM6-227
Membrane Forcep
With Delicate
Ring Tip
10 cm



SIMZ-YM6-228
Debakey Atraumatic
Tissue Forceps
Straight
Serrated Teeth
160mm



SIMZ-YM6-229
Debakey Atraumatic
Tissue Forceps Curved
Serrated Teeth
160mm

TISSUE MEMBRANE FORCEPS



SIMZ-YM6-230
Russian Long
Tissue Forceps
Serrated Tip
200mm



SIMZ-YM6-231
Russian
Tissue Forceps
Serrated Tip
150mm



SIMZ-YM6-232
Tissue Punch Forceps
To Remove Tissue
After Using Punch
150mm



SIMZ-YM6-233
Narrow
135mm



SIMZ-YM6-234
105mm



SIMZ-YM6-235
wide
135mm



SIMZ-YM6-236
Bishop Harman
Micro Tissue Forcep
8.6cm 1x2 Teeth
- 0.3mm Tips



SIMZ-YM6-237
Pin-Tweezers
length 150mm
with lock



SIMZ-YM6-238
without lock
150mm



SIMZ-YM6-239
Haltepinzette with lock
155mm, for pins, towards
the tips, etc



SIMZ-YM6-240
Carry Forceps
150mm




SIMZ-YM6-241
Carry Forceps
150mm



SIMZ-YM6-242
Carry Forceps
150mm


OSTEOTOME CONVEXES / CONCAVE / CONIQUES



SIMZ-YM6-243
SET OF 5 PIECES
Convexes

SIMZ-YM6-244
SET OF 5 PIECES
Concaves

SIMZ-YM6-245
SET OF 5 PIECES
Coniques



	CONVEXES STRAIGHT	SIMZ-YM6-243	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONCAVES STRAIGHT	SIMZ-YM6-244	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONIQUES STRAIGHT	SIMZ-YM6-245	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm


OSTEOTOME CONVEXES / CONCAVE / CONIQUES



SIMZ-YM6-246
SET OF 5 PIECES
Convexes Curved

SIMZ-YM6-247
SET OF 5 PIECES
Concaves Curved

SIMZ-YM6-248
SET OF 5 PIECES
Coniques Curved



	CONVEXES CURVED	SIMZ-YM6-246	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONCAVES CURVED	SIMZ-YM6-247	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONIQUES CURVED	SIMZ-YM6-248	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm

TREPHINE

Burs with depth markings are suited for precise and simple transplantation of cylindrical bone grafts as well as for all other bone-extraction tasks.



SIMZ-YM6-249 Ø Int. 2 mm Ø Ext. 3 mm	SIMZ-YM6-250 Ø Int. 3 mm Ø Ext. 3.5 mm	SIMZ-YM6-251 Ø Int. 4 mm Ø Ext. 5 mm	SIMZ-YM6-252 Ø Int. 5 mm Ø Ext. 6 mm	SIMZ-YM6-253 Ø Int. 6 mm Ø Ext. 7 mm	SIMZ-YM6-254 Ø Int. 7 mm Ø Ext. 8 mm	SIMZ-YM6-255 Ø Int. 8 mm Ø Ext. 9 mm	SIMZ-YM6-256 Ø Int. 9 mm Ø Ext. 10 mm	SIMZ-YM6-257 Ø Int. 10 mm Ø Ext. 11 mm
---	---	---	---	---	---	---	--	---

Ø4.0 Ø5.0 Ø6.0 Ø7.0



SIMZ-YM6-258
Tissue Punch with pin
Set of 4 Pcs



SIMZ-YM6-259
Tissue Punch
Set of 5 Pcs



SIMZ-YM6-260
Tissue Punch
Set of 5 Pcs



SIMZ-YM6-261
Trephine Drill
Set of 8 Pcs



Ø2.0mm Ø2.5mm Ø2.8mm Ø3.2mm Ø3.65mm Ø4.2mm Ø5.2mm

SIMZ-YM6-262
Short Implant Drills
Set of 3 Pcs

SIMZ-YM6-263
Long Implant Drills
Set of 4 Pcs



SIMZ-YM6-264
Trephine Drill
Set of 8 Pcs

BONE SCRAPER STRAIGHT & CURVED



SIMZ-YM6-265
Bone Scraper
Straight

SIMZ-YM6-266
Bone Scraper
Curved

TISSUE PUNCH STRAIGHT & CURVED



Tissue Punch
Straight

Tissue Punch
Curved



SIMZ-YM6-279
Bone Mill



SIMZ-YM6-280
Basins For Bone



SIMZ-YM6-281
Bone Morselizer Mill Crusher
With Titanium Plates



SIMZ-YM6-282
Dental Bur Cleaning
Brush Copper

KEYS FOR IMPLANT SCREW



SIMZ-YM6-283



SIMZ-YM6-284



SIMZ-YM6-285



SIMZ-YM6-286

TORQUE WRENCH & RATCHET



SIMZ-YM6-287 Star Torque Wrench 6.35mm Universal



SIMZ-YM6-288 Round Torque Wrench 4.0mm Universal



SIMZ-YM6-289 Hex Driver With 4.0mm Driver



SIMZ-YM6-290 Star Torque Wrench 6.35mm Universal



SIMZ-YM6-291 Hex Torque Wrench 4.00mm Universal



SIMZ-YM6-292 Star Torque Wrench 6.35mm Universal



SIMZ-YM6-293 Star Torque Wrench 6.35mm Universal



SIMZ-YM6-294 Star Torque Wrench 6.35mm Universal



SIMZ-YM6-295 Star Torque Wrench 6.35mm Universal



SIMZ-YM6-296 Star Torque Wrench 6.35mm Universal



SIMZ-YM6-297 Star Ratchet Wrench 6.35mm Universal



SIMZ-YM6-298 Star Ratchet Wrench 6.35mm Universal



SIMZ-YM6-299 Star Ratchet Wrench 6.35mm Universal



SIMZ-YM6-300 Star Ratchet Wrench 4.0mm Universal



SIMZ-YM6-301 Round Ratchet Wrench 4.0mm Universal



SIMZ-YM6-302 Round Ratchet Wrench 3.35mm Universal




SIMZ-YM6-303 Star Ratchet Wrench 4.0mm Universal



SIMZ-YM6-304 Star Ratchet Wrench 6.35mm Universal

- Инструменты для имплантологии
 - Implantologie
 - Implantologie
- Інструменти для імплантології
 - Implantoloia
 - زراعة الأسنان



BONE EXPANDER SCREW & SURGICAL KITS



SIMZ-YM6-310
 Bone Universal Kit
 Drill: 2mm, 2.8mm, 3.2mm, 3.65mm, 4.2mm, 4.5mm, 4.8mm, 5.2mm
 Trephines Bit:
 3.0mm, 4.0mm, 5.0mm, 6.0mm, 7.0mm, 8.0mm, 9.0mm, 10mm.
 Driver:
 2.5mm, 2.5mm Hand Pcs
 1.25mm, 1.25mm Hand Pcs
 Short Driver 2.5mm
 Long Driver 2.5mm
 1.25mm, Short Hex Driver Hand Pcs
 1.25mm, Long Hex Driver Hand Pcs
 Torque Wrench 4.0mm



SIMZ-YM6-311
 Bone Expander Kit
 Expander: 2.6mm, 3.0mm, 3.4mm, 3.8mm, 4.3mm
 Ratchet Driver
 Parallel Pin
 Pilot Drill
 Drill 1.4mm
 Drill Driver 4mm
 Hand Driver 4mm
 Ratchet 4.0mm



SIMZ-YM6-312
 Bone Universal Kit
 Drill: 2mm, 2.8mm, 3.2mm, 3.65mm, 4.2mm, 4.5mm, 4.8mm, 5.2mm
 Trephines Bit:
 3.0mm, 4.0mm, 5.0mm, 6.0mm, 7.0mm, 8.0mm, 9.0mm, 10mm.

Bone Compression:
 2mm, 3.3mm, 2.8mm, 4.2mm, 3mm, 3.5mm
 Driver:
 2.5mm, 2.5mm Hand Pcs
 1.25mm, 1.25mm Hand Pcs
 Short Driver 2.5mm
 Long Driver 2.5mm
 1.25mm, Short Hex Driver Hand Pcs
 1.25mm, Long Hex Driver Hand Pcs
 Torque Wrench 6.35mm 1 Pcs



SIMZ-YM6-313
 Conical Drills: 20 Pcs
 Tissue Punches:
 3.0mm, 3.50mm, 4.0mm, 4.5mm, 5.0mm, 6.0mm, 6.50mm
 Bone Expander:
 1.6/3.8mm, 1.9/3.0mm, 2.1/3.6mm, 2.6/3.8mm, 2.9/4.3mm
 Drill Adapter: 1 Pcs
 1.25mm, 2.50mm Hex Driver Hand Pcs

Pilot Drill: 1.6mm
 Drill Bit: 1.4mm, 1.6mm, 1.8mm, 2.0mm, 2.2mm
 Ratchet Driver: 4mm/6.35, 1.25/2.50
 Rimer: 3.2/4.7mm, 3.7/4.7mm
 Parallel Pin: 4.5mm, 2.5mm
 Double step Drill: 2.4/2.8mm
 Depth Gauge: 1 Pcs
 Torque Wrench 6.35mm Universal: 1 Pcs



SIMZ-YM6-314
Drills kit (Black Coated)
 The Kit Includes:
Dental Drills Implant Kit 8 Pcs
Drill Sizes:
 2.0mm, 2.5mm, 2.8mm, 3.2mm,
 3.65mm, 4.2mm, 5.2mm, 5.5mm
 Marking Depths on the drills
 for accurate drilling.
 Sizes 6.0mm, 8.0mm, 10.0mm,
 11.5mm, 13.0mm & 16.0mm.



SIMZ-YM6-315
Drill Kit (Polished)
 The Kit Includes:
Dental Drills Implant Kit 8 Pcs
Drill Sizes:
 2.0mm, 2.5mm, 2.8mm, 3.2mm,
 3.65mm, 4.2mm, 5.2mm, 5.5mm
 Marking Depths on the drills for
 accurate drilling.
 Sizes 6.0mm, 8.0mm, 10.0mm,
 11.5mm, 13.0mm & 16.0mm.

The Kit Includes:

- 1.8mm Pilot Drill
- 1525 (2.0mm)
- 2535 (3.0mm)
- 3545 (4.0mm)
- 4555 (5.0mm)
- 1828 (2.3mm)
- 2838 (3.3mm)
- 3848 (4.3mm)
- 4858 (5.3mm)
- 2228 (2.5mm)
- 3238 (3.5mm)
- 4248 (4.5mm)
- 5258 (5.5mm)



SIMZ-YM6-316
Osseodensification Burs Kits

- Инструменты для имплантологии
- Implantologie
- Implantologie

- Инструменти для имплантології
- Implantoloia
- زراعة الأسنان



SIMZ-YM6-317

Compression Expander Screws Kit

The Kit Includes:

Kit Includes

3x Bone Compression Screws

1x Sinus Elevation Screw (3.5mm)

Screws are compatible with any brands 4.0mm Torque Wrench / Ratchet Wrench / Adapter's

Kit Advantages

Screws are Laser marked to indicate the depth of 8mm, 10mm & 13mm



SIMZ-YM6-318

Hex Drivers Adapter Kit

The Kit Includes:

4 Pcs Hex Drivers Kit Fits Ratchet / Torque Wrench

2 x 25mm & 2x 2.5mm Long & Short Driver

Manufactured from High Grade Stainless Steel - Highly Polished



SIMZ-YM6-319

Bone Expander Kit

The Kit Includes:

1) Expander Screws

1.6 / 2.6 mm, 1.9 / 3.0 mm, 2.1 / 3.4 mm

2.6 / 3.8 mm, 2.9 / 4.3 mm

With Laser Depth Markings on 8mm, 10mm & 11.5mm.

2) Saw Disks

7.0 mm, 10.0 mm

3) Pilot Drill 1.6mm

4) Hand piece Adapter

5) Ratchet Extender

6) Hand Driver

7) Ratchet Wrench 4.0mm

8) Sterilization Cassette / Box.



SIMZ-YM6-320

Bone Compact Kit

The Kit Includes:

1) Twist Drills

2.0mm, 2.5mm, 2.8mm 3.2mm,

3.65mm, 4.2mm, 5.2mm

With Laser Depth Markings on 6.0mm, 8.0mm, 10.0mm,

11.5mm, 13.0mm & 16.0mm

2) Ratchet / Wrench Hex Drivers

(Color Coded)

1.25mm Short / Long, 2.42mm Short / Long

3) Contra Angle Hex Drivers

1.25mm, 2.42mm

4) Parallel Pins Long & Short

5) Drill Extender

6) Marking Drill

7) Ratchet Wrench 6.35mm

8) Sterilization Cassette / Box / Organizer



SIMZ-YM6-321

Trimmer Kit

The Kit Includes:

- 1) Trephine Drills**
3.0mm, 4.0mm
5.0mm, 6.0mm
7.0mm

- 2) Tissue Punch**

3.5mm, 4.0mm
4.5mm, 5.0mm
5.5mm

- 3) Guide Punch**

4.0mm, 5.0mm

- 4) Saw Disks**

7.0mm, 10.0mm

- 5) Sterilization Cassette.**



SIMZ-YM6-322

Bone Compression Kit

The Kit Includes:


- Convex Compression Screws**
1.6/2.7mm, 2.0/3.3mm
2.4/3.6mm, 3.1/4.1mm

- Cancave Compression Screws**

3.0mm, 3.5mm, 4.0mm

- Spade Marking Drill**

Motor Mount Drivers Long & Short
Hand Driver
Ratchet Adapter
Ratchet Wrench



SIMZ-YM6-323

Surgical Drill Kit

The Kit Includes:

- 1) Twist Drills**

2.0mm, 2.8mm
3.2mm, 3.65mm
4.2mm, 4.8mm, 5.2mm

With Laser Depth Markings on 6.0mm,
8.0mm, 10.0mm, 11.5mm, 13.0mm & 16.0mm

- 2) Ratchet / Wrench Hex Drivers (Color Coded)**

1.25mm Short / Long, 2.42mm Short / Long

- 3) Parallel Pins**

Long & Short

- 4) Motor Mount Hex**

1.25mm, 2.5mm

- 5) Ratchet Wrench 4mm**

- 6) Sterilization Cassette / Box / Organizer**

- Инструменты для имплантологии
- Implantologie
- Implantologie

- Инструменты для имплантолѳи
- Implantoloia
- زراعة الاسنان



**SIMZ-YM6-324
Universal Kit 50 PCS**

The Kit Includes:

**16pcs Conical Drills
with Stoppers**

- 1.8mm / 2.2mm - S1, 1.8mm / 2.2mm - S3
- 1.8mm / 2.2mm - S5, 1.8mm / 2.2mm - S7
- 2.2mm / 2.7mm - S1, 2.2mm / 2.7mm - S3
- 2.2mm / 2.7mm - S5, 2.2mm / 2.7mm - S7
- 2.3mm / 3.3mm - S1, 2.3mm / 3.3mm - S3
- 2.3mm / 3.3mm - S5, 2.3mm / 3.3mm - S7
- 2.8mm / 4.5mm - S1, 2.8mm / 4.5mm - S3
- 2.8mm / 4.5mm - S5, 2.8mm / 4.5mm - S7

S = Stopper Length

Bone Expander Screws

- 1.6mm / 2.6mm, 1.9mm / 3.0mm
- 2.9mm / 3.4mm, 2.3mm / 3.8mm
- 2.9mm / 4.3mm,

Titanium Coated Drills

- Drill Bit 2.0mm L 18.50, Drill Bit 2.8mm L 18.50
- Drill Bit 3.65mm L 18.50, Drill Bit 4.2mm L 18.50
- Drill Bit 4.5mm L 18.50, Drill Bit 4.8mm L 18.50
- Drill Bit 5.2mm L 18.50

8pcs Trephine Drills

- Int: 2.0mm, Ext: 3.0mm, Int: 3.0mm, Ext: 4.0mm
- Int: 4.0mm, Ext: 5.0mm, Int: 5.0mm, Ext: 6.0mm
- Int: 6.0mm, Ext: 7.0mm, Int: 7.0mm, Ext: 8.0mm
- Int: 8.0mm, Ext: 9.0mm, Int: 9.0mm, Ext: 10.0mm
- Compression Screw 2.8 / 4.0mm, Compression Screw 4.0mm
- Driver 2.5mm
- Compression Screw 2.0 / 2.8mm, Compression Screw 2.8 / 3.4mm
- Tissue Punch 4.5mm, Tissue Punch 5.8mm
- Pilot Drill 2.0mm
- Drill Extender
- Ratchet Driver 2.5mm, Hand Driver 2.5mm
- Torque Wrench 6.35mm (10-40Ncm)
- Depth Gauge
- Sterilization Cassette / Box



**SIMZ-YM6-325
Conical Implant Drill Kit**

The Kit Includes:

Drills Sizes:

- 1.8mm / 2.2mm, 2.2mm / 2.7mm
- 2.3mm / 3.3mm, 2.7mm / 3.7mm
- 2.7mm / 4.5mm


Stopper Sizes:

- 1.0mm / L 13.0, 2.0mm / L 11.5
- 3.0mm / L 10.0, 5.0mm / L 8.0
- 7.0mm / L 6.0

Overall Drills Length:

- 13.0mm (32mm with Shaft)

BONE INJECTORS, COLLECTORS & SUCTION TUBE



1.5mm



SIMZ-YM6-339
Suction
Tubes
180mm

3mm



SIMZ-YM6-340
Suction
Tubes
180mm

5mm



SIMZ-YM6-341
Suction
Tubes
180mm



SIMZ-YM6-342
Fig. 6FR
Frazier
150mm - Ø 2 mm

SIMZ-YM6-343
Ø 3 mm
length 180mm
SIMZ-YM6-344
Ø 4 mm
length 180mm
SIMZ-YM6-345
Ø 5 mm
length 180mm



SIMZ-YM6-346
Suction Cannula Kit
with 3 Needles



Suction Pipe



Suction Needle
1mm



Suction Needle
1.5mm




Suction Needle
2mm



- Оттисные ложки
- Porte-empreinte
- Abdrucklöffel
- Відбиткові ложки
- Cubetas para impresiones
- صواني الانطباع

PERFORATED RIM LOCK UPPER



All the above impression trays are available in these 5 sizes

UPPER		U5	U4	U3	U2	U1
		XS	S	M	L	XL
A	A	52	53	54	60	62
	B	68	69	73	75	85

LOWER		L5	L4	L3	L2	L1
		XS	S	M	L	XL
A	A	52	53	54	60	62
	B	68	69	73	75	85

All Dimension are in mm

SOLID RIM LOCK UPPER




IMPRESSION TRAYS

- Оттисные ложки
- Відбиткові ложки
- Porte-empreinte
- Cubetas para impresiones
- Abdrucklöffel
- صواني الاطباع

PERFORATED FULL DENTURE UPPER



LOWER




All the above impression trays are available in these 4 sizes

UPPER	A	S	M	L	XL	LOWER	A	S	M	L	XL
		B	60	64	67			70	B	57	60
		70	73	76	79			70	70	73	86

All Dimension are in mm

SOLID FULL DENTURE UPPER



PERFORATED FULL DENTURE UPPER

- | | |
|----------------------|----------------------|
| SET OF 06 PCS | SET OF 08 PCS |
| Upper: S, M, L | Upper: S, M, L, XL |
| Lower: S, M, L | Lower: S, M, L, XL |
| SIMZ-YM7-037 | SIMZ-YM7-038 |


SOLID FULL DENTURE UPPER

- | | |
|----------------------|----------------------|
| SET OF 06 PCS | SET OF 08 PCS |
| Upper: S, M, L | Upper: S, M, L, XL |
| Lower: S, M, L | Lower: S, M, L, XL |
| SIMZ-YM7-039 | SIMZ-YM7-040 |

- Оттисные ложки
- Відбиткові ложки
- Porte-empreinte
- Cubetas para impresiones
- صواني الانطباع
- Abdrucklöffel

SOLID FULL DENTURE

UPPER



LOWER




SET OF 06 PCS

Upper: S, M, L
Lower: S, M, L
SIMZ-YM7-049

SET OF 08 PCS

Upper: S, M, L, XL
Lower: S, M, L, XL
SIMZ-YM7-050

All the above impression trays are available in these 4 sizes



UPPER		S	M	L	XL
	A	60	64	67	70
B	70	73	76	79	

LOWER		S	S	M	L
	A	57	60	64	64
B	70	70	73	86	





All Dimension are in mm

PERFORATED FULL DENTURE

UPPER



LOWER



SET OF 06 PCS

Upper: S, M, L
Lower: S, M, L
SIMZ-YM7-059

SET OF 08 PCS

Upper: S, M, L, XL
Lower: S, M, L, XL
SIMZ-YM7-060

IMPRESSION TRAYS

- Оттисные ложки
- Відбиткові ложки
- Porte-empreinte
- Cubetas para impresiones
- Abdrucklöffel
- صواني الاطباع

UPPER




LOWER




SET OF 06 PCS

Upper: S, M, L
Lower: S, M, L
SIMZ-YM7-067

All the above impression trays are available in these 3 sizes



UPPER		U3	U2	U1
		S	M	L
	A	57	64	67
	B	70	71	81



LOWER		L3	L2	L1
		S	M	L
	A	57	57	60
	B	76	79	81

All Dimension are in mm

UPPER



LOWER




SET OF 06 PCS

Upper: S, M, L
Lower: S, M, L
SIMZ-YM7-074

- Оттисковые ложки
- Видбиткові ложки
- Porte-empreinte
- Cubetas para impresiones
- Abdrucklöffel
- صواني الانطباع

English Pattern Mesh Type



SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
SIMZ-YM7-083


SET OF 08 PCS
Upper: S, M, L, XL
Lower: S, M, L, XL
SIMZ-YM7-084

All the above impression trays are available in these 4 sizes

UPPER	U 1	U 2	U 3	U 4	LOWER	L 1	L 2	L 3	L 4
	S	M	L	XL		S	M	L	XL
A	12	13	15	15	A	14	15	15	17
B	50	51	55	58	B	50	51	55	60
C	67	72	78	80	C	70	77	79	80

All Dimension are in mm

English Pattern Punched Type



SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
SIMZ-YM7-093

SET OF 08 PCS
Upper: S, M, L, XL
Lower: S, M, L, XL
SIMZ-YM7-094

IMPRESSION TRAYS

- Оттисные ложки
- Відбиткові ложки
- Porte-empreinte
- Cubetas para impresiones
- Abdrucklöffel
- صواني الاطباع

ENGLISH PATTERN SOLID TYPE

UPPER




LOWER



SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
SIMZ-YM7-103


SET OF 08 PCS
Upper: S, M, L, XL
Lower: S, M, L, XL
SIMZ-YM7-104

All the above impression trays are available in these 4 sizes



UPPER		U 1	U 2	U 3	U 4
		S	M	L	XL
A		12	13	15	15
B		50	51	55	58
C		67	72	78	80

LOWER		L 1	L 2	L 3	L 4
		S	M	L	XL
A		14	15	15	17
B		50	51	55	60
C		70	77	79	80



All Dimension are in mm

- Оттисные ложки
- Porte-empreinte
- Abdrucklöffel
- Відбиткові ложки
- Cubetas para impresiones
- صواني الانطباع

SOLID FULL DENTURE

UPPER



LOWER



SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
SIMZ-YM7-113

SET OF 08 PCS
Upper: S, M, L, XL
Lower: S, M, L, XL
SIMZ-YM7-114

PERFORATED FULL DENTURE

UPPER



LOWER



SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
SIMZ-YM7-123

SET OF 08 PCS
Upper: S, M, L, XL
Lower: S, M, L, XL
SIMZ-YM7-124

IMPRESSION TRAYS

- Оттисковые ложки
- Відбиткові ложки
- Porte-empreinte
- Cubetas para impresiones
- Abdrucklöffel
- صواني الاطباع

ORTHODONTIAL SOLID

UPPER



S
SIMZ-YM7-125



M
SIMZ-YM7-126



L
SIMZ-YM7-127



XL
SIMZ-YM7-128

LOWER



S
SIMZ-YM7-129



M
SIMZ-YM7-130



L
SIMZ-YM7-131



XL
SIMZ-YM7-132

SET OF 06 PCS

Upper: S, M, L
Lower: S, M, L
SIMZ-YM7-133

SET OF 08 PCS

Upper: S, M, L, XL
Lower: S, M, L, XL
SIMZ-YM7-134

ORTHODONTIAL PERFORATED

UPPER



S
SIMZ-YM7-135



M
SIMZ-YM7-136



L
SIMZ-YM7-137



XL
SIMZ-YM7-138

LOWER



S
SIMZ-YM7-139



M
SIMZ-YM7-140



L
SIMZ-YM7-141



XL
SIMZ-YM7-142

SET OF 06 PCS

Upper: S, M, L
Lower: S, M, L
SIMZ-YM7-143


SET OF 08 PCS

Upper: S, M, L, XL
Lower: S, M, L, XL
SIMZ-YM7-144


- Оттисные ложки
 - Porte-empreinte
 - Abdrucklöffel
- Відбиткові ложки
 - Cubetas para impresiones
 - صواني الانطباع

PERFORATED WITHOUT RIM LOCK

UPPER



LOWER



SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
SIMZ-YM7-157

SET OF 08 PCS
Upper: XS, S, M, L
Lower: XS, S, M, L
SIMZ-YM7-158

SET OF 08 PCS
Upper: S, M, L, XL
Lower: S, M, L, XL
SIMZ-YM7-159

SET OF 10 PCS
Upper: XS, S, M, L, XL
Lower: XS, S, M, L, XL
SIMZ-YM7-160

SET OF 12 PCS
Upper: XXS, XS, S, M, L, XL
Lower: XXS, XS, S, M, L, XL
SIMZ-YM7-161

All the above impression trays are available in these 6 sizes

UPPER		U1	U2	U3	U4	U5	U6
		XXS	XS	S	M	L	XL
	A	61	64	66	69	71	73
B	47	51	53	57	60	62	

LOWER		L1	L2	L3	L4	L5	L6
		XXS	XS	S	M	L	XL
	A	55	60	72	74	76	78
B	47	49	54	57	59	61	


All Dimension are in mm

IMPRESSION TRAYS

- Оттисковые ложки
- Porte-empreinte
- Abdrucklöffel
- Відбиткові ложки
- Cubetas para impresiones
- صواني الانطباع

SOLID WITHOUT RIM LOCK

UPPER



LOWER



SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
SIMZ-YM7-174

SET OF 08 PCS
Upper: XS, S, M, L
Lower: XS, S, M, L
SIMZ-YM7-175

SET OF 08 PCS
Upper: S, M, L, XL
Lower: S, M, L, XL
SIMZ-YM7-176

SET OF 10 PCS
Upper: XS, S, M, L, XL
Lower: XS, S, M, L, XL
SIMZ-YM7-177

SET OF 12 PCS
Upper: XXS, XS, S, M, L, XL
Lower: XXS, XS, S, M, L, XL
SIMZ-YM7-178

All the above impression trays are available in these 6 sizes

UPPER	U1	U2	U3	U4	U5	U6
	XXS	XS	S	M	L	XL
A	61	64	66	69	71	73
B	47	51	53	57	60	62

LOWER	L1	L2	L3	L4	L5	L6
	XXS	XS	S	M	L	XL
A	55	60	72	74	76	78
B	47	49	54	57	59	61

All Dimension are in mm

- Оттисные ложки
- Porte-empreinte
- Abdrucklöffel
- Відбиткові ложки
- Cubetas para impresiones
- صواني الانطباع

PARTIAL SOLID



SIMZ-YM7-179 44mm
SIMZ-YM7-180 55mm
SIMZ-YM7-181 60mm
SIMZ-YM7-182 65mm



SIMZ-YM7-183



-YM7-184 44mm
-YM7-185 55mm
-YM7-186 60mm
-YM7-187 65mm



SIMZ-YM7-188



SIMZ-YM7-189

Set of 5 Pcs SIMZ-YM7-190

PARTIAL PERFORATED



SIMZ-YM7-191 44mm
SIMZ-YM7-192 55mm
SIMZ-YM7-193 60mm
SIMZ-YM7-194 65mm



SIMZ-YM7-195



SIMZ-YM7-196 44mm
SIMZ-YM7-197 55mm
SIMZ-YM7-198 60mm
SIMZ-YM7-199 65mm



SIMZ-YM7-200



SIMZ-YM7-201



SIMZ-YM7-202




SIMZ-YM7-203

Set of 7 Pcs. SIMZ-YM7-204

IMPRESSION TRAYS

- Оттисные ложки
- Відбиткові ложки
- Porte-empreinte
- Cubetas para impresiones
- Abdrucklöffel
- صواني الاعطايح



Perforated

1/2L

1/2R

SIMZ-YM7-205 1/2L/S
SIMZ-YM7-206 1/2L/M
SIMZ-YM7-207 1/2L/L

SIMZ-YM7-208 1/2R/S
SIMZ-YM7-209 1/2R/M
SIMZ-YM7-210 1/2R/L



non-Perforated

1/2L

1/2R

SIMZ-YM7-211 1/2L/S
SIMZ-YM7-212 1/2L/M
SIMZ-YM7-213 1/2L/L

SIMZ-YM7-214 1/2R/S
SIMZ-YM7-215 1/2R/M
SIMZ-YM7-216 1/2R/L



Universal

Back

Universal

SIMZ-YM7-217 UN/1
Front

SIMZ-YM7-218 UN/1
Front

with swivel handle and fixing screw



Back

SIMZ-YM7-219 1
Front
with swivel handle

SIMZ-YM7-220 1
Front



Stolley

SIMZ-YM7-221
ST/1

Partial occlusal stainless steel impression trays

They allow you to take an impression of two semi-arches antagonistic



with swivel handle and solid

SIMZ-YM7-222
OCP



SIMZ-YM7-223
OCG

IMPRESSION TRAYS

- Оттисные ложки
- Porte-empreinte
- Abdrucklöffel
- Відбиткові ложки
- Cubetas para impresiones
- صواني الانطباع



	A mm.	B mm.
U 0	59	52
U 1	66	52
U 2	71	62
U 3	71	64
U 4	80	70

UPPER



	A mm.	B mm.
L 0	65	46
L 1	67	50
L 2	76	54
L 3	79	56
L 4	79	61

LOWER

SET OF 10 PCS
Upper: XS, S, M, L, XL
Lower: XS, S, M, L, XL
SIMZ-YM7-234

IMPRESSION TRAYS

- Оттисные ложки
- Porte-empreinte
- Abdrucklöffel
- Відбиткові ложки
- Cubetas para impresiones
- صواني الاطباع

Ehricke Stainless Steel Impression Trays

UPPER



S
SIMZ-YM7-235
U1



M
SIMZ-YM7-236
U2



L
SIMZ-YM7-237
U3

LOWER



S
SIMZ-YM7-238
L1



M
SIMZ-YM7-239
L2



L
SIMZ-YM7-240
L3

SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
SIMZ-YM7-241

	A mm.	B mm.	UPPER
U 1	64	46	
U 2	68	53	
U 3	79	61	

	A mm.	B mm.	LOWER
L 1	70	48	
L 2	70	53	
L 3	72	57	

- Оттисные ложки
 - Porte-empainte
 - Abdrucklöffel
- Відбиткові ложки
 - Cubetas para impresiones
 - صواني الانطباع

UPPER



S
SIMZ-YM7-242
U1



M
SIMZ-YM7-243
U2



L
SIMZ-YM7-244
U3

LOWER



S
SIMZ-YM7-245
L1



M
SIMZ-YM7-246
L2



L
SIMZ-YM7-247
L3

SET OF 06 PCS

Upper: S, M, L
Lower: S, M, L
SIMZ-YM7-248

	A mm.	B mm.	UPPER
U 1	62	50	
U 2	70	56	
U 3	78	62	

	A mm.	B mm.	LOWER
L 1	70	55	
L 2	78	59	
L 3	86	63	






Scalers Tray
SIMZ-YM8-001 200 x 100 mm



Scalers Trays
SIMZ-YM8-002 205 x 115 x 17 mm
SIMZ-YM8-003 205 x 152 x 17 mm
SIMZ-YM8-004 290 x 115 x 17 mm



Instruments Trays
SIMZ-YM8-005 200 x 110 x 30 mm
SIMZ-YM8-006 215 x 115 x 35 mm
SIMZ-YM8-007 220 x 130 x 30 mm
SIMZ-YM8-008 230 x 125 x 35 mm
SIMZ-YM8-009 245 x 140 x 40 mm
SIMZ-YM8-010 285 x 160 x 45 mm
SIMZ-YM8-011 315 x 175 x 50 mm
SIMZ-YM8-012 330 x 185 x 50 mm



Instruments Trays
SIMZ-YM8-013 150 x 100 x 50 mm
SIMZ-YM8-014 200 x 100 x 50 mm
SIMZ-YM8-015 220 x 120 x 50 mm
SIMZ-YM8-016 240 x 140 x 50 mm
SIMZ-YM8-017 200 x 150 x 50 mm
SIMZ-YM8-018 250 x 200 x 50 mm
SIMZ-YM8-019 300 x 250 x 50 mm
SIMZ-YM8-020 300 x 200 x 50 mm
SIMZ-YM8-021 350 x 250 x 50 mm
SIMZ-YM8-022 400 x 300 x 50 mm
SIMZ-YM8-023 450 x 300 x 50 mm
SIMZ-YM8-024 450 x 350 x 50 mm



Instruments Trays With Lid
SIMZ-YM8-025 150 x 100 x 50 mm
SIMZ-YM8-026 200 x 100 x 50 mm
SIMZ-YM8-027 220 x 120 x 50 mm
SIMZ-YM8-028 240 x 140 x 50 mm
SIMZ-YM8-029 200 x 150 x 50 mm
SIMZ-YM8-030 250 x 200 x 50 mm
SIMZ-YM8-031 300 x 250 x 50 mm
SIMZ-YM8-032 300 x 200 x 50 mm
SIMZ-YM8-033 350 x 250 x 50 mm
SIMZ-YM8-034 400 x 300 x 50 mm
SIMZ-YM8-035 450 x 300 x 50 mm
SIMZ-YM8-036 450 x 350 x 50 mm



Mayo Tray Perforated
SIMZ-YM8-037 150 x 100 x 25 mm
SIMZ-YM8-038 200 x 100 x 25 mm
SIMZ-YM8-039 220 x 120 x 25 mm
SIMZ-YM8-040 240 x 140 x 25 mm
SIMZ-YM8-041 200 x 150 x 25 mm
SIMZ-YM8-042 250 x 200 x 25 mm
SIMZ-YM8-043 300 x 250 x 25 mm
SIMZ-YM8-044 300 x 200 x 25 mm
SIMZ-YM8-045 350 x 250 x 25 mm
SIMZ-YM8-046 400 x 300 x 25 mm
SIMZ-YM8-047 450 x 300 x 25 mm
SIMZ-YM8-048 450 x 350 x 25 mm



Scaler Tray
SIMZ-YM8-049 200 x 115 mm
SIMZ-YM8-050 205 x 150 mm
SIMZ-YM8-051 290 x 150 mm
SIMZ-YM8-052 300 x 200 mm



Instruments Tray Perforated
SIMZ-YM8-053 150 x 100 x 50 mm
SIMZ-YM8-054 200 x 100 x 50 mm
SIMZ-YM8-055 220 x 120 x 50 mm
SIMZ-YM8-056 240 x 140 x 50 mm
SIMZ-YM8-057 200 x 150 x 50 mm
SIMZ-YM8-058 250 x 200 x 50 mm
SIMZ-YM8-059 300 x 250 x 50 mm
SIMZ-YM8-060 300 x 200 x 50 mm
SIMZ-YM8-061 350 x 250 x 50 mm
SIMZ-YM8-062 400 x 300 x 50 mm
SIMZ-YM8-063 450 x 300 x 50 mm
SIMZ-YM8-064 450 x 350 x 50 mm

Sterilizing boxes & instruments Trays



- Instruments Box with Hooks**
- SIMZ-YM8-065 170 x 70 x 30 mm
 - SIMZ-YM8-066 180 x 80 x 40 mm
 - SIMZ-YM8-067 200 x 100 x 40 mm
 - SIMZ-YM8-068 200 x 100 x 50 mm
 - SIMZ-YM8-069 220 x 120 x 50 mm
 - SIMZ-YM8-070 250 x 150 x 50 mm
 - SIMZ-YM8-071 250 x 200 x 50 mm
 - SIMZ-YM8-072 300 x 150 x 60 mm
 - SIMZ-YM8-073 300 x 200 x 50 mm
 - SIMZ-YM8-074 350 x 250 x 50 mm
 - SIMZ-YM8-075 400 x 200 x 100 mm
 - SIMZ-YM8-076 450 x 250 x 50 mm



- Instruments Box**
- SIMZ-YM8-077 170 x 70 x 30 mm
 - SIMZ-YM8-078 180 x 80 x 40 mm
 - SIMZ-YM8-079 200 x 100 x 40 mm
 - SIMZ-YM8-080 200 x 100 x 50 mm
 - SIMZ-YM8-081 220 x 120 x 50 mm
 - SIMZ-YM8-082 250 x 150 x 50 mm
 - SIMZ-YM8-083 250 x 200 x 50 mm
 - SIMZ-YM8-084 300 x 150 x 60 mm
 - SIMZ-YM8-085 300 x 200 x 50 mm
 - SIMZ-YM8-086 350 x 250 x 50 mm
 - SIMZ-YM8-087 400 x 200 x 100 mm
 - SIMZ-YM8-088 450 x 250 x 50 mm



- Instruments Box Perforated**
- SIMZ-YM8-089 170 x 70 x 30 mm
 - SIMZ-YM8-090 180 x 80 x 40 mm
 - SIMZ-YM8-091 200 x 100 x 40 mm
 - SIMZ-YM8-092 200 x 100 x 50 mm
 - SIMZ-YM8-093 220 x 120 x 50 mm
 - SIMZ-YM8-094 250 x 150 x 50 mm
 - SIMZ-YM8-095 250 x 200 x 50 mm
 - SIMZ-YM8-096 300 x 150 x 60 mm
 - SIMZ-YM8-097 300 x 200 x 50 mm
 - SIMZ-YM8-098 350 x 250 x 50 mm
 - SIMZ-YM8-099 400 x 200 x 100 mm
 - SIMZ-YM8-100 450 x 250 x 50 mm



- Gloves Box**
SIMZ-YM8-101 300 x 150 mm



- Catheters Tray with lid**
- SIMZ-YM8-102 200 x 60 x 40 mm
 - SIMZ-YM8-103 200 x 90 x 40 mm
 - SIMZ-YM8-104 320 x 95 x 50 mm
 - SIMZ-YM8-105 328 x 90 x 50 mm
 - SIMZ-YM8-106 365 x 90 x 50 mm
 - SIMZ-YM8-107 440 x 115 x 50 mm



- Box With Filter**
- SIMZ-YM8-108 100 x 50 mm
 - SIMZ-YM8-109 250 x 150 mm
 - SIMZ-YM8-110 150 x 65 mm
 - SIMZ-YM8-111 200 x 100 mm
 - SIMZ-YM8-112 300 x 200 mm



- Instruments Trays With Cover**
- SIMZ-YM8-113 210 x 110 x 50 mm
 - SIMZ-YM8-114 240 x 160 x 50 mm
 - SIMZ-YM8-115 300 x 180 x 50 mm
 - SIMZ-YM8-116 360 x 210 x 50 mm
 - SIMZ-YM8-117 400 x 230 x 55 mm

- Instruments Trays Without Cover**
- SIMZ-YM8-118 210 x 110 x 50 mm
 - SIMZ-YM8-119 240 x 160 x 50 mm
 - SIMZ-YM8-120 300 x 180 x 50 mm
 - SIMZ-YM8-121 360 x 210 x 50 mm
 - SIMZ-YM8-122 400 x 230 x 55 mm



- Instruments Box with lid**
- SIMZ-YM8-123 170 x 70 x 30 mm
 - SIMZ-YM8-124 180 x 80 x 40 mm
 - SIMZ-YM8-125 200 x 100 x 40 mm
 - SIMZ-YM8-126 200 x 100 x 50 mm
 - SIMZ-YM8-127 220 x 120 x 50 mm
 - SIMZ-YM8-128 250 x 150 x 50 mm
 - SIMZ-YM8-129 250 x 200 x 50 mm
 - SIMZ-YM8-130 300 x 150 x 60 mm
 - SIMZ-YM8-131 300 x 200 x 50 mm
 - SIMZ-YM8-132 350 x 250 x 50 mm
 - SIMZ-YM8-133 400 x 200 x 100 mm
 - SIMZ-YM8-134 450 x 250 x 50 mm



Kidney Tray
 SIMZ-YM8-135 160 mm
 SIMZ-YM8-136 200 mm
 SIMZ-YM8-137 250 mm
 SIMZ-YM8-138 280 mm
 SIMZ-YM8-139 300 mm



Kidney Tray Shallow
 SIMZ-YM8-140 170 mm
 SIMZ-YM8-141 220 mm
 SIMZ-YM8-142 240 mm



Kidney Tray With Cover
 SIMZ-YM8-143 160 mm
 SIMZ-YM8-144 200 mm
 SIMZ-YM8-145 250 mm
 SIMZ-YM8-146 280 mm
 SIMZ-YM8-147 300 mm



Kidney Tray
 SIMZ-YM8-148 160 mm
 SIMZ-YM8-149 200 mm
 SIMZ-YM8-150 250 mm
 SIMZ-YM8-151 280 mm
 SIMZ-YM8-152 300 mm



Forceps Jars
 SIMZ-YM8-153 30 x 90 mm
 SIMZ-YM8-154 50 x 100 mm
 SIMZ-YM8-155 50 x 130 mm
 SIMZ-YM8-156 50 x 175 mm
 SIMZ-YM8-157 75 x 100 mm
 SIMZ-YM8-158 75 x 130 mm
 SIMZ-YM8-159 75 x 175 mm
 SIMZ-YM8-160 100 x 175 mm
 SIMZ-YM8-161 100 x 200 mm



Methot Jar Without Lid
 SIMZ-YM8-162



Dressing Jar / Cotton Container
 SIMZ-YM8-163 80 x 80 mm
 SIMZ-YM8-164 100 x 60 mm
 SIMZ-YM8-165 100 x 80 mm
 SIMZ-YM8-166 100 x 100 mm
 SIMZ-YM8-167 120 x 80 mm
 SIMZ-YM8-168 120 x 100 mm
 SIMZ-YM8-169 150 x 100 mm
 SIMZ-YM8-170 150 x 150 mm
 SIMZ-YM8-171 180 x 150 mm
 SIMZ-YM8-172 180 x 180 mm
 SIMZ-YM8-173 200 x 200 mm
 SIMZ-YM8-174 240 x 240 mm



Sputum Cup
 SIMZ-YM8-175



Dressing Jar With Lid & Knob
 SIMZ-YM8-176 80 x 80 mm
 SIMZ-YM8-177 90 x 90 mm
 SIMZ-YM8-178 100 x 100 mm
 SIMZ-YM8-179 120 x 120 mm
 SIMZ-YM8-180 100 x 130 mm
 SIMZ-YM8-181 150 x 150 mm
 SIMZ-YM8-182 125 x 170 mm



Round Bowls
 SIMZ-YM8-183 80 x 35 mm
 SIMZ-YM8-184 125 x 60 mm



**Spirit Lamp
SIMZ-YM8-185**



**Spirit Lamp
SIMZ-YM8-186**



**Spirit Lamp
SIMZ-YM8-187**



**Cotton Holding Jar
SIMZ-YM8-188**



**Cotton Holder
SIMZ-YM8-189**



**Cotton Holder
SIMZ-YM8-190**



**Cotton Wool Dispenser Acc.
to Methol
SIMZ-YM8-191**



**Needle Sterilization Tub
SIMZ-YM8-192**



**Round Endodontic Boxes
complete with lid
SIMZ-YM8-193 Ø mm. 94 - h. mm. 70**



**Endodontic Box
SIMZ-YM8-194**



SIMZ-YM8-195
100 x 180 x 11mm
Mini size



SIMZ-YM8-196
100 x 180 x 11mm
Mini size



SIMZ-YM8-197
100 x 180 x 11mm
Mini size



SIMZ-YM8-198
100 x 180 x 11mm
Mini size



SIMZ-YM8-199
100 x 180 x 11mm
Mini size



Lid
Coperchio

SIMZ-YM8-200
103 x 182 x 10mm
Mini size



SIMZ-YM8-201
Silver



SIMZ-YM8-202 Blue
SIMZ-YM8-203 Red



SIMZ-YM8-204 Blue
SIMZ-YM8-205 Red



SIMZ-YM8-206
Silver



SIMZ-YM8-207
Green



SIMZ-YM8-209
Blue



SIMZ-YM8-208
Golden



SIMZ-YM8-210
Red




SIMZ-YM8-211
Silver



SIMZ-YM8-212
Midi
125 x 105 x 54mm



SIMZ-YM8-213
Maxi
204 x 105 x 54mm



SIMZ-YM8-214
Maxi
204 x 105 x 54mm



Round Sterilizer Drum

- SIMZ-YM8-215 mm. Ø 80 x 80 H
- SIMZ-YM8-216 mm. Ø 125 x 80 H
- SIMZ-YM8-217 mm. Ø 125 x 100 H
- SIMZ-YM8-218 mm. Ø 125 x 125 H
- SIMZ-YM8-219 mm. Ø 150 x 80 H
- SIMZ-YM8-220 mm. Ø 150 x 100 H
- SIMZ-YM8-221 mm. Ø 150 x 130 H
- SIMZ-YM8-222 mm. Ø 150 x 150 H
- SIMZ-YM8-223 mm. Ø 150 x 190 H
- SIMZ-YM8-224 mm. Ø 190 x 125 H
- SIMZ-YM8-225 mm. Ø 190 x 145 H
- SIMZ-YM8-226 mm. Ø 190 x 165 H
- SIMZ-YM8-227 mm. Ø 190 x 190 H
- SIMZ-YM8-228 mm. Ø 190 x 240 H
- SIMZ-YM8-229 mm. Ø 190 x 290 H
- SIMZ-YM8-230 mm. Ø 240 x 120 H
- SIMZ-YM8-231 mm. Ø 240 x 145 H
- SIMZ-YM8-232 mm. Ø 240 x 165 H
- SIMZ-YM8-233 mm. Ø 240 x 190 H
- SIMZ-YM8-234 mm. Ø 240 x 240 H
- SIMZ-YM8-235 mm. Ø 240 x 290 H
- SIMZ-YM8-236 mm. Ø 240 x 340 H
- SIMZ-YM8-237 mm. Ø 240 x 390 H
- SIMZ-YM8-238 mm. Ø 290 x 145 H
- SIMZ-YM8-239 mm. Ø 290 x 165 H
- SIMZ-YM8-240 mm. Ø 290 x 190 H
- SIMZ-YM8-241 mm. Ø 290 x 240 H
- SIMZ-YM8-242 mm. Ø 290 x 290 H
- SIMZ-YM8-243 mm. Ø 290 x 340 H
- SIMZ-YM8-244 mm. Ø 290 x 390 H
- SIMZ-YM8-245 mm. Ø 340 x 145 H
- SIMZ-YM8-246 mm. Ø 340 x 165 H
- SIMZ-YM8-247 mm. Ø 340 x 190 H
- SIMZ-YM8-248 mm. Ø 340 x 240 H
- SIMZ-YM8-249 mm. Ø 340 x 290 H
- SIMZ-YM8-250 mm. Ø 340 x 340 H
- SIMZ-YM8-251 mm. Ø 340 x 390 H
- SIMZ-YM8-252 mm. Ø 390 x 140 H
- SIMZ-YM8-253 mm. Ø 390 x 165 H
- SIMZ-YM8-254 mm. Ø 390 x 190 H
- SIMZ-YM8-255 mm. Ø 390 x 240 H
- SIMZ-YM8-256 mm. Ø 390 x 290 H
- SIMZ-YM8-257 mm. Ø 390 x 340 H
- SIMZ-YM8-258 mm. Ø 390 x 390 H
- SIMZ-YM8-259 mm. Ø 390 x 490 H



Round Sterilizer Drums "Filter"

- SIMZ-YM8-260 mm. Ø 80 x 80 H
- SIMZ-YM8-261 mm. Ø 125 x 80 H
- SIMZ-YM8-262 mm. Ø 125 x 100 H
- SIMZ-YM8-263 mm. Ø 125 x 125 H
- SIMZ-YM8-264 mm. Ø 150 x 100 H
- SIMZ-YM8-265 mm. Ø 150 x 130 H
- SIMZ-YM8-266 mm. Ø 150 x 150 H
- SIMZ-YM8-267 mm. Ø 150 x 190 H
- SIMZ-YM8-268 mm. Ø 190 x 145 H
- SIMZ-YM8-269 mm. Ø 190 x 165 H
- SIMZ-YM8-270 mm. Ø 190 x 190 H
- SIMZ-YM8-271 mm. Ø 190 x 240 H
- SIMZ-YM8-272 mm. Ø 240 x 145 H
- SIMZ-YM8-273 mm. Ø 240 x 165 H
- SIMZ-YM8-274 mm. Ø 240 x 190 H
- SIMZ-YM8-275 mm. Ø 240 x 240 H
- SIMZ-YM8-276 mm. Ø 240 x 290 H
- SIMZ-YM8-277 mm. Ø 240 x 340 H
- SIMZ-YM8-278 mm. Ø 240 x 390 H
- SIMZ-YM8-279 mm. Ø 290 x 145 H
- SIMZ-YM8-280 mm. Ø 290 x 165 H
- SIMZ-YM8-281 mm. Ø 290 x 190 H
- SIMZ-YM8-282 mm. Ø 290 x 240 H
- SIMZ-YM8-283 mm. Ø 290 x 290 H
- SIMZ-YM8-284 mm. Ø 290 x 340 H
- SIMZ-YM8-285 mm. Ø 290 x 390 H
- SIMZ-YM8-286 mm. Ø 340 x 145 H
- SIMZ-YM8-287 mm. Ø 340 x 165 H
- SIMZ-YM8-288 mm. Ø 340 x 190 H
- SIMZ-YM8-289 mm. Ø 340 x 240 H
- SIMZ-YM8-290 mm. Ø 340 x 290 H
- SIMZ-YM8-291 mm. Ø 390 x 140 H
- SIMZ-YM8-292 mm. Ø 390 x 165 H
- SIMZ-YM8-293 mm. Ø 390 x 190 H
- SIMZ-YM8-294 mm. Ø 390 x 240 H
- SIMZ-YM8-295 mm. Ø 390 x 290 H
- SIMZ-YM8-296 mm. Ø 390 x 340 H
- SIMZ-YM8-297 mm. Ø 390 x 490 H




Round Bowls With Cover

- SIMZ-YM8-298 Ø 40 mm ml.30
- SIMZ-YM8-299 Ø 60 mm ml.50
- SIMZ-YM8-300 Ø 80 mm ml.150
- SIMZ-YM8-301 Ø 100 mm ml.250
- SIMZ-YM8-302 Ø 120 mm ml.450
- SIMZ-YM8-303 Ø 150 mm ml.750
- SIMZ-YM8-304 Ø 180 mm ml.1500
- SIMZ-YM8-305 Ø 200 mm ml.2100
- SIMZ-YM8-306 Ø 250 mm ml.3500


Round Bowls With Out Cover

- SIMZ-YM8-307 Ø 40 mm ml.30
- SIMZ-YM8-308 Ø 60 mm ml.50
- SIMZ-YM8-309 Ø 80 mm ml.150
- SIMZ-YM8-310 Ø 100 mm ml.250
- SIMZ-YM8-311 Ø 120 mm ml.450
- SIMZ-YM8-312 Ø 150 mm ml.750
- SIMZ-YM8-313 Ø 180 mm ml.1500
- SIMZ-YM8-314 Ø 200 mm ml.2100
- SIMZ-YM8-315 Ø 250 mm ml.3500

INSTRUMENTS MANAGEMENT SYSTEM



Cassette For 4 Instruments
SIMZ-YM8-316



Cassette For 5 Instruments
SIMZ-YM8-317



Cassette For 6 Instruments
SIMZ-YM8-318



Cassette For 7 Instruments
SIMZ-YM8-319



Cassette For 10 Instruments
SIMZ-YM8-320



Cassette For 20 Instruments
SIMZ-YM8-321



Cassette For Extraction Forceps
SIMZ-YM8-322



Cassette For Root Elevator
SIMZ-YM8-323



SIMZ-YM8-324
Cassette Tray
15 Pcs




SIMZ-YM8-325
Cassette Tray
15 Pcs

STERILIZING MESH TRAY

Different sizes & dimensions are available at buyer's demand



MESH TRAYS
SIMZ-YM8-326
Different sizes & dimensions are available at buyer's demand



MESH TRAYS
SIMZ-YM8-327
05 Pcs



MESH TRAYS
SIMZ-YM8-328
07 Pcs



MESH TRAYS
SIMZ-YM8-329
10 Pcs



MESH TRAYS
SIMZ-YM8-330
20 Pcs



SIMZ-YM8-331
Mesh Cassettes
5 Pieces, 7 Pieces, 10 Pieces,
15 Pieces & 20 Pieces



MESH TRAY
SIMZ-YM8-332 240x130x50mm
SIMZ-YM8-333 240x130x60mm
SIMZ-YM8-334 240x250x50mm
SIMZ-YM8-335 240x250x60mm
SIMZ-YM8-336 480x250x50mm
SIMZ-YM8-337 480x250x60mm
SIMZ-YM8-338 480x250x70mm



MESH TRAY
SIMZ-YM8-339 240x250x 50mm
SIMZ-YM8-340 240x250x 60mm
SIMZ-YM8-341 240x250x 70mm
SIMZ-YM8-342 480x250x 50mm
SIMZ-YM8-343 480x250x 60mm
SIMZ-YM8-344 480x250x 70mm
SIMZ-YM8-345 480x250x100mm



MESH TRAY
SIMZ-YM8-346 240x250x50mm
SIMZ-YM8-347 240x250x60mm
SIMZ-YM8-348 240x250x70mm
SIMZ-YM8-349 280x250x50mm
SIMZ-YM8-350 280x250x60mm



MESH TRAY
SIMZ-YM8-351 280x250x 70mm
SIMZ-YM8-352 280x250x100mm
SIMZ-YM8-353 280x250x 50mm
SIMZ-YM8-354 280x280x 60mm
SIMZ-YM8-355 280x280x 70mm
SIMZ-YM8-356 300x280x 60mm



SIMZ-YM8-357
18×7×3.5cm
5 Instruments Cassette
Detachable Lids
+
5 Instruments Cassette Hinged
Lids with button Lock



SIMZ-YM8-358
29×18.5×3.5cm
20 Instruments Cassette Hinged
Lids with surgical kit fitting,
button lock



Instruments Cassette
SIMZ-YM8-359 5 Pieces Cassette
SIMZ-YM8-360 7 Pieces Cassette
SIMZ-YM8-361 10 Pieces Cassette



SIMZ-YM8-362
18×11×3.5cm
Extra large cassette for 6
Instruments, with compartment
for small parts hinged lids



SIMZ-YM8-363
29×18.5×3.5cm
20 Instruments Cassette
Detachable Lids
+
20 Instruments Cassette Hinged
Lids with button Lock



SIMZ-YM8-364
10×13×2.2cm
10 Instruments Cassette Hinged
Lids, Click Lock Available in 2
designs



SIMZ-YM8-365
29×18.5×3.5cm
20 Instruments Cassette Detachable
Lids with elevator fitting
+
20 Instruments Cassette Hinged
Lids with elevator fitting . button lock



SIMZ-YM8-366
18×9×2.2cm
7 Instruments Cassette Hinged
Lids, Click Lock, Available in 2
designs



SIMZ-YM8-367
21×14×3.5cm
15 Instruments cassette,
Detachable lids, button lock



SIMZ-YM8-368
21×14×3.5cm
15 Instruments cassette,
detachable lids, button
lock, with fitting of micro
surgery kit



SIMZ-YM8-369
21×14×4.5cm
8 pcs luxator/elevator
cassette, Detach0
ablelids, button lock



SIMZ-YM8-370
15x9x5cm
Box for dental implant tools and
surgical kits, Capacity:8
tools, 18 tools, 24 Tools, 58 Tools
+ removable compartment



Instruments Tray Perforated
SIMZ-YM8-371 288 x 187 x 40 mm



Instruments Tray Perforated Base
SIMZ-YM8-372 288 x 187 x 40 mm



Instruments Tray Solid Base
SIMZ-YM8-373 288 x 187 x 40 mm

DESCRIPTION	PAGE#
A	
ABOU-RASS EXPLORERS	125
ACUTE IMPLANTOLOGY	307
ADAMS PLIER	275-288
ADAMS PLIERS FOR BENDING WIRE ARCHES	275
ADERER CLAMP THREE JAW WIRE BENDING PLIER	288
ADERER FACEBOW LOOP FORMING PLIERS	279
ADERER LARGE WIRE BENDING AND LOOP FORMING PLIERS	274
ADERER PLIER	274
ADHESIVE REMOVING PLIERS	283 & 286
ADSON / MICRO-ADSON TWEEZERS	118
ADSON 7X7 CURVED TISSUE FORCEPS	323
ADSON 7X7 TISSUE FORCEPS	323
ADSON BROWN TWEEZERS	117
ADSON CARBIDE TIPPED TISSUE FORCEPS	323
ADSON NEEDLE HOLDERS	242
ADSON T.C TWEEZERS	118-119
ADSON TWEEZERS	117 & 119
ADSON WITH TEETH TWEEZERS	116
ADSON-BROWN TWEEZERS	118
ADSON-TC TWEEZERS	117
ADULT BUCCAL DENTAL PHOTOGRAPHIC MIRRORS	217
ADULT PALATAL DENTAL PHOTOGRAPHIC MIRRORS	217
ALBRECHT EXTRACTING FORCEPS (FOR CHILDREN)	31
ALGINATE & PLASTER SPATULAS	211
ALLEN PERIODONTAL INSTRUMENTS	259
ALLIS STERILIZING & TISSUE CLAMPS	238
ALLIS-BABY STERILIZING & TISSUE CLAMPS	238
ALL-PURPOSE PLIERS TC GOLD	241
ALLUMINIUM CALIPER	298
AMALGAM & PLASTIC FILLING INSTRUMENTS	136-138
AMALGAM CARRIERS	255
AMALGAM GUN	336
AMALGAM INSTRUMENTS	185 & 254-255
AMALGAM INSTRUMENTS AND BURNISHERS FILLING INSTRUMENTS	139-140
AMALGAM WELL	252 & 255
AMSTERDAM TWEEZERS	114
ANATOMIC HANDLE EXPLORERS	124
ANGLE PLIER	274
ANGLED ROOT FORCEPS 45°	41
ANGLED ROOT FORCEPS 90°	41
ANISWORTH RUBBER DAM CLAMPS FORCEPS	244
ANTERIOR	148
ANTERIOR SINGLE ENDED	134
ANTKOVIK WAX KNIVES	210
ANVIL DENTAL LABORATORY INSTRUMENT	216
APEX FORCEPS	37-38
APEXO ERGO DYNAMIC ROOT ELEVATOR	86
APEXO ROOT ELEVATORS	71 & 78-79 86 & 88
APICAL RETENTION FORCEPS	66
APCAL ROOT ELEVATORS	73 & 77 & 100
APPLICATOR PERIODONTAL INSTRUMENTS	157
ARCH BENDING PLIERS	284 & 288
ARCH COUNTING PLIER	276
ARROW CLASP FORMING PLIERS	280
ARTICULATORS	214-216
ATRAUMATIC EXTRACTION BAYONET ROOT ELEVATORS	95-96
A-TRAUMATIC EXTRACTION FORCEPS	35
AUSTIN TISSUE RETRACTORS	224
B	
BABY TOOTH EXTRACTION FORCEPS	26
BACKHAUS HAEMOSTATIC FORCEPS & TOWEL CLAMPS	237
BACKHAUS NOT-PERFORATING AND NOT-DAMAGING	237
BALL FILLING INSTRUMENTS	139 & 143-144
BALL MARKING EXPLORERS	127
BALL SCALERS	150
BAND REMOVING PLIERS	277 & 286
BAND SEATER AUTO CLEAVABLE	273
BAND SEATING & REMOVING INSTRUMENTS	273
BAND SEATING & REMOVING SILICON GRIP	273
BAND SEATING PLIERS	283
BARNHART CURETTES	154
BARRAQUER NEEDLE HOLDERS TC GOLD	241
BASINS FOR BONE	328
BAUMGARTNER NEEDLE HOLDERS TC GOLD	239
BAUMGARTNER NEEDLE HOLDERS	242
BAYONET ELEVATOR	96
BEALE CEMENT SPATULAS	206
BEEBEE CROWN SCISSORS	321
BEEBEE GUM SCISSORS TC GOLD	233
BEERENDONK MINI GAUGE	298
BEIN ROOTS ELEVATOR	69 & 72 & 91-94
BENDING PLIER	288
BENNET FILLING INSTRUMENTS	142
BERNARD BAYONET ELEVATOR	96
BERTEN ROOT ELEVATOR	69 & 92
BEYER BONE ROUNGER	265-266
BIRD BEAK PLIERS	283 & 288
BIRD'S BEAK CLAW WITH CUTTER	281
BISHOP HARMAN MICRO TISSUE FORCEP	324
BLACK TITANIUM COLOUR COATED INSTRUMENTS	300-304
BLADE AND PUSHER	316
BLUE PLASMA COATING & OPERATIVE INSTRUMENTS	162
BLUMENTHAL BONE ROUNGER	264
BOGHOSIAN IMPLANTOLOGY	301
BOHM PLIERS FOR TELESCOPE CROWNS	267

DESCRIPTION	PAGE#
BOLEY GAUGE	298
BONE BRACKET GAUGE	269
BONE CARRIER TITANIUM	314
BONE COMPACT KIT	333
BONE COMPRESSION KIT	334
BONE CRUSHER + MALLET	330
BONE CRUSHER	330
BONE CRUSHER BONE GRINDING	330
BONE CRUSHER TITANIUM METAL CUTTER SET OF 2 PCS	330
BONE CURETTES	201
BONE EXPANDER KIT	331 & 333
BONE EXPANDER SCREW & SURGICAL KITS	331
BONE FILES	199 & 256
BONE GRAFT/CARRIER TITANIUM	318
BONE GRAFTING INSTRUMENTS	315-317
BONE GRAFTING INSTRUMENTS COMPACTOR /TAMPERS	315
BONE INJECTORS, COLLECTORS & SUCTION TUBE	336
BONE MILL	328
BONE MORSELIZER MILL CRUSHER	328
BONE ROUNGER	264-266
BONE SCRAPER	328
BONE SCRAPER STRAIGHT & CURVED	328
BONE SPLITTING CHISELS	257
BONE SPOON/ PACKER	315
BONE THICKNESS GAUGE	269 & 297
BONE UNIVERSAL KIT	331
BOX FOR DENTAL IMPLANT TOOLS AND SURGICAL KITS	361
BOX WITH FILTER	354
BOYNTON NEEDLE HOLDERS	243
BRACKET HOLDER TWEEZERS	270
BRACKET REMOVING PLIER	278
BRACKET REMOVING PLIERS ANGLED	286
BRACKET REMOVING PLIERS STRAIGHT	286
BRACKET REMOVING PLIERS WITH PLASTIC TIP	286
BRACKET TWEEZERS SLOT ALIGNER	270
BRADLEY ROOT ELEVATOR	82
BREWER RUBBER DAM CLAMPS FORCEPS	244
BROWN TWEEZERS	118
BRUNINGS CHEEK RETRACTOR	223
BUCCAL DENTAL PHOTOGRAPHIC MIRRORS	217
BUCCAL/LINGUAL CURETTES	152
BUCCAL-SMALL DENTAL PHOTOGRAPHIC MIRRORS	217
BUCK PERIODONTAL FILES	176
BUCK PERIODONTAL INSTRUMENTS	177 & 196
BUFFALO WAX KNIVES	209
BUNTING CURETTES	154
BUNTING CURETTE-CHISEL COMBINATION	154
BURNISHER PERIODONTAL INSTRUMENTS	157
BURNISHERS	302
BUSCH EXTRACTING FORCEPS (ENGLISH PATTERN)	24
BUSER PERIODONTAL INSTRUMENTS	259
C	
CALIPER / CALIPERS RULE	296
CARBIDE TIPPED ROOT TIP FORCEPS	39
CARRY FORCEPS	324
CARVERS INSTRUMENTS	163-165 & 207
CASFROVIEJO NEEDLE HOLDERS TC GOLD	241
CASSETTE FOR 10 INSTRUMENTS	359
CASSETTE FOR EXTRACTION FORCEPS	359
CASSETTE FOR ROOT ELEVATOR	359
CASSETTE TRAY 15 PCS	359
CASTROVIEJO	296
CASTROVIEJO CALIPERS	296
CASTROVIEJO GUM SCISSORS	230
CASTROVIEJO MEASURING LOOP	297
CASTROVIEJO MICRO SCISSORS	234
CASTROVIEJO NEEDLE HOLDERS	242
CASTROVIEJO NEEDLE HOLDERS TC GOLD	240
CATHETERS TRAY WITH LID	354
CATTONI PROBE	126
CATTONI SET	153
CEMENT SPATULAS	173-175 & 205-206
CERAMIC COATED SCISSORS T.C	235
CERI SAW	268
CERI SAW CRUCIBLE TONGS	268
CHALAZION MUCOCLELE FORCEPS	323
CHEATLE STERILIZING & TISSUE CLAMPS	238
CHEEK RETRACTOR	222-224
CHEEK RETRACTOR RETRACTORS COMBINATION	224
CHILD BUCCAL DENTAL PHOTOGRAPHIC MIRRORS	217
CHILDREN PATTERN	22
CHIP BLOWER	254
CHISELS, BONE FILES & SPOONE CURETTE	256-257
CHISELS	165 & 167 & 257
CHOMPRET ROOT ELEVATOR	102
CHOMPRET SYNDESMOTOME	102
CINCH BACK / TIE BACK DISTAL END BENDING PLIERS	287
CLEOIDS FILLING INSTRUMENTS	165
CLEVELAND BONE ROUNGER	265
COARSE SUCTION TUBES	336
COGSWELL ELEVATOR ERGO DYNAMIC	86
COGSWELL ROOT ELEVATOR	82
COHEN TISSUE NIPPER	266
COHEN UNIVERSAL FORCEPS	53
COLEMAN ROOT ELEVATOR	80
COLLEGE TWEEZERS	115
COLLINS TONGUE FORCEPS	237

DESCRIPTION	PAGE#	DESCRIPTION	PAGE#
COLUMBIA CURETTES	153	DETAILING STEP PLIERS / BAYONET BENDING PLIERS	287
COLUMBIA PERIODONTAL INSTRUMENTS	160	DIAMOND COATED EXTRACTION FORCEPS - BETTER GRIP & TACTILE	40
COLUMBIA RETRACTORS	223	FEEL - LESS TRAUMA	252-253
COLUMBIA SET	130	DIDNER MATRIX RETAINERS	219-220
COMBINATION CURETTES	152	DIEFFENBACH SCALPEL KNIVES	270
COMELLAS FORCEPS ANATOMICAL HANDLE	67	DIRECT BOND BRACKET HOLDER TWEEZERS	270
COMPACTOR	317	DIRECT BONDING BRACKET HOLDERS	291
COMPASSO DE WILLIS ORTHODONTIC INSTRUMENTS	269	DIRECT BRACKET REMOVER TC. GOLD	262
COMPOSITE GOLD COATED INSTRUMENTS	143	DISC SYSTEM DENTAL SYRINGES	255
COMPOSITE GOLD TIP INSTRUMENTS	144	DISTAL AMALGAM INSTRUMENTS	273
COMPOSITE INSTRUMENTS	159	DISTAL BENDER SILICON GRIP	282
COMPOSITE INSTRUMENTS BLUE TITANIUM WITH CASSETTES	204	DISTAL END CUTTER FLUSH CUT T.C	290
COMPOSITE INSTRUMENTS TITANIUM BLUE ENDODONTIC	204	DISTAL END CUTTER LONG HANDLE TC. GOLD	285
COMPOSITE PACKING INSTRUMENT	303	DISTAL END CUTTER T.C	290
COMPRESSION EXPANDER SCREWS KIT	333	DISTAL END CUTTER TC. GOLD	290
COMPRESSOR DOUBLE FLASK	212	DISTAL END CUTTER, SAFETY HOLD PLIER TC. GOLD	88
CONCAVES SET OF 5 PIECES	325	DISTAL ROOT ELEVATORS	295
CONDENSER IMPLANTOLOGY	308	DOUBLE ENDED OSTEOTOMY DEPTH GAUGE	272
CONE-SOCKET FRONT SURFACE SET OF 12 PCS	112	DOUBLE ENDED TAMPERS	288
CONE-SOCKET MAGNIFYING MOUTH MIRRORS	111	DOUBLE ROUNDED JAW BENDING PLIERS	125
CONE-SOCKET MIRROR HANDLE	113	DOUBLE-ENDED EXPLORERS	196
CONE-SOCKET PLAIN SET OF 12 PCS	111	DR. ALLEN END-CUTTING INTRASULCULAR KNIFE	117
CONE-SOCKET RHODIUM SET OF 12 PCS	112	DRESSING FORCEPS	355
CONE-SOCKET SOLID	113	DRESSING JAR / COTTON CONTAINER	355
CONE-SOCKET SOLID OCTAGONAL	113	DRESSING JAR WITH LID & KNOB	332
CONICAL IMPLANT DRILL KIT	335	DRILL KIT	213
CONTOURING INSTRUMENT	301	DYERLINE FLASK, SILVER	299
CONTOURING PLIER	276 & 288	DYNAMOMETER	
CONVEXES SET OF 5 PIECES	326	E	
COSMETIC CONTOURING IMPLANTOLOGY	303	EASY ACCESS TWEEZERS	270
COTTON HOLDER	356	EHRICKE STAINLESS STEEL IMPRESSION TRAYS	273
COTTON HOLDING JAR	356	ELASTIC INSTRUMENT SILICON GRIP	271
COTTON MERIAM TWEezer	175	ELASTIC INSTRUMENTS	244
COTTON WOOL DISPENSER ACC. TO METHOT	356	ELASTIC SEPARATING & LIGATURE PLIER	292
COUPLAND GOUGE ROOT ELEVATOR	85	ELASTIC SEPARATING PLIER	292
COUPLAND ROOT ELEVATOR	73 & 101	ELASTOMERIC LIGATURE PLACING FORCEPS	281
COW HORN EXTRACTING FORCEPS	65	ELASTOMERIC SEPARATOR PLACING PLIERS	253
CRANA KAPLAN PERIODONTAL INSTRUMENTS	160 & 196	ELLIOT MATRIX RETAINERS	166
CRANE KAPLAN TWEEZERS	117	ENDED CHEEK RETRACTORS	269
CRANE PICK ROOT ELEVATOR	82	ENDO FINGER RULER SPAN MEASURE SCALE	356
CRANE-KAPLAN CURETTES	150	ENDODONTIC BOX	169
CRANE-KAPLAN PERIODONTAL INSTRUMENTS	189 & 192	ENDODONTIC EXCAVATORS	134-135
CRILE-RYDER NEEDLE HOLDERS TC GOLD	239	ENDODONTICS INSTRUMENTS	114
CRILE-WOOD NEEDLE HOLDERS	242	ENDO-LOCKING GROOVED TWEEZERS	85
CRILE-WOOD NEEDLE HOLDERS TC GOLD	239	ENGLISH PATTERNS ROOT ELEVATOR	86
CROSS-BAR HANDLE	102	ERGO DYNAMIC ELEVATORS	210
CROWN & SHELL PLIER	277	ERMERT WAX KNIVES	261
CROWN PLIER WITH PLASTIC INSERTS	267	EURO BREECH LOADING DENTAL SYRINGES	118
CROWN REMOVER	248-250	EWALD / HUDSON TWEEZERS	122-123
CROWN REMOVER TURKISH STYLE	250	EXAMINATION SET	129 & 173-175
CROWN SCISSORS	321	EXAVATORS	169-171
CROWN SPREADER	248-249	EXCAVATORS	158
CROWN SPREADERS WITH SPRING	249	EXCAVATORS INSTRUMENTS	270
CRAYER EXTRACTING FORCEPS (AMERICAN PATTERN)	50	EJECT BRACKET HOLDER TWEEZERS	124-132 & 172,185
CRAYER ROOT ELEVATOR	75-92	EXPLORERS	132
CRAYER ROOT OR FRAGMENT ELEVATOR	81	EXPLORERS INSTRUMENTS	361
CRAYER ROOTS ELEVATOR	91	EXTRA LARGE CASSETTE FOR 6 INSTRUMENTS	65-66
CUMINE CURETTES	150	EXTRACTING FORCEPS (AMERICAN PATTERN HOLE HANDLE)	45-55
CURETTE	190 & 193	EXTRACTING FORCEPS (ANATOMICAL HANDLE)	56-64
CURETTES	145-156 & 178-202 & 322	EXTRACTING FORCEPS (ENGLISH PATTERN HOLE HANDLE)	35-38
CUSHING TAYLOR TWEEZERS	116	EXTRACTING FORCEPS (ENGLISH PATTERN)	01-34 & 67
CUSHING TWEEZERS	119	EXTRACTING FORCEPS (FOR CHILDREN)	26-27 & 31-32
CUTTER TC	282	EXTRACTING FORCEPS KITS	68
CUTTER, UTILITY PLIERS & WIRE FORMING PLIERS	283	EXTRACTING FORCEPS KITS GOLD COATED SET OF 25 PCS	68
		EXTRACTING FORCEPS KITS POLISH SET OF 25 PCS	68
		EXTRACTING FORCEPS KITS SET OF 37 PCS	68
		EXTRA-FLEX IMPLANTOLOGY	300 & 304
D		F	
DE BACKY ATRAUMATIC TWEEZERS	115	FAHNENSTOCK WAX KNIVES	208
DE LA ROSA PLIER	276	FALCON BEAK CONTOURING PLIERS	275 & 288
DEAN DISSECTING SCISSORS	235	FARABEUF RETRACTOR	225
DEAN GUM SCISSORS	228	FEDI PERIODONTAL CHISEL	198
DEAN GUM SCISSORS TC GOLD	233	FELDMAN RETRACTOR	225
DEBAKEY ATRAUMATIC TISSUE FORCEPS	323	FILLING INSTRUMENTS	129-144 & 163-173 & 186 & 203-204
DE-BAKEY WITH TUNGSTEN CARBIDE (TC) TWEEZERS	116	FIRST AID EXAMINATION KIT SET OF 3 PIECES	122-123
DEEP GRIPE EXTRACTING FORCEPS	35	FIRST AID EXAMINATION SET OF 3 PIECES	122-123
DEEP GRIPPING ROOT FORCEPS	27 & 34	FIRST AID EXAMINATION SET OF 4 PIECES	122-123
DEEP ROOT PICK FORCEPS	36-38	FLAGG TWEEZERS	114
DELICATE NEEDLE HOLDERS TC GOLD	239	FLASKS & COMPRESSORS	212-213
DELUX BOLEY GAUGE	298	FLAT NOSE PLIERS	275-280
DELUXE BRACKET HOLDER	270	FLEXI-THIN IMPLANTOLOGY	300
DELUXE BRACKET HOLDER TWEEZERS	270	FLOHR ROOT ELEVATORS	73 & 77, 91-92
DENHARD RETRACTORS - MOUTH GAGS	221	FLUCH-CUT DISTAL AND CUTTING PLIERS T.C	278
DENTAL BUR CLEANING BRUSH COPPER	328	FLUSH-CUT DISTAL END CUTTER T.C	285
DENTAL FLASK 4 PARTS, BRASS	213	FOERSTER DRESSING FORCEPS	237
DENTAL FLASK, 4 PARTS SILVER	213	FORCEPS FOR PLACING DENTALASTICS	293
DENTAL MIRROR MAGNIFY	175	FORCEPS JARS	355
DENTAL PHOTOGRAPHIC MIRRORS	217	FRAHM CARVERS INSTRUMENTS	164
DENTAL SCISSORS	231	FRAZIER SUCTION TUBES	337
DENTAL SINGLE FLASK, BRASS	213	FREEDMAN "DUCKHEAD" IMPLANTOLOGY	302
DENTAL SYRINGE 3-RING	261	FREEDMAN BURNISHERS	302
DENTAL SYRINGES	261-263	FREEDMAN LARGE CONTACT FORMING	302
DENTAL SYRINGES 2-RING	261	FREER PERIODONTAL INSTRUMENTS	260
DENTAURUM RULER	299	FREER SPOONE CURETTE	256
DENTINS REMOVALS	169		
DENTURE FLASK LOWER, BRASS	212		
DENTURE FLASK UPPER, BRASS	212		
DEPTH & GINGIVAL GAUGES	295		
DERF NEEDLE HOLDERS	242		
DERF NEEDLE HOLDERS TC GOLD	239		

DESCRIPTION	PAGE#
FRET SAW	268
FRIEDMAN BONE ROUNDER	265
FRIEDMAN CRAIG GOUGE ROOT ELEVATOR	85
FRIEDMANN ROOTS ELEVATOR	76
FURRER COPPER RING REMOVING PLIER	267
G	
GARTNER ROOTS ELEVATOR	93
GERALD TWEEZERS	115-117
GILLES TWEEZERS	115
GILLIES NEEDLE HOLDER WITH SCISSORS	115
GILLIES NEEDLE HOLDER WITH SCISSORS TC GOLD	240
GINGIVAL RETRACTORS	161
GLICK BLADE / PLUGGER	135
GLOVES BOX	354
GOLDMAN CURETTES	145
GOLDMAN FOX EXPLORERS	127-128
GOLDMAN FOX PERIODONTAL INSTRUMENTS	160
GOLDMAN FOX SCISSORS	234
GOLDMAN FOX UNIVERSITY OF MICHIGAN EXPLORERS	128
GOLDMAN FOX WILLIAMS EXPLORERS	128
GOLDMAN PERIODONTAL INSTRUMENTS	188-191
GOLDMAN-FOX CURETTE SET	152
GOLDMAN-FOX CURETTES	151
GOLDMAN-FOX GUM SCISSORS	227
GOLDMAN-FOX GUM SCISSORS T.C GOLD	232
GOLDMAN-FOX PERIODONTAL INSTRUMENTS	260
GOLDMAN-FOX PERIODONTAL KNIVES	197
GOLDMAN-FOX SET	153
GOSLEE FLAT NOSE PLIERS	280
GRACEY PERIODONTAL INSTRUMENTS	160
GRACEY CURETTE (BLACK & GOLD COATING) SET OF 7 PCS	180
GRACEY CURETTES	156 & 178-179 & 180-184
GRACEY CURETTES SET OF 7 PCS	180-184
GRACEY MINI FIVE	187
GRACEY PROPHY CURETTES	156
GRADED HANDLE HYLIN CARVER	207
GRAFT MATERIAL SPOON CONCAVE TIP FOR LARGER GRAFTS	315
GRAFT PACKER	315
GREGG IMPLANTOLOGY	301
GRITMAN CEMENT SPATULAS	206
GROOVED FINE PLIER	278
GROOVED LUXATING ELEVATOR	95
GROSS MAIER DRESSING FORCEPS	237
GUM SCISSORS	226-230 & 234-235
GUM SCISSORS T.C GOLD	232-233
H	
HAEMOSTATIC FORCEPS & TOWEL CLAMPS	237
HAEMOSTATIC FORCEPS	236-237
HALSEY NEEDLE HOLDERS	242
HALSEY NEEDLE HOLDERS TC GOLD	239
HALSTEAD MOSQUITO HAEMOSTATIC FORCEPS	236
HALTEPINZETTE WITH LOCK	324
HAMMER HEAD DISTAL END BENDING PLIERS	287
HAMMER WITH WOODEN HANDLE	268
HANDLES FOR MOUTH MIRRORS	113
HARD WIRE CUTTER 15 T.C	282
HARD WIRE CUTTER STRAIGHT T.C	282
HARD WIRE CUTTER T.C	285
HARRIS EXTRACTING FORCEPS (AMERICAN PATTERN)	44
HARTMANN (BABY MOSQUITO) HAEMOSTATIC FORCEPS	236
HATCHES FILLING INSTRUMENTS	166
HEAT CARRIERS	133
HEBER BAYONET ROOTS ELEVATOR	94
HEBER ROOTS ELEVATOR	94
HEIDBRINK ROOT APPEX ELEVATOR	90
HEIDBRINK ROOT ELEVATOR	100-101
HEIDBRINK ROOT ELEVATORS	76-78 & 98-99
HEIDBRINK UNIVERSAL MINI-ELEVATOR FOR FRAGMENTS	100
HEIDEBRINK ROOT TIP PICKS	98
HEIDEMANN FILLING INSTRUMENTS	143-144
HEIDEMANN SPATULA	157 & 168
HEISS HAEMOSTATIC FORCEPS	236
HEMINGWAY BONE CURETTES	201
HEMINGWAY TWEEZERS	117
HENAHAN EXTRACTING FORCEPS (AMERICAN PATTERN)	53
HERBERT PERIODONTAL INSTRUMENTS	258
HERTZ EXTRACTING FORCEPS (ENGLISH PATTERN)	20
HEX DRIVERS ADAPTER KIT	333
HIRSCHFELD PERIODONTAL FILES	176
HIRSCHFELD PERIODONTAL INSTRUMENTS	177
HMG HEISTER CHEEK RETRACTOR	222
HOE CURETTES	152
HOE/CHISEL CURETTES	153
HOLLENBACK CARVERS INSTRUMENTS	163
HOLLENBACK FILLING INSTRUMENTS	136
HOLLENBACK IMPLANTOLOGY	304
HOLLOW CHOP CONTOURING PLIERS	284
HOLLOW MIRROR HANDLE ROUND	113
HOLLOW MIRROR HEXAGONAL HANDLE	113
HOLLOWARE INSTRUMENTS	351-360
HOOK CRIMPING PLIER ANGLED TC. GOLD.	290
HOOK KNIVES	210
HOURLIGAN ROOT ELEVATOR	85
HOW ANGLED UTILITY PLIERS	286
HOW PLIER	275 & 283
HOW PLIERS OFFSET	283

DESCRIPTION	PAGE#
HOW STRAIGHT UTILITY PLIERS	286
HOWARD BONE FILE	199
HOWARD PERIODONTAL INSTRUMENTS	259
HOWARD ROOT ELEVATOR	100
HULL EXTRACTING FORCEPS (AMERICAN PATTERN)	49
HYGIENIST CURETTES	145
HYGIENIST JACQUETTE PERIODONTAL INSTRUMENTS	189
HYGIENIST PERIODONTAL INSTRUMENTS	188 & 191
HYGIENIST-JACQUETTE CURETTES	145
HYGIENIST-TOWNER CURETTES	145
HYGIENIST-TOWNER PERIODONTAL INSTRUMENTS	188 & 191
HYLIN CARVERS INSTRUMENTS	207
HYLIN ROOT ELEVATOR	80 & 90
HYLIN-CERAMICA CARVERS INSTRUMENTS	207
HYLIN-LE CRON CARVERS INSTRUMENTS	207
IMPLANT	313
IMPLANT DEPTH GAUGE & KEY	295
IMPLANT DEPTH GAUGE	295
IMPLANT GAGES	295-338
IMPRESSION TRAY	339-352
IMPRESSION TRAY SET OF 06 PCS	341-343 & 351-353
INDIANA UNIVERSITY ANTERIOR	154
INDIANA UNIVERSITY POSTERIOR	154
INSTRUMENTS BOX	354
INSTRUMENTS BOX PERFORATED	354
INSTRUMENTS BOX WITH HOOKS	354
INSTRUMENTS BOX WITH LID	354
INSTRUMENTS CASSETTE 10 PIECES CASSETTE	361
INSTRUMENTS CASSETTE 5 PIECES CASSETTE	361
INSTRUMENTS CASSETTE 7 PIECES CASSETTE	361
INSTRUMENTS MANAGEMENT SYSTEM	359
INSTRUMENTS TRAY PERFORATED	353 & 361
INSTRUMENTS TRAY PERFORATED BASE	361
INSTRUMENTS TRAY SOLID BASE	361
INSTRUMENTS TRAY WITHOUT COVER	361
INSTRUMENTS TRAYS	353
INSTRUMENTS TRAYS WITH LID	353
INTERPROXIMAL CARVERS INSTRUMENTS	164
INTERPROXIMAL CARVER	301
INTERPROXIMAL CARVER, LONG	301
INTERPROXIMAL OFF-ANGLE CARVERS INSTRUMENTS	164
INTRALIGAMENTAL ANESTHESIA SYRINGE	262
IRIS DOUBLE BEND GUM SCISSORS T.C GOLD	232
IRIS GUM SCISSORS	226
IRIS GUM SCISSORS T.C GOLD	232
IRIS SCISSORS	234
IRONING CRESCENT	175
IVORY MATRIX RETAINERS	252
IVORY RUBBER DAM CLAMPS FORCEPS	244
IVORY WITHOUT SHIELDS MATRIX RETAINERS	253
IWANSON SPRING CALIPER FOR METAL AND PORCELAIN	297
IWANSON SPRING CALIPER FOR WAX	297
J	
JACQUETTE PERIODONTAL INSTRUMENTS	189 & 192
JACQUETTE SCALERS	146
JACQUETTE SCALERS PERIODONTAL INSTRUMENTS	189 & 192
JACQUETTE PERIODONTAL INSTRUMENTS	160
JARABAK PILERS	279 & 284
JECT DENTAL SYRINGES	262
JECT ULTIMA INJECTS	262
JENNINGS RETRACTORS - MOUTH GAGS	221
JOHNSON CONTOURING PLIERS	275
JOHNSON CURETTES	147
JOHNSON PLIERS	288
JONES HAEMOSTATIC FORCEPS & TOWEL CLAMPS	237
JOSEPH GUM SCISSORS	227-228
K	
KALT NEEDLE HOLDERS TC GOLD	241
KELLS EXTRACTING FORCEPS (AMERICAN PATTERN)	49
KELLY GUM SCISSORS	228
KELLY GUM SCISSORS TC GOLD	233
KELLY HAEMOSTATIC FORCEPS	236
KELLY-RANKIN HAEMOSTATIC FORCEPS	236
KERRISON BONE ROUNDER	264
KEYS FOR IMPLANT SCREW	329
KIDNEY TRAY	355
KIDNEY TRAY SHALLOW	355
KIDNEY TRAY WITH COVER	355
KILNER (RAGNELL) GUM SCISSORS	227
KIRKLAND PERIODONTAL INSTRUMENTS	160
KIRKLAND PERIODONTAL KNIVES	197
KOCHER RETRACTOR	225
KOEHLER LUXATION FORCEPS	31
KOPP ROOTS ELEVATOR	93-94
KRAMER NEVINS SINUS LIFT	312
KRAMER-NEVINS PERIODONTAL CHISELS	200
KRAMER-NEVINS BONE CHISEL	256
L	
LA GRANGE NEEDLE HOLDERS TC GOLD	241
LA GRANGE GUM SCISSORS	227
LA GRANGE GUM SCISSORS T.C GOLD	232
LADMORE FILLING INSTRUMENTS	140
LADMORE IMPLANTOLOGY	302
LAGRANGE SCISSORS	234

DESCRIPTION	PAGE#
LANGENBECK RETRACTOR	225
LAP & L-KEY JOINT PLIERS T.C	290-291
LARGE PLACING/ CONDENSING	304
LAXUTOR ROOT ELEVATOR	74
LAXUTOR ROOT ELEVATORS SET OF 9 PCS	87
LE CRON CARVERS INSTRUMENTS	207
LE CRON CEMENT SPATULAS	206
LECLUSE BAYONET ROOT ELEVATOR	97
LEMPERT DELICATE BONE ROUNGER	265
LESSMANN WAX KNIVES	208
LICHTENBERG-RYDER NEEDLE HOLDERS TC GOLD	241
LIGATURE CUTTERS T.C	277
LIGATURE INSTRUMENT SILICON GRIP	273
LIGATURE INSTRUMENTS	271
LIGATURE PLIERS	293
LIGATURE PLIERS T.C GOLD	293
LIGHT WORE BIRD BEAK PLIERS WITH GROOVE	284
LINDO LEVIAN ROOT ELEVATORS	76
LINDO LEVIAN EXTRACTING FORCEPS (AMERICAN PATTERN)	41
LINGUAL / PALATAL ARCH FORMING PLIES	287
LINGUAL ARCH FORMING PLIERS	284
LINGUAL DENTAL PHOTOGRAPHIC MIRRORS	217
LINGUAL LARGE DENTAL PHOTOGRAPHIC MIRRORS	217
LINGUAL LOWER DENTAL PHOTOGRAPHIC MIRRORS	217
LINGUAL UPPER DENTAL PHOTOGRAPHIC MIRRORS	217
LISTER BANDAGE SCISSORS	231
LISTON BONE ROUNGER	266
L-KEY JOINT PLIERS T.C	285-289
LOCKLIN DISSECTING SCISSORS	235
LOCKLIN GUM SCISSORS	228
LONDON COLLEGE TWEEZERS	114 & 116
LONDON HOSPITAL ROOT ELEVATOR	80
LONG IMPLANT DRILLS SET OF 4 PCS	327
LOOP FORMING PLIER	288
LORNA STERILIZING & TISSUE CLAMPS	238
LUCAS CURETTES	160 & 186,202 & 322
LUCAS SPOON SHAPE	202
LUCAS SURGICAL CURETTE	202
LUER BONE ROUNGER	264
LUKS PERIODONTAL PROBES	134
LUXACUT KIT	104
LUXACUT LUXATION INSTRUMENT	104
LUXACUT ROOT ELEVATOR	104
LUXATING ELEVATOR / PERIOTOME	96
LUXATING ELEVATOR	86
LUXATING ELEVATOR STRAIGHT	95
LUXATING ELEVATOR WITH PERIOTOME TIP	96
LUXATING ELEVATORS SET OF 08 PCS - BLUE TITANIUM COATED TIPS	106
LUXATING ELEVATORS SET OF 10 PCS - BLUE TITANIUM COATED TIPS	105
LUXATING ELEVATORS SET OF 10 PCS - GOLD COATED TIPS	106
LUXATING ELEVATORS SET OF 10 PCS	105
LUXATING ROOT ELEVATOR	86
LUXATOR ROOT ELEVATOR	71 & 79
LUXATOR/ELEVATOR CASSETTE 8 PCS	361
LUXIER ROOTS ELEVATOR	89
M	
MACMILLAN GOUGE ROOT ELEVATOR	85
MAGNIFYING MOUTH MIRRORS	111
MALLETS - PLASTER SHEARS - SAWS & BLADES	268
MANUAL ASPIRATING SYRINGE	261
MARGINAL RIDGE & EMBRASURE SHAPING INSTRUMENT	303
MARGINTRIMMERS	166
MARKLEY FILLING INSTRUMENTS	136 & 143-144
MARQUETTE FILLING INSTRUMENTS	136
MARTIN BONE CURETTES	201
MATHIEU LIGATURE PLIERS	293
MATHIEU NEEDLE HOLDER TC GOLD	240
MATHIEU NEEDLE HOLDERS	243
MATHIEU-KOCHER NEEDLE HOLDER TC GOLD	240
MATHIEU-KOCHER NEEDLE HOLDERS	243
MATRIX IVORY MATRIX RETAINERS	253
MATRIX RETAINERS	252-253
MAXI INSTRUMENTS TRAY	357
MAYO DISSECTING SCISSORS	235
MAYO GUM SCISSORS T.C GOLD	232
MAYO STANDARD OPERATING SCISSORS	229
MAYO TRAY PERFORATED	353
MAYO-HEGAR NEEDLE HOLDERS	242
MAYO-HEGAR NEEDLE HOLDERS TC GOLD	239
MC CALL SET	130 & 151
MEAD EXTRACTING FORCEPS (ENGLISH PATTERN)	19
MEAD PERIODONTAL INSTRUMENTS	258
MEASURING LOOP	297
MEASURING PROBE	295
MEAT PERIODONTAL INSTRUMENTS	260
MEDAN ROOTS ELEVATOR	91
MEDIUM PLACING/ CONDENSING	304
MEDIUM/LARGE CONTACT FORMING	304
MEISSNER EXTRACTING FORCEPS	67
MEMBRANE 3 POSITIONING IMPLANTOLOGY	308
MEMBRANE FORCEP	323
MEMBRANE TWEEZERS KUEHNE	324
MERCURY CARVERS INSTRUMENTS	164
MERIAM TWEEZERS	114-116
MERRIT A EXPLORERS	128
MERSHON BAND PUSHER	273
MESH CASSETTES 5 PIECES, 7 PIECES, 10 PIECES, 15PIECES & 20 PIECES	360

DESCRIPTION	PAGE#
MESH TRAYS 05 PCS	360
MESH TRAYS 07 PCS	360
MESH TRAYS 10 PCS	360
MESH TRAYS 20 PCS	360
MESH TRAYS	360
MESH TYPE IMPRESSION TRAY	342
MESH TYPE IMPRESSION TRAY SET OF 06 PCS	342
MESH TYPE IMPRESSION TRAY SET OF 08 PCS	342
MESIAL & DISTAL ROOT ELEVATORS	88
MESIAL ROOT ELEVATORS	88
METHIEU PLIER	293
METHOT JAR WITHOUT LID	355
METZENBAUM GUM SCISSORS	226
METZENBAUM GUM SCISSORS TC GOLD	233
METZENBAUM STANDARD SCISSORS	235
MICOR ADSON TISSUE FORCEP	323
MICRO BOLEY GAUGE	299
MICRO BOLEY GAUGE TWISTED	299
MICRO FRIEDMAN BONE ROUNGER	265
MICRO IRIS ANGULAR GUM SCISSORS	226
MICRO IRIS GUM SCISSORS	226
MICRO LUER BONE ROUNGER	264
MICRO MASQUITO STERILIZING & TISSUE CLAMPS	238
MICRO MINI CURETTES	156
MICRO SURGERY INSTRUMENT	120
MICRO-ADSON TWEEZERS	119
MICRO-CUSHING TWEEZERS	119
MICRO-MINI IMPLANTOLOGY	300 & 304
MICRO-MINI PIN & LIGATURE CUTTER T.C	282
MIDDELDORFF RETRACTOR	225
MIDI INSTRUMENTS TRAY	357
MILLER APEXO ROOT ELEVATOR	83
MILLER BONE FILES	199 & 256
MILLER CURETTES	322
MILLER ROOT ELEVATORS	75
MILLER SPOONE CURETTE	256
MILLER TWEEZERS & SPLINTER FORCEPS	121
MINI 1 IMPLANTOLOGY	300-304
MINI 3 EXTRA-FLEX IMPLANTOLOGY	300
MINI 4 EXTRA-FLEX IMPLANTOLOGY	300
MINI DISTAL AMALGAM INSTRUMENTS	255
MINI FIVE CURETTES	149
MINI LIGATURE CUTTERS T.C	277-278
MINI PIN & LIGATURE CUTTER T.C	282
MINI PIN & LIGATURE CUTTER TC.GOLD	290
MINI RETRACTOR	224
MINI SIZE COPERCHIO LID	357
MINI-LANGENBECK RETRACTOR	224
MINNESOTA RETRACTORS	223
MIRROR HANDLE AUTOCLAVABLE MATERIAL	113
MIRROR HANDLE SILICON GRIP	113
MIRROR HANDLE SOLID	113
MIRROR HANDLE SOLID OCTAGONAL	113
MISCELLANEOUS PRODUCTS	251
MISCH SOCKET SCOOP	70
MIXING SPATULA WITH SHOVEL	314
MODELLING INSTRUMENTS	174
MODIFIED MCCALL SCALER	190 & 193
MODIFIED MCCALL CURETTES	151
MODIFIED TG WITH OCHSENDEIN CUTOUTS	200
MOLARS ROOTS ELEVATOR	90
MOLT PERIODONTAL INSTRUMENTS	260
MORSE CURETTES	148
MOUTH GAG MOLT (ADULTS)	222
MOUTH GAG MOLT (CHILDREN)	222
MOUTH MIRRORS	111-112 & 129
MOUTH MIRRORS SET OF 12 PCS	111
MOUTH PROPS	251
MOVEABLE HEIGHT GAUGE	269
MULTIFUNCTIONAL METER	295
MULTI-PURPOSE PLIERS	287
N	
NABERS PROBE	126
NANCE LOOP FORMING PLIERS	276-287
NANCE PLIER TC. GOLD	291
NANCE PLIERS	276
NAPKIN HOLDER WITH METAL CHAIN	251
NAPKIN HOLDER WITH PLASTIC CHAIN	251
NEBRASKA CURETTES	149
NEEDLE HOLDERS	239-243
NEEDLE STERILIZATION TUB	356
NEUMANN GUM SCISSORS	227
NEVI ANTERIOR SCALERS	146
NEVI POSTERIOR SCALERS	146
NEVIUS EXTRACTING FORCEPS (AMERICAN PATTERN)	48-53
NOYES GUM SCISSORS	230
NYSTROM FILLING INSTRUMENTS	165
NYSTROM MATRIX RETAINERS	253
O	
OCCLUSAL ANATOMY INSTRUMENT	303
OCHSENBIN PERIODONTAL INSTRUMENTS	258
ODONT/JACQUETTE	149
OHL PERIODONTAL INSTRUMENTS	259
OLSEN-HEGAR NEEDLE HOLDER WITH SCISSORS	240
OLSEN-HEGAR NEEDLE HOLDER WITH SCISSORS TC GOLD	240
OMEGA LOOP FORMING PLIERS	283

DESCRIPTION	PAGE#
ORBAN KNIFE PERIODONTAL INSTRUMENTS	160
ORBAN PERIODONTAL KNIVES	197
ORINGER MEDIUM RETRACTORS	223
ORTHO BAND PUSHER SCALER	272
ORTHO INSTRUMENT	271
ORTHO LIGATURE DIRECTOR SINGLE ENDED	272
ORTHO LIGHT BAND PUSHER SCALER	273
ORTHO LIGHT DOUBLE ENDED LIGATURE DIRECTOR	272
ORTHO LIGHT DOUBLE ENDED LIGATURE DIRECTOR SCALER	273
ORTHO LIGHT LIGATURE DIRECTOR DOUBLE ENDED	272
ORTHO LIGHT SCALER	271
ORTHO PLIERS	289
ORTHO RULER / CONE SOCKET	299
ORTHODONTAL PERFORATED SET OF 06 PCS	345
ORTHODONTAL PERFORATED SET OF 08 PCS	345
ORTHODONTAL SOLID SET OF 06 PCS	345
ORTHODONTAL SOLID SET OF 08 PCS	345
ORTHODONTIC INSTRUMENTS	269-273
ORTHOLITE LIGATURE DIRECTOR SCALER	272
OSSEODENSIFICATION BURS KITS	332
OSTEOMETER CLAMP	297
OSTEOMETER MEASURING FORCEPS	269
OSTEOTOME CONVEXES / CONCAVE / CONIQUES	325-326
OVAL BONE CURETTES	201
OVAL SHAPE EXCAVATORS	170
OVAL SPOON	169
P	
P.K. THOMAS WAX CARVERS SET OF 5 PCS	203
P.K. THOMAS WAXING INSTRUMENTS	203
PALATAL-ADULT DENTAL PHOTOGRAPHIC MIRRORS	217
PALATAL-ADULT LARGE DENTAL PHOTOGRAPHIC MIRRORS	217
PALATAL-CHILD DENTAL PHOTOGRAPHIC MIRRORS	217
PALMER RUBBER DAM CLAMPS FORCEPS	244
PARTLY EXTRACTING FORCEPS (AMERICAN PATTERN)	46
PARTIAL PERFORATED SET OF 7 PCS	348
PARTIAL SOLID SET OF 5 PCS	348
PARTSON PERIODONTAL INSTRUMENTS	258
PDL ELEVATOR SET OF 7 PCS	108
PDL ELEVATOR SET OF 8 PCS	110
PDL ELEVATORS SET OF 10 PCS	109
PEDIATRIC TOOTH EXTRACTING FORCEPS	26
PEDO FORCEPS	40
PEESO WIRE BENDING PLIERS	280
PEET SPLINTER FORCEPS TC GOLD	241
PEET TWEEZERS & SPLINTER FORCEPS	121
PERFORATED FULL DENTURE SET OF 06 PCS	340-344
PERFORATED FULL DENTURE SET OF 08 PCS	340-344
PERFORATED FULL DENTURE UPPER	339
PERFORATED FULL DENTURE UPPER SET OF 06 PCS	339
PERFORATED FULL DENTURE UPPER SET OF 08 PCS	339
PERFORATED RIM LOCK UPPER IMPRESSION TRAY	338
PERFORATED WITHOUT RIM LOCK	346
PERFORATED WITHOUT RIM LOCK SET OF 06 PCS	346
PERFORATED WITHOUT RIM LOCK SET OF 08 PCS	346
PERFORATED WITHOUT RIM LOCK SET OF 10 PCS	346
PERFORATED WITHOUT RIM LOCK SET OF 12 PCS	346
PERIODONTAL CHISELS	198 & 200
PERIODONTAL FILES	176-177
PERIODONTAL INSTRUMENTS	129-130 & 157-168 & 173-177 & 188-200
PERIODONTAL INSTRUMENTS SET	175
PERIODONTAL INSTRUMENTS SET OF 10 PIECES	175
PERIODONTAL KNIVES	196-197
PERIODONTAL PROBES	133-134-135
PERIODONTAL SET OF 10 PIECES	129
PERIODONTAL SURGERY	312
PERIOSTEALS ELEVATOR	319-320-321
PERIOTOME GOLD TIP ROOT ELEVATOR	107
PERIOTOME LUXATING ELEVATORS SET OF 04 PCS - GOLD COATED TIPS	107
PERIOTOMES ROOT ELEVATOR	98 & 101
PERRY TWEEZERS	114-116
PHYSICK EXTRACTING FORCEPS (AMERICAN PATTERN)	42
PIN & LIGATURE CUTTER T.C.	282
PIN & LIGATURE CUTTERS TC. GOLD.	290
PIN HOLDING PLIER	291
PIN-TWEEZERS	324
PISTOL-JECT INJECTS	262
PLAIN MOUTH MIRRORS	111
PLASTER KNIFE	209
PLASTER KNIVES	209
PLASTER SHEARS	268
PLASTIC FILLING INSTRUMENTS	141-144
PLIER STAND PLASTIC	294
PLIER STAND STAINLESS STEEL	294
PLIERS FOR ORTHODONTICS & PROSTHETICS	274-294
PLIERS STAND STAINLESS ROD	294
PLIERS TO PLACE MODULES	293
PLUGGERS	175,132 & 317
POINTED CURETTE	151,190 & 193
POLOKOFF BONE FILE	199
POSIGRIP ARTICULATING PAPER FORCEPS	121
POSITIVE POSITIONER	269
POSTERIOR	148
POSTERIOR BAND REMOVING PLIERS - LONG	283
POSTERIOR BAND REMOVING PLIERS - SHORT	283
POSTERIOR JACQUETTE PERIODONTAL INSTRUMENTS	189 & 192
POSTERIOR SCALERS	146
POSTERIOR SINGLE ENDED	135

DESCRIPTION	PAGE#
POTT ROOTS ELEVATOR	97
POTTS "T-HANDLE" ELEVATOR-LEFT	70
POTTS "T-HANDLE" ELEVATOR-RIGHT	70
POTTS ELEVATOR-LEFT	70
POTTS ELEVATOR-RIGHT	70
POTTS ROOT ELEVATOR	83
POTTS-SMITH DRESSING FORCEPS TWEEZERS	116
POTTS-SMITH TISSUE FORCEPS 1X2 TEETH TWEEZERS	116
PRICHARD PERIODONTAL INSTRUMENTS	259
PROBE	126 & 172-175
PROJECTION SYSTEM DENTAL SYRINGES	262
PRONG 3 CLASP ADJUSTABLE PLIER TC. GOLD.	291
PUNCH BONE ROUNGER	266
PUNCHED TYPE IMPRESSION TRAY SET OF 06 PCS	342
PUNCHED TYPE IMPRESSION TRAY SET OF 08 PCS	342
Q	
QUATRO DENTAL SYRINGES	261
R	
RAMUS RETRACTOR	225
RASP. MOLT	321
RASPATOR, BENQUE	321
RASPATOR, BONE SCRAPER	319
RASPATOR, FREER	319,320
RASPATOR, MICRO FREER	320
RASPATOR, MICRO RASPATOR / CHISEL	320
RASPATOR, MICRO-FREER	319
RASPATOR, MOLT	320
RASPATOR, NORDBERG	321
RASPATOR, PAPILLA	319
RASPATOR, PRICHARD	314
RASPATOR, RASPATORY	259 & 319
RASPATOR, WAUTIER	321
RASPATOR, WILLIGER	259
RASPATORY CHIPPING WORKERS	321
RASPATORY MICRO MOLT	321
RASPATORY MOLT WORK	321
RECTANGULAR ARCH BENDING PLIERS	284
REDDY THIN PERIODONTAL INSTRUMENTS	196
REMINGTON CURETTES	148
REMINGTON HOE PERIODONTAL CHISELS	200
RETRACTOR AND TONGUE PUSHERS BISHOP	223
RETRACTOR COMBINATION CHEEK RETRACTOR	223
RETRACTOR WITH ERGO GRIP	222
RETRACTORS - MOUTH GAGS	221-222
RETRACTORS	223-225
RETRO FILLING PLUGGERS	135
RHODES BACK-ACTION PERIODONTAL CHISELS	200
RHODIUM FRONT SURFACE SET OF 12 PCS	112
RIDGE MAPPING DOUBLE CALIPER	296
ROACH CEMENT SPATULAS	206
ROCHESTER-KOCHER HAEMOSTATIC FORCEPS	236
ROCHESTER-PEAN HAEMOSTATIC FORCEPS	236
ROGGE STERILIZING & TISSUE CLAMPS	238
ROMEROWSKI FILLING INSTRUMENTS	140
ROOT CANAL PLUGGER	132
ROOT CANAL SPREADER	185
ROOT ELEVATOR	69-109
ROOT FRAGMENT FORCEPS	28-29
ROOT PLANNING CHISELS & BONE FILES	256
ROOTS SPLINTER FORCEPS	30
ROSER-KONIG RETRACTORS - MOUTH GAGS	221
ROUND BONE CURETTES	201
ROUND BOWLS	355
ROUND DENTAL FLASK, BRASS	212
ROUND DENTAL FLASK, SILVER	213
ROUND DRUMS WITH COVER	358
ROUND DRUMS WITHOUT COVER	358
ROUND ENDODONTIC BOXES	356
ROUND STERILIZER DRUMS	358
ROUND STERILIZER DRUMS "FILTER"	358
ROUTURIER WISDOM EXTRACTING FORCEPS	23
RUBBER DAM CLAMPS	247
RUBBER DAM CLAMPS FORCEPS	244
RUBBER DAM CLAMPS KITS	245-246
RUBBER DAM CLAMPS KITS WITH 08 CLAMPS	246
RUBBER DAM CLAMPS KITS WITH 11 CLAMPS	246
RUBBER DAM CLAMPS KITS WITH 12 CLAMPS	246
RUBBER DAM CLAMPS KITS WITHOUT CLAMPS	246
RUBBER DAM FRAME	244
RUBBER DAM KITS	245
RUBBER DAM TEMPLATE	245
RULER	299
RULER WEDGE	299
RUSKIN BONE ROUNGER	266
RUSSIAN LONG TISSUE FORCEPS	324
RUSSIAN TISSUE FORCEPS	324
S	
SACHS WAX KNIVES	210
SANDERS PERIODONTAL INSTRUMENTS	196
SAW BLADES	268
SCALER / CURETTE SETS	151-153
SCALER CURETTE	190 & 193
SCALERS	146 & 190,193-195
SCALERS TRAY	353
SCALPEL BLADE & SCALPEL KNIVES	219

DESCRIPTION	PAGE#
SCALPEL HANDLE 6 WAY ADJUSTABLE ANGLED	218
SCALPEL HANDLES	218
SCALPEL KNIVES	219-220
SCHLEMMER ROOTS ELEVATOR	90
SCHLUGER PERIODONTAL FILES	176
SCHLUGER PERIODONTAL INSTRUMENTS	177
SCHMECKEBIER APEXO ROOT ELEVATOR	70&73&83
SCHURE BAND SEATER	271
SCHWARZ PLIER	280
SELDIN PERIODONTAL INSTRUMENTS	259
SELDIN RETRACTOR	225
SELDIN ROOT ELEVATOR	71-72-77-82-86-78-88
SELDIN ROOT OR FRAGMENT ELEVATOR	82
SELDIN ROOTS ELEVATOR	91
SELF LOCKING GROOVED	323
SELF-LOCKING TWEEZERS	114
SEMKEN TAYLOR TWEEZERS	114
SENN MILLER RETRACTOR	224
SEPARATING FORCEPS	42
SEPARATING PLIERS FOR SEPARATORS	292
SHORT IMPLANT DRILLS SET OF 3 PCS	327
SICKEL-SCALER POSTERIOR CURETTES	147
SICKLE SCALERS	145-148
SICKLE SCALERS INSTRUMENTS	188 & 191
SILICONE HANDLES	185
SINGLE ENDED PROBE HOLLOW HANDLE	124
SINGLE ENDED PROBE SOLID HANDLE	124
SINUS IMPLANTOLOGY	307
SINUS LIFT ELEVATOR	305-306-309
SINUS LIFT ELEVATORS SET OF 6 PCS	311
SINUS LIFT INSTRUMENT	309-311
SINUS LIFT INSTRUMENTS BLUE PLASMA COATING	161
SINUS WINDOW BONE TAMPERS	317
SINUSLIFT IMPLANTOLOGY	308
SLIM WEINGART PLIERS	283
SLIPPER TWEEZERS	270
SMAHA ELASTOMERIC LIGATURE PLACING FORCEPS	292
SMALL FLASK WITH COMPRESSOR, BRASS	212
SOCHSENBEIN PERIODONTAL CHISEL	198
SOFT WIRE AND LIGATURE CUTTER T.C	285
SOLDERING PLIERS - TWEEZERS & SPLINTER FORCEPS	121
SOLDERING PLIERS	121
SOLID FULL DENTURE SET OF 06 PCS	340
SOLID FULL DENTURE SET OF 06 PCS	344
SOLID FULL DENTURE SET OF 08 PCS	340
SOLID FULL DENTURE SET OF 08 PCS	344
SOLID FULL DENTURE UPPER	339
SOLID FULL DENTURE UPPER SET OF 06 PCS	339
SOLID FULL DENTURE UPPER SET OF 08 PCS	339
SOLID HANDLE EXPLORERS	124
SOLID RIM LOCK UPPER IMPRESSION TRAY	338
SOLID WITHOUT RIM LOCK	347
SOLID WITHOUT RIM LOCK SET OF 06 PCS	347
SOLID WITHOUT RIM LOCK SET OF 08 PCS	347
SOLID WITHOUT RIM LOCK SET OF 10 PCS	347
SOLID WITHOUT RIM LOCK SET OF 12 PCS	347
SOLT PERIODONTAL CHISELS	200
SOLT PERIODONTAL INSTRUMENTS	196
SPATULAS	168
SPECIAL INOX WAX KNIVES	210
SPENCER GUM SCISSORS	228
SPIRIT LAMP	356
SPLINTER FORCEPS	30
SPOON AND PUSHER (TRINOVO)	318
SPOON AND PUSHER	318
SPOON AND PUSHER WITH SAFETY STOP	318
SPOON EXCAVATORS	169-171
SPREADERS	133
SPRINGFIELD ANTERIOR CURETTES	147
SPRINGFIELD CURETTES	147
SPUTUM CUP	355
STANDARD ARTICULATING PAPER FORCEPS	121
STANDARD EXTRACTING FORCEPS (AMERICAN PATTERN)	53
STANDARD LIGATURE CUTTER T.C. GOLD	290
STANDARD OPERATING SCISSORS	229
STERIL SCALPEL BLADES	219
STERILIZING & TISSUE CLAMPS	238
STERILIZING BOXES & INSTRUMENTS TRAYS	354
STERILIZING MESH TRAY	360
STERNBERG RETRACTORS	223
STIEGLITZ ROOT FRAGMENTS	30
STIEGLITZ SPLINTER FORCEPS	30
STIKSA EXTRACTING FORCEPS (ENGLISH PATTERN)	20
STOCKES RUBBER DAM CLAMPS FORCEPS	244
STOLLEY IMPRESSION TRAY	349
STOP PLIERS V-BEED	287
STRAIGHT DIRECT BOND BRACKET REMOVER T.C	291
SUCTION CANNULA KIT WITH 3 NEEDLES	337
SUCTION NEEDLE	337
SUCTION PIPE	337
SUCTION TUBES	336-337
SUGARMAN FILE	176
SUGARMAN FILE SET	176-177
SUPER CUT GOLDMAN FOX SCISSORS	234
SUPER CUT IRIS SCISSORS	234
SUPER CUT LAGRANGE SCISSORS	234
SURGICAL BALL HOOK CRIMPING PLIERS	287
SURGICAL BALL HOOK CRIMPING PLIERS, STRAIGHT	287

DESCRIPTION	PAGE#
SURGICAL CHISELS	312-313
SURGICAL DRILL KIT	334
SUTURE FORCEPS	119
SUTURE PLIERS SUTURE FORCEPS	323
SWEDISH TYPE DENTAL SYRINGES	262
SWING EXTRACTING FORCEPS (AMERICAN PATTERN)	52
SYNDESMOTOMES ROOTS ELEVATOR	98
T	
TAMPER FOR BONE MATERIAL	316
TANNER FILLING INSTRUMENTS	136
TANNER IMPLANTOLOGY	304
TARTAR REMOVER INSTRUMENTS	195
TAYLOR CURETTES	148
TAYLOR TWEEZERS	115-119
TC GRIPPER CARBIDE TIPPED TISSUE FORCEPS	323
TEFLON INTERCHANGEABLE INSERTS	258
THERMO FORMING PLIERS	289
THICK GROOVED PLIER	281
THREE JAW PLIERS	284
TISSUE FORCEPS	117
TISSUE GRAFT FORCEP	323
TISSUE MEMBRANE FORCEPS	323
TISSUE PUNCH	328
TISSUE PUNCH FORCEPS	324
TISSUE PUNCH SET OF 5 PCS	327
TISSUE PUNCH STRAIGHT & CURVED	328
TISSUE PUNCH WITH PIN SET OF 4 PCS	327
TITANIUM GRAFT PACKER/PLACEMENT SITE DEPTH GAUGE	315
TITANIUM METAL LOCKING PLIER	323
TITANIUM PLASMA COATED POINT	313
TITANIUM PLASMA COLOUR COATED IMPLANT SINUS LIFT INSTRUMENTS	307
TITANIUM PLASMA COLOUR COATED INSTRUMENTS	309-310-311-314
TITANIUM PLASMA COLOUR COATED SINUS LIFT INSTRUMENTS	305-306-308
TOFFLEMIRE RUBBER DAM CLAMPS	248
TOFFLEMIRE-JUNIOR MATRIX RETAINERS	252-253
TOFFLEMIRE-SENIOR MATRIX RETAINERS	252
TOFFLEMIRE-UNIVERSAL MATRIX RETAINERS	252-253
TOMAS EXTRACTING FORCEPS (AMERICAN PATTERN)	48
TOMAS EXTRACTING FORCEPS (FOR CHILDREN)	32
TOMATI BONE ROUNDER	266
TORQUE WRENCH & RATCHET	329
TORQUING PLIERS	287
TOWNER CURETTES	145
TOWNER PERIODONTAL INSTRUMENTS	188-191
TOWNER-HYGIENIST CURETTES	145
TOWNER-JACQUETTE PERIODONTAL INSTRUMENTS	188-191
TRELAT PERIODONTAL INSTRUMENTS	259
TREPINE	327
TREPINE DRILL SET OF 8 PCS	327
TRIMMER KIT	334
TRIO DENTAL SYRINGES	261
TUNNELIERUNGS INSTRUMENT	316
TWEED / ANGEL RECTANGULAR ARCH FORMING AND TORQUING PLIERS	287
TWEED / ANGLE RECTANGULAR ARCH FORMING & TORQUING PLIERS	279
TWEED LOOP FORMING PLIERS	276-284
TWEED LOOP PLIER	277
TWEEZERS	114-120
TWIEL-ON ELASTIC MODULE INSTRUMENT	271
TYPE A DENTAL SYRINGES	261
TYPE C DENTAL SYRINGES	261
TYPE EURO DENTAL SYRINGES	261
U	
UNIQUE RETRACTOR	224
UNIVERSAL CROWN SCISSORS	321
UNIVERSAL CURETTES	154
UNIVERSAL CUT & HOLD DISTAL END CUTTER T.C	282
UNIVERSAL DOLPHIN STYLE WIRE BENDING & CUTTING PLIERS	280
UNIVERSAL EXTRACTING FORCEPS (FOR CHILDREN)	32
UNIVERSAL GUM SCISSORS TC GOLD	233
UNIVERSAL KIT 50 PCS	335
UNIVERSAL MINI-ELEVATOR FOR FRAGMENTS	101
UNIVERSAL PLIERS	267-279-281
UNIVERSAL SCALER	152
UNIVERSITY OF MICHIGAN EXPLORERS	127-128
UNIVERSITY OF NORTH CAROLINA EXPLORERS	127
USC TOWNER PERIODONTAL INSTRUMENTS	196
UTILITY PLIER WEINGART T.C. GOLD	291
V	
V BENDING PLIERS	279
VEHE CEMENT SPATULAS	206
VERNIER GAUGE	298
VICKERS PERIODONTAL INSTRUMENTS	258
W	
W. BARRY ROOTS ELEVATOR	97
WAHL BONE FILES	199-256
WAKAIHEAT CARRIERS	133
WAKEFIELD PERIODONTAL INSTRUMENTS	258
WALDSACHS WIRE BENDING AND CUTTING PLIERS	280
WANGEN HOLDER WANGEN HOLDER WITH DUAL	222
WANGEN HOLDER	222
WARD CARVERS INSTRUMENTS	163
WARD FILLING INSTRUMENTS	138
WARWICK JAMES ROOT ELEVATOR	99-101
WATER SYRINGE	254

DESCRIPTION	PAGE#
WAX AND MODELLING CARVERS	206
WAX KNIVES	208-209-210
WAX SPATULAS	211
WAX SPATULAS WOODEN HANDLE	211
WEDELSTAEDT FILLING INSTRUMENTS	185
WEINGART PLIER	279-283
WEINGART UTILITY PLIERS	286
WEINGART-MINI UTILITY PLIERS	286
WEITLANDER RETRACTORS - MOUTH GAGS	221
WEITLANER RETRACTORS - MOUTH GAGS	221
WESTCOTT FILLING INSTRUMENTS	143-140-144
WESTCOTT SCISSORS	230
WESTON CEMENT SPATULAS	205
WHITE FILLING INSTRUMENTS	142
WHITE ROOT ELEVATOR	69-92
WHITESIDE CURETTES	149
WILAND CARVERS INSTRUMENTS	164
WILKINS/TUFTS EXPLORER	125
WILLIAM PROBE SET DOUBLE ENDED	131
WILLIAMS EXPLORERS	127-128
WILLIAMS PROBE	126
WILLIGER CURETTE	256
WINTER EXTRACTING FORCEPS (AMERICAN PATTERN)	42
WINTER ROOT ELEVATOR	197-103
WIRE BENDING AND LOOP FORMING PLIERS	286
WISDOM EXTRACTING FORCEPS	23
WISDOM ROOTS ELEVATOR	89
WITZEL EXTRACTING FORCEPS (AMERICAN PATTERN)	41
WITZEL EXTRACTING FORCEPS (ENGLISH PATTERN)	24
WITZEL FORCEPS WITH GOLD SET OF 3 PCS	67
WITZEL FORCEPS WITH HOLE AND SCREW PATTERN SET OF 3 PCS	67
WITZEL ROOT FRAGMENT FORCEPS	28
WOODSON FILLING INSTRUMENTS	142
WOODWARD EXTRACTING FORCEPS (AMERICAN PATTERN)	54
WOODWARD HAND PRESSURE CHISEL	85
WOODWARD ROOT ELEVATOR	84
WRAP AROUND EXTRACTION FORCEPS	39
WRAP AROUND FORCEPS	39
WYNMAN CROWN GRIPPER	248
X	
X RAY FILM HANGERS	251
Y	
YOUNG TONGUE FORCEPS	237
YOUNG WIRE BENDING & LOOP FORMING PLIERS	277
YOUNGER-GOOD CURETTES	149
YOUNGER-GOOD UNIVERSAL CURETTES	154
Z	
ZAHLE CEMENT SPATULAS	206

DESCRIPTION	PAGE#
-------------	-------



SIMZ MEDICAL SUPPLIES LLC

Milwaukee, Wisconsin, USA



+1 (414) 514-3542

+1 (732) 513-4737



info@simzmedicalsupplies.com



www.simzmedicalsupplies.com